



Qaelum



DOSE

PATIENT RADIATION DOSE MONITORING SOFTWARE

USER MANUAL V24

KINGDOM OF SAUDI ARABIA

QAEUM NV
Kolonel Begaultlaan 1B
3012 Leuven (Wilsela)
Belgium

QAEUM INC.
530 7th Avenue
Suite 902
New York, NY 10018
United States

Website: www.qaelum.com
E-Mail: info@qaelum.com
Phone: +32 (0) 16 585157,
+1 347 579 5756
VAT: 0840.884.387
Bank account: BNP Paribas Fortis
IBAN: BE79 0017 8626 8033
BIC: GEBABEBB

CONTENT

1.	First access to the system.....	8
1.1.	Login.....	8
1.2.	Create an account	8
1.3.	Edit account information	9
2.	User management	11
2.1.	Introduction	11
2.2.	User management v2.....	11
2.2.1.	Users tab	12
2.2.2.	Roles tab.....	13
2.2.3.	Teams tab.....	16
2.3.	User management [deprecated].....	19
2.4.	User overview	19
2.5.	Team management.....	20
2.6.	Team structure	20
2.7.	Create / Edit a team	22
2.8.	Create a filter	22
2.8.1.	Edit a filter	23
2.9.	User-Team.....	24
2.10.	Role Management	25
2.10.1.	Assigning and removing user roles.....	25
2.11.	Restrictions	27
2.12.	Default roles and their functionalities.....	29
2.12.1.	Admin	29
2.12.2.	Radiologist.....	29
2.12.3.	Radiographer	29
2.12.4.	Physicist.....	30
2.12.5.	Quality Manager	30
3.	Quick overview	31
4.	Start portal.....	33
5.	Study level.....	34
5.1.	Study details page	36
5.1.1.	General study information (Overview)	36

- 5.1.2. Series information41
- 5.1.3. Peak Skin Dose (PSD).....45
- 5.1.4. Organ doses.....55
- 5.1.5. History56
- 5.1.6. Activity stream56
- 5.1.7. Actions58
- 6. Patient level59
 - 6.1. Patient radiation dose passport.....59
 - 6.2. Study overview60
 - 6.2.1. Color coding of studies61
- 7. Notifications.....62
 - 7.1. Notifications from Study Details62
 - 7.2. Notifications Module.....64
 - 7.3. Notification Center69
 - 7.3.1. Description of rules74
 - 7.3.2. E-mail notifications94
 - 7.3.3. Import and export of notification rules95
- 8. Device level100
 - 8.1. Site section.....100
 - 8.1.1. Add site.....100
 - 8.1.2. Edit site102
 - 8.1.3. Remove site.....102
 - 8.1.4. Dashboard.....102
 - 8.1.5. Device list102
 - 8.1.6. Dynamic reports102
 - 8.2. Device section.....106
 - 8.2.1. Add Devices.....106
 - 8.2.2. Add/Edit station name107
 - 8.2.3. Edit device107
 - 8.2.4. Remove device110
 - 8.3. Configuration.....110
 - 8.3.1. DICOM configuration.....110
 - 8.3.2. Correction factors/parameters.....112

- 8.3.3. Bodypart mapping.....115
- 8.4. Device dashboards.....121
 - 8.4.1. Import of dashboards126
 - 8.4.2. Export127
 - 8.4.3. Advanced filter129
 - 8.4.4. Manual entry of an XA study132
 - 8.4.5. Dosimetry overview.....133
 - 8.4.6. Tools dashboard.....136
 - 8.4.7. Export/Save data146
 - 8.4.8. Favorite dashboard146
 - 8.4.9. Compliance monitoring147
 - 8.4.10. Optimization Gain Prediction widget.....169
- 8.5. Optional features170
 - 8.5.1. Reject analysis.....170
 - 8.5.2. Contrast Administration Monitoring.....172
 - 8.5.3. Virtual Dosimetry.....174
 - 8.5.4. Global Noise Level180
 - 8.5.5. Nuclear Medicine.....183
- 9. Live dashboard201
 - 9.1. Notifications from the Live Dashboard203
- 10. Modality level.....205
 - 10.1. Modality comparison205
 - 10.1.1. List of configurations.....208
 - 10.2. Modality Export215
 - 10.2.1. Modality export.....215
 - 10.2.2. SSM Export217
 - 10.2.3. Full export218
- 11. Location level.....220
- 12. Logbook221
 - 12.1. Device level222
 - 12.2. Hospital level226
 - 12.3. Configuration of e-mail sending in Logbook.....226
- 13. Webservice integration.....228
 - 13.1. PACS: Patient Dose passport.....228

- 13.2. Proactive warning.....231
- 13.3. Clinical quality evaluation.....232
- 14. Settings.....238
 - 14.1. Warnings238
 - 14.1.1. Unconfigured data.....238
 - 14.1.2. Scheduler240
 - 14.2. Engine monitor241
 - 14.2.1. Engine Monitor241
 - 14.2.2. Study analysis241
 - 14.2.3. Log file browser241
 - 14.3. Data tables.....242
 - 14.3.1. Conversion factors242
 - 14.3.2. Study Groups247
 - 14.3.3. Dicom bodypart mapping.....248
 - 14.4. Data management248
 - 14.4.1. Remove data248
 - 14.4.2. Analysis configuration.....248
 - 14.4.3. Notification Recalculation Settings.....248
 - 14.5. Tools250
 - 14.5.1. User Action log.....250
 - 14.5.2. DICOM Tools251
 - 14.5.3. Notification configuration.....252
- 15. Integrations.....255
- 16. About.....256
 - 16.1. Online Training Center.....257
 - 16.1.1. Sign in258
 - 16.1.2. How to find a video259
 - 16.2. Service desk262
 - 16.2.1. Invitation for Service Desk262
 - 16.2.2. Configure your login details.....263
 - 16.2.3. Reporting via Service Desk.....263
 - 16.2.4. All requests overview269
 - 16.2.5. Profile269

16.3. Symbols.....270

17. Supported browsers271

18. Minimal hardware requirements.....271

19. Functional description272

 19.1. Background.....272

 19.2. Description of the solution.....272

 19.3. Supported modalities273

 19.4. DOSE: Individual patient benefits273

 19.5. DOSE residual risks275

 19.6. DOSE intended user and training.....276

 19.6.1. Intended Use.....276

 19.6.2. Medical conditions277

 19.6.3. Contraindications278

 19.6.4. Performance of the device.....278

20. Incident reporting281



21. Glossary281

SOFTWARE IDENTIFICATION

Project ID
Customer project ID

Software identification
DOSE - Patient radiation dose monitoring software
VERSION: v24

Please read the instructions for use carefully before using this application.

 [User manual](#)
 [Technical requirements](#)





Link to Online Instructions for Use
<https://qaelum.com/media-center/documentation>

To access our on-line Training Center and watch educational videos, please visit:
knowledge.qaelum.com

To report a bug or to give feedback, please use our ServiceDesk:
Qaelum ServiceDesk
support@qaelum.com

CE 1639

[Cancel](#)

Qaelum NV
Kolonel Begaultlaan 1B
3012 Leuven (Wisele)
Belgium

2024-03-25
(YYYY-MM-DD)

HRI =
(01)05430003132
012(8012)v24

USER MANUAL

DOSE v24

DOCUMENT NAME	DOCUMENT DESCRIPTION	DOCUMENT REVISION NUMBER
Qaelum - DOSE V24 - User Manual – EN - KSA	User manual v24 – English version - KSA	Version 1

GENERAL INFO	
Target regulatory jurisdiction:	KINGDOM OF SAUDI ARABIA
Notice to the user and/or patient:	Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the Saudi Food and Drug Authority


Date of issue: 25MAR2024 (initial version)

	WHO	Application date
Prepared by:	Mahta Mazloumi, Robin Stapleford and Ana Dolcet	08MAR2024
Revised by:	Tom Van Herpe	12MAR2024
Translated by:	N/A	N/A
Controlled by:	N/A	N/A
Approved by:	Jurgen Jacobs	13MAR2024

1. First access to the system

1.1. Login


The application can only be accessed by users who are logged in. The user is asked to input his/her login credentials on the start page. New accounts can be created by clicking on **Create account** below the blue 'Log in' button. The user can select the preferred language from a list of languages as defined in the relevant regulations by clicking on the language button in the top right corner of the page. The background of the login page can be customized if the user provides Qaelum with the desired image.



Login page

1.2. Create an account

A user can create an account by clicking on the **Create account** button. A guest account will be created, and the administrator (this can be Qaelum or a user with the appropriate role) can assign new users the right team, role and permissions. The user can select the desired language.



Create account

Enter all account details

Username

First name

Last name

E-mail address

Password

Repeat password

Language
 English


Cancel Create account

Create account

Administrators can also create new user accounts inside the DOSE application by clicking on **Actions->Register Account** in **Settings/User Management**.


1.3. Edit account information

After logging in to the application, a user can change the account information and password by clicking on the **Profile** button, situated at the top right corner of the page.



Profile button


The user can select the default Start Portal setting and choose between DOSE and COACH (if available). Additionally, the start screen preference 'Quick overview', 'Expert mode', 'Live Dashboard', or 'Live Dashboard Modality' can be selected here.



The screenshot shows a 'Profile editor' dialog box with a close button (X) in the top right corner. Below the title bar, there is a user icon and the text 'Configure the parameters of this profile' followed by a blue link 'Update password'. The main area contains several input fields and dropdown menus: 'Real Name' with the value 'qaelum', 'E-mail address' with 'support@qaelum.com', 'Phone number' (empty), 'Language' with a dropdown showing 'en', 'Start portal preference' with a dropdown showing 'DOSE', and 'Start screen preference' with a dropdown showing 'Expert mode'. At the bottom, there are two buttons: a red 'Cancel' button and a blue 'Update' button. A red warning message at the bottom reads: '* This application will be refreshed when you update the parameters of your profile'.

Edit profile


We strongly recommend users to change passwords immediately in cases when the credentials are created by an Administrator. For safety and privacy reasons it is vital to choose a strong password and to change it frequently.




The screenshot shows a 'Change Password' dialog box with a close button (X) in the top right corner. It contains three input fields: 'Current Password', 'New Password', and 'Confirm New Password'. At the bottom, there are two buttons: a blue 'Commit' button and a red 'Discard' button.

Change password

It is also possible to change the color of the background. By clicking on **Switch Lights**, the background will change from dark to light. The dark background is the default setting, and is more user friendly.



Switch lights



The light color version of DOSE

2. User management

2.1. Introduction

It is possible to manage access to specific sites and devices by using **Teams** and it is possible to manage the functionalities the user can access by using **Roles**. This chapter will explain how to manage teams and roles within DOSE.


From v24, there is a new interface for user management called "User management v2" that will coexist with the previous "User management" for some releases, but it will eventually replace it. For new users we recommend using directly "User management v2".

2.2. User management v2


In the **Settings** of DOSE, a section called **User Management** can be found by any user with "User management" rights. Upon clicking on it, 3 tabs can be accessed: Users, Roles and Teams.

2.2.1. Users tab


On the first tab, the user has the possibility to create new users or to manage the existing ones.



By clicking on a single user in the table, all details of the user are shown. The details are greyed up as this is a read-only function.



By clicking on the user edit button (the blue button on the right end of the table), the details become editable.



In edit mode, it is possible to edit the user details as well as assign the user to teams and roles.

In the current version, it is not possible to edit multiple users at the same time in the user tab, however, it is possible to assign multiple users to a role/team in the following tabs.

2.2.2. Roles tab

The interface of the Roles tab is coherent with the interface of the Users tab. By clicking on a single role in the table, all details of the role are shown. It is possible to search for a user belonging to the role in this read-only mode.

Roles Overview		Manage users of a role	+ Add a new role	Functionality Key	Functionality Description
Name	Description			AGENT_MANAGEMENT	Agent Management
admin	Show All	Show All		CLINICAL_QUALITY	Clinical quality
Admin	ADMIN			COMPLIANCE	Compliance Viewer
Admin_role_test	Admin_role_test			COMPLIANCE_MANAGEMENT	Compliance management
ROLE_ADMINISTRATOR	ROLE_ADMINISTRATOR			DASHBOARD_MANAGEMENT	Dashboard management
Series Data Admin	Edit series data			DATA_MANAGEMENT	Data management
Super Administrator	SUPER_ADMINISTRATOR			DEVICE_INTEGRATION_MANAGEMENT5	Device integration manager
				DEVICE_LEVEL	Device level
Username	Real Name				
Show All	Show All				
Ana	Ana Dolcet				
AnaTest	Ana Dolcet				
AnaTestGerman	Ana TestGerman				
Bruno	Bruno Borges				
Daniel2	Daniel				
FirstLogin	FirstLogin FirstLogin				
GuillaumeD	Guillaume Deschamps				

By clicking on the role edit button (the blue button on the right end of the table), a dialog is shown to allow editing the name of the role, its description and the list of functionalities assigned to it.

Modify role entry

Modify the parameters for the current role entry Remove


Role Name
Admin

Role Description
ADMIN


Functionality
<input checked="" type="checkbox"/> Agent management
<input checked="" type="checkbox"/> Clinical quality
<input type="checkbox"/> COACH management
<input type="checkbox"/> COACH read
<input checked="" type="checkbox"/> Compliance
<input checked="" type="checkbox"/> Compliance management
<input checked="" type="checkbox"/> Dashboard management
<input checked="" type="checkbox"/> Data management
<input checked="" type="checkbox"/> Device integration manager
<input checked="" type="checkbox"/> Device level
<input checked="" type="checkbox"/> Device management
<input checked="" type="checkbox"/> Dose management
<input checked="" type="checkbox"/> Export management
<input checked="" type="checkbox"/> Live Dashboard
<input checked="" type="checkbox"/> Live Dashboard Modality
<input checked="" type="checkbox"/> Location level
<input checked="" type="checkbox"/> Logbook
<input checked="" type="checkbox"/> Logbook Administrator
<input checked="" type="checkbox"/> Modality level

Cancel
 Save

The same dialog is shown when clicking on the Add a new role button.




Users can be assigned to the role on the User tab, editing each user. As an alternative, the button *Manage users of a role* of selected role allows assigning several users at a time to a role.



2.2.3. Teams tab

The interface of the Roles tab is coherent with the other tabs. By clicking on a single team in the table, all details of the team are shown. In this read-only mode it is possible to search for a device or a user belonging to the team.



By clicking on the team edit button (the blue button on the right end of the table), a dialog is shown to allow editing the name of the team, its description, the parent team and the list of devices belonging to it.

Update Team [Close]

Modify the parameters for current team entry Remove


Parent team
root_team

Team Name: Hospital 1 Team Description: only hospital 1


Add devices for the current team Add Devices

Site Name	Modality	Device Name	
Demo Hospital 1	RF	Room 07 - Siemens Fluorospot Compact FD	
Demo Hospital 1	CT	SIEMENS SOMATOM Definition AS+	
Demo Hospital 1	CT	SIEMENS SOMATOM Definition Flash 2	
Demo Hospital 1	CT	SIEMENS SOMATOM Definition Flash	
Demo Hospital 1	CT	GE MEDICAL SYSTEMS BrightSpeed	
Demo Hospital 1	CT	SIEMENS SOMATOM Force Data An	
Demo Hospital 1	CT	Room 05 - Siemens Somatom Definition Flash	
Demo Hospital 1	NM	Room 08 - NM Hotlab	
Demo Hospital 1	CBCT	Room 09 - NIM NT5G	
Demo Hospital 1	CT	Room 04 - Toshiba Aquilion 16	

The red button on the right of each row in the table allows removing the devices, while clicking on Add devices allows selecting devices to add to the team.



Users can be assigned to the team on the User tab, editing each user. As an alternative, the button *Manage users of a team* of selected team allows assigning several users at a time to a team.



2.3. User management [deprecated]

In the **Settings** of DOSE, a section called **User Management** can be found. Upon clicking on it, 4 tabs can be accessed:

- User Overview: an overview of all users, possibility to create, edit, or lock a user
- Team Management: used to create/edit/view the teams
- User-Team: to assign and check user teams
- User-Role: to assign and check user roles

2.4. User overview


User Overview contains a table displaying an overview of all registered users. Users can be created or edited by using the “Actions” button in the bottom left corner.

USERNAME	E-MAIL ADDRESS	REAL NAME	LOCKED	LANGUAGE
agfa	agfa@agfa.agfa	Agfa	NO	en
evgenia.boldyreva	evgenia.boldyreva@qaelum.com	Evgénia Boldyreva	NO	en
filippo.miniati	filippo.miniati@qaelum.com	Filippo Miniati	NO	en
qaelum	jurgen.jacobs@qaelum.com	Jurgen Jacobs	NO	en
roleAdmin	roleAdmin@roleAdmin.com	roleAdmin	NO	en
sam.lox	sam.lox@qaelum.com	Sam Lox	NO	en

User overview tab

It is possible to activate or deactivate the “Toggle Filter Bar”, allowing to filter by the Username, E-mail address, Real Name, Locked/Unlocked status and Language. With the


dropdown arrow additional columns can be shown (e.g. Registration date, Reason For Lock, etc.). The toggle bar is deactivated by default.



Toggle filter bar

2.5. Team management

Teams are used to allow access to specific sites and devices. The **Team Management** tab is used to manage the teams and team's filters (*see Team Structure below*).




Team management tab

2.6. Team structure


A team is defined by one or more filters that contain a list of specific sites and devices. A single or multiple teams can be assigned to a user. In order to explain how it is managed, the scheme “Team and filter structure” will be used. In this example, User1 and User2 are part of the Team1. Because Team1 is defined by Filter1 and Filter2, User1 and User2 will only have access to the specific sites and devices contained in Filter1 and Filter2. User4 is part of Team2 and will have access to the sites and devices specified in Filter3. User3 is part of both Team1 and Team2 and will have access to all sites and devices specified in Filter1, Filter2 and Filter3.

It is also possible to give access to a site but not to its devices. In this case, the user will be able to see all the information about the site in DOSE, but will not have access to any devices of this site.

In order to have a certain hierarchy in the user management, every team has a parent team, as shown in the scheme “Team relationship”. In this example, User1 is part of Team1, which is the parent of Team2. In this case, User2 will only have access to the sites and devices defined for Team2, but User1 will have access to all sites and devices from Team1 and Team2, even if Team2 is not directly assigned to him.



Scheme users filter_1




Scheme users_2

2.7. Create / Edit a team

To create or edit a team, click on **Create/Edit Teams**.

Select a **Parent team** that can be an existing team or standalone (**root_team**), enter a name and a description for the new team, and click on the **Create** button. All users from the selected parent team will have access to the new team.



To edit a team:

- select the team to be edited
- modify the parent team, if necessary
- edit the team name and/or the description, if necessary
- save changes by clicking on the **Edit** button

It is also possible to delete a team. To do so, select the desired team and click on the **Delete** button.

!!WARNING!!

If you delete a team, all sub-teams (children) of this team will also be deleted. This means that specific access to the site/devices that was linked to the team will not be possible anymore by users of the deleted team.


Create/Edit team

2.8. Create a filter

Creation of a new filter is performed using the **Create/Edit Teams** window

- Enter a name and description of the new filter
- Select the team that needs to be filtered (Team List)
- Select the sites to be added to the new filter (Site List)
- Devices can be filtered by *Modality, Vendor, Model, and Dose information* in the **Filter** section
- Click on **Submit Filter** to get the list of devices

- Select the devices to be added to the new filter by ticking the checkbox next to the devices
- Click on **Save Filter**




When a team is selected, all the filters associated with this team are shown at the bottom of the window.

For more information, refer to the video *How to create a new team and customize user access* in our online training center.

2.8.1. Edit a filter

To edit a filter:

- Click on the blue **Edit Filters** button at the bottom left
- Select the team that the filter is attached to
- Select the filter
- Once the filter is selected, the list of sites and devices within the filter are displayed
- The name/description, Site List, and Equipment List within the filter can be modified by selecting/unselecting the desired sites and devices
- Click on the **Edit** button to save the changes



Filters can also be deleted by selecting the filter and clicking on the **Delete** button.


For more information, refer to the video *How to modify an existing team* in our online training center.

2.9. User-Team

When a team is created, it is possible to assign this team to users within the **User - Team** tab. To check the current team(s) that users belong to, select the user(s) and click on the blue button **Check User Teams**. To assign a new team to a user:

- Select the user
- Select the team(s)
- Click on the blue **Save** button

It is also possible to remove a user from a team by selecting the user and the team, and clicking on the red **Remove Selected Team(s)** button.



User team tab

2.10. Role Management


Roles are used to grant specific permissions to users. A dedicated customized list of functionalities can be assigned to each role. Each functionality allows access to certain features of DOSE. This means that every user can be easily assigned a suitable function according to her/his tasks in the hospital. The **User – Role** tab is used to manage the roles. The User List and Role List are displayed, along with an overview of all existing roles and associated functionalities below. From this window it is possible to:

- Check user roles by selecting users and clicking on the blue **Check User Role** button
- Assign a role to a user
- Remove a role from a user
- Create or edit a role

2.10.1. Assigning and removing user roles

To assign a role to a user, select the user and the desired role and click on the blue **Save** button.

To remove a role from a user, select the user and the role to be removed and click on the red **Remove Selected Role(s)** button.



For more information, refer to the video *How to create a new user and assign teams and roles* in our online training center.


To create or edit a role, click on the blue **Create/Edit Roles** button.

To create a new role:

- Enter a name and a description
- Check all desired functionalities for this role
- Click on the blue “Create” button

To edit a role:

- Select the role to be edited
- In the left panel, all functionalities that are already set for the selected role will be displayed
- In the right panel, all the functionalities assigned to this role are shown as checked
- The name and/or the description can be edited
- The desired functionalities can be checked to add and/or unchecked to remove
- Click on the “Edit” button to save all changes



Create/Edit roles window

For more information, refer to the video *How to modify an existing role* in our online training center.

2.11. Restrictions

In order to avoid potential mistakes in using the User Management tool, some restrictions are in place. A list of all restrictions can be found below.

General restrictions

- Only the SuperAdministrator, RoleManager and users with the “User Management” functionality can access User Management in Settings.

Team management restrictions

- The SuperAdministrator and RoleManager can see everything, including all existing teams, sites and devices.
- The SuperAdministrator and the RoleManager only can delete a team.
- Users with the "User management" permission can only access teams they belong to.
- Users with the "User management" permission can only see the sites and devices they have been assigned access to, depending on the teams that they are part of. If a user is part of a team that can only see site A, she/he will only be able to create or edit teams within site A.

User – Team

- SuperAdministrator can see all users and all teams.
- RoleManager can see all users and all teams except the "root_team".
- Other users with the "User management" permission can see all users but can only see the teams they are part of.

User – Role

- SuperAdministrator can see all users and all roles.
- RoleManager can see all users and all roles except the roles "Role_administrator" and "Super Administrator".
- SuperAdministrator and RoleManager only can delete a role.
- Other users with the "User management" permission can only create or edit roles that have been assigned to them.
- Users with the "User management" permission can only assign the roles they hold to other users.

2.12. Default roles and their functionalities

2.12.1. Admin

Access to:

- Agent Management
- Clinical Quality
- Compliance Management
- Dashboard Management
- Data Management
- Device Level
- Device Management
- Dose Management
- Export Management
- Location Level
- Logbook Administrator
- Modality Level
- Notification Administrator
- Patient Level
- Physical parameters Management
- Scheduler Management
- Settings Management
- Site Management
- Study management
- User Management

2.12.2. Radiologist

Access to:

- Clinical Quality
- Compliance
- Device Level
- Export management
- Logbook Administrator
- Modality Level
- Notification
- Patient Level

2.12.3. Radiographer

Access to:

- Compliance
- Device Level
- Export Management

- Logbook
- Modality Level
- Notification
- Patient Level

2.12.4. Physicist

Access to:

- Compliance Management
- Device Level
- Device Management
- Dose Management
- Export Management
- Location Level
- Logbook Administrator
- Modality Level
- Notification Administrator
- Patient Level
- Physical parameters management

2.12.5. Quality Manager


Access to:

- Clinical Quality
- Data Management
- Device Level
- Export Management
- Location Level
- Logbook Administrator
- Modality Level
- Notification Administrator
- Patient Level
- Scheduler Management

For more information about these permissions or functionalities, please consult the *Description of Roles* document.


3. Quick overview

This start page gives users an 'easy landing' page that shows a simple KPI analysis of the department. When users log in to the system, 'Quick Overview' will be the first page that will be shown. There is a possibility to change to 'Expert mode' via the drop-down menu in the upper right corner.



Quick overview

The question mark button in the right corner explains all the components of the Quick overview.



Users have the possibility to choose the KPI overview for all modalities or only for one modality by selecting and clicking on "Update"

Choice of modality for KPI overview

Refer to the figure below showing each numbered field from the Quick overview.


1. From this menu the extended version of DOSE can be accessed.
2. Percentage of studies in the 'achievable' range.
3. Number of studies performed.
4. The red bars show the average value of the previous four weeks for each weekday.
5. Number of patients examined.
6. The color of the arrow shows whether the trend is an improvement.
7. Average study duration (from the moment the study was started to the moment the study was closed).
8. The percentage shows the difference between the current value and the average of the four comparison points.
9. The current weekday, up until the current time, is compared against the same weekday, up until the current time, of the previous four weeks.
10. Percentage of studies in the acceptable range.
11. The blue line is the average of the four other bars.
12. Percentage of studies outside the acceptable range.



Quick overview


4. Start portal

After the login to DOSE, the user arrives in the Start Portal: this is a personal page that every user can customize with charts and tabs the user needs the most.




Start portal

The Start portal tab is divided into 2 tabs: the 'Overview', with the favorite graphs, and 'Favourites', where the user can manage the previously starred graphs. In Favourites, the user can remove a component from the Overview page. Preferably, a maximum of 5 components should be added to the Overview page.



Favourites


Please refer to the specific section on how to mark (star) charts and tabs as favorites and how to set them as part of the user's personal Start Portal.




Start Portal - overview

5. Study level

On **Study** level users have the possibility to look up studies based on search values (i.e. Accession Number, Study ID or Study UID).




Search values



ALERTS	STUDY DATE	PATIENT GENDER	PATIENT AGE	MODALITY	DEVICE NAME	STUDY DESCRIPTION	BODYPART EXAMINED	EFFECTIVE DOSE	PARAMETER TYPE
3	2015-10-23 11:09:32	F	29 Y	CT	Room 00 - Philips Ingenuity	RAD ct abd 22	ABDOMEN	7.64 mSv	424.7 mGy.cm

Study level filtered by accession number

The **Parameter Type** shown in the list can be changed by clicking on the **Filter** button in the upper-right corner.



Change the parameter type shown on the list

A study can be created manually by clicking on the **Actions** button in the lower-right corner and selecting **Add study**. The study details can be entered here.

Study details

Accession Number *	<input type="text"/>
Study ID *	<input type="text"/>
Study UID *	<input type="text"/>
Study description	<input type="text"/>
Body Part Examined *	HEAD
Study Date *	11/11/19 03:12:51.844 PM
patientId *	<input type="text"/>
Gender *	M
Age Type *	Y
Age *	0
Weight (kg) *	0
Length (cm) *	0
Collection ID *	<input type="text"/>
Modality *	<input type="text"/>
Parameter type *	<input type="text"/>
Total Dose *	0.0

Manual addition of a study

When weight and height are entered, the body mass index will be calculated automatically.


For interventional procedures and NM examinations, it is possible to create studies manually with specific templates (see [Manual entry of an XA study](#) and [Study Overview section in the NM chapter](#)).

5.1. Study details page

This page opens when a study is selected (by double clicking).

5.1.1. General study information (Overview)

Here all general information about the study and the patient can be found. If weight and height are collected, the BMI is automatically calculated and displayed.



Study Details
Patient Name: Unknown

Overview Series information Contrast details Organ doses History Activity stream

Study date
23/10/2015 11:09

Study description
RAD ct abd 22

Category
Abdomen

DLP Limit(s)
> ACH (upper): 300.0 mGy.cm
> ACC (upper): 700.0 mGy.cm

Max CTDIvol [mGy]
8.1

Body part examined
ABDOMEN

Manufacturer
Philips

Study ID
123

Operator
Operator 1

Patient information
F | 29 Y |

Accession Number

Total mAs
2.76E5

Modality
CT

Modelname
Ingenuity Core

Study UID

Performing physician
Performing physician 1

DLP [mGy.cm]
424.7

Effective Dose ICRP103 [mSv]
7.64

Total scanlength (mm)
477

Event #
3

Software version
4.1

Device Name
Room 00 - Philips Ingenuity Core

Requesting physician
Requesting physician 3

Hide notifications

STATIC LIMIT ▲ DLP

Study details overview

If a certain study/series is linked to a study group with static limits, it is possible to proactively see whether this study is compliant to the limits.

ACH = Achievable level

ACC= Acceptable level

The compliance bar color shown below the thumbnail follows this logic:



- Default – gray (not linked to a study group)



- Below ACH level - green
- Between ACH and ACC level - orange
- Above ACC level - red


This compliance color bar will always show the color of the highest compliance notification. For example, if a CT study has an orange DLP notification and a red CTDI notification, the resulting compliance bar will be red.

An exception to this behavior in DOSE is the Compliance Monitoring table, in which case the compliance bar will be the color corresponding to the selected parameter (DLP, CTDI, etc.) for optimization purposes.

If the series contained in a study are linked to a study group limit, their compliance color will be shown in Series Details and the compliance bar at study level will show the color of the series with the highest color.

The compliance bar does not take into account other notifications that are not directly generated by the study group configuration (i.e. Notification Center rules, dynamic alerts and manually created notifications).

To see the study groups that a study/series is linked to, the user can click on *Show study groups* below the thumbnail.



The screenshot shows a study details page with the following information:

- Study date:** [M] 62 Y
- Study description:** RAD ct lwz, [12-7-2-88]
- Protocol name:** 02_LWZ
- Hide study details:** [button]
- Hide Study groups:** [button]
- Study group table:**

Study group	Parameter and level	Lower ACC	Lower ACH	Upper ACH	Upper ACC
Lumbar spine	CT_EFFECTIVE_DOSE_103	0.0	0.0	0.0	8.9
	mGy cm (STUDY)				
	DLP mGy.cm (STUDY)	0.0	0.0	0.0	592.0

Study groups linked to a study



Radiation triangles


The radiation triangles that are shown next to/below the thumbnail show a relative representation of the effective dose (according to American College of Radiology (ACR) Appropriateness Criteria). Please note that the effective dose is based on the BodyPart, and thus if it is missing, the triangles are empty. Below you can find the explanation of the Relative Radiation Level.

Table 1. Relative radiation level designations along with common example examinations for each classification

Relative Radiation Level*	Adult Effective Dose Estimate Range	Pediatric Effective Dose Estimate Range	Example Examinations
○	0 mSv	0 mSv	Ultrasound; MRI
☢	<0.1 mSv	<0.03 mSv	Chest radiographs; Hand radiographs
☢☢	0.1-1 mSv	0.03-0.3 mSv	Pelvis radiographs; Mammography
☢☢☢	1-10 mSv	0.3-3 mSv	Abdomen CT, Nuclear medicine bone scan
☢☢☢☢	10-30 mSv	3-10 mSv	Abdomen CT without and with contrast; Whole body PET
☢☢☢☢☢	30-100 mSv	10-30 mSv	CTA chest abdomen and pelvis with contrast; Transjugular intrahepatic portosystemic shunt placement

*The RRL assignments for some of the examinations cannot be made, because the actual patient doses in these procedures vary as a function of factors (eg, the region of the body exposed to ionizing radiation, the imaging guidance that is used, etc.). The RRLs for these examinations are designated as "Varies."


Please note that dose is shown in scientific notation when its value exceeds 1000 or is inferior to 0.01. The point is used to separate units from decimals.



Dose notation


5.1.1.1. NOTIFICATIONS

Alerts are created when study/series dose parameters are outside acceptable or achievable limits. All study alerts are shown in the **Notifications** tab. Notifications can be for the dynamic or static limits (according to the study group and limit configuration) or for a Notification Center rule.



Visual limits

Hovering over these alerts will display alert details and double-clicking will open the **Notifications** tab.



Notification tooltip

Study Details

Patient Name: Unknown

Cancel Actions

Overview Series information Organ doses History Activity stream Notifications

Study Date	Origin	Limit type	Parameter key	Type	Study	Message	Value	Difference	Range	
07/07/2017 08:19:53	Series	STATIC	CTDI	RED	Test Ana SERIES - 83 Y M	Outside acceptable range Value: 12.84 [0,10]	12.84	28 %	[0.0; 10.0]	Handle Details
07/07/2017 08:19:53	Series	STATIC	DLP	RED	Test Ana SERIES - 83 Y M	Outside acceptable range Value: 610.88 [0,500]	610.88	22 %	[0.0; 500.0]	Handle Details
07/07/2017 08:19:53	Study	STATIC	EFFECTIVEDOSEICRP103	GREEN	83 Y M	achievable Effective dose ICRP 103	0	-0 %	[0.0; 0.0]	Archived Details
07/07/2017 08:19:53	Series	STATIC	CTDI	RED	COMPLIANCE TEST - 83 Y M	Outside acceptable range Value: 12.84 [0,12]	12.84	7 %	[0.0; 12.0]	Archived Details

Notifications tab

If archived notifications are present (i.e. notifications that were saved after a recalculation), they can be visualized by selecting the **Show archived notifications** checkbox.



Archived notifications

For more information about notifications, please refer to the **Notifications** section in this User Manual.

5.1.1.2. CLINICAL QUALITY EVALUATION

Clinical quality evaluations (CIQs) make it possible to perform a subjective scoring of the clinical image quality. This is possible from every **Actions** menu in the Study Details and from the **Live dashboard**. On the CIQ evaluation page the user (with the appropriate role assigned in User Management) can score the image contrast, image sharpness and patient positioning. Additional comments can also be entered.

This feature is mostly used as PACS integration (Webservice). For more information see the **Clinical Quality Evaluation** section of **Webservice integration**.



Clinical quality evaluation scores

If no evaluations are found, the overview is collapsed by default. If data is present, the evaluations are visible.


5.1.1.3. STUDY COMPARISONS

The current study can be compared to all studies from the same device or modality, or to studies with the same study description, examined body part, study composition, protocol name. A histogram is available for DLP, CTDI, scanlength, number of series, series analysis and study composition. Every histogram can be expanded to full window, exported, saved and printed. The histogram data can also be exported.



Study comparison

5.1.2. Series information



DATE TIME	ACQ. PROTOCOL	CTDI (mGy)	DLP (mGy*cm)	PEAK VOLTAGE (kVp)	TUBE CURRENT / REF. TUBE CURRENT (mAs)	PHANTOM	SCAN LENGTH (mm)	ACQ. TYPE	PATIENT POSITION
2013-08-01 11:21:27	Topogram AP	0.0	0.0	120	35		768	Constant Angle Acquisition	
2013-08-01 11:21:27	Aorta z.k.	12.19	736.41	120	446	IEC Body Dosimetry Phantom	604	Spiral Acquisition	
2013-08-01 11:21:27	PreMonitoring	2.74	2.83	120	54	IEC Body Dosimetry Phantom	9	Sequenced Acquisition	
2013-08-01 11:21:27	Monitoring	16.43	15.77	120	54	IEC Body Dosimetry Phantom	9	Sequenced Acquisition	
2013-08-01 11:21:27	Aorta airt	13.72	820.22	120	572	IEC Body Dosimetry Phantom	598	Spiral Acquisition	
2013-08-01 11:21:27	Aorta ven	18.64	1107.14	120	512	IEC Body Dosimetry Phantom	594	Spiral Acquisition	
2013-08-01 11:21:27	Thorax	13.72	432.06	120	572	IEC Body Dosimetry Phantom	315	Spiral Acquisition	


Series information tab

This tab contains all the information about the series/instances that comprise the selected study.

By selecting the spiral/axial CT acquisition protocols, Advanced CT Analysis (if applicable/available) can be accessed.

5.1.2.1. WEIGHTED CTDI_{vol}

Weighted CTDI_{vol} is a weighting factor of the acquisitions (excluding localizers) taking into account the scan length of the acquisition. This length is calculated using the ratio of DLP/CTDI_{vol} instead of the DICOM scan length, as it is more consistent between different scanner vendors. This ratio is shown as a column in the series information.



Bodypart	Effective Dose ICRP 103	DLP / CTDI _{vol} [mm]
HEAD	0	N/A
HEAD	0.1	175.16

Ratio DLP/CTDI_{vol} for calculation of Weighted CTDI_{vol}

$$WeightedCTDI_{vol} = \frac{TotalDLP}{TotalScanlength}$$

where

$$TotalDLP = \sum_i^n DLP_i$$


$$Totalscanlength = \sum_i^n \frac{DLP_i}{CTDI_{vol_i}}$$

Where *i* is each series (excluding topograms) and *n* the total number of series for that study.

Weighted CTDI_{vol} is included as a parameter for the DRL configurations in **Settings/Study groups** and as a parameter for the automatic notifications. For more information see [Compliance Configuration](#) and [Notification Center](#).


5.1.2.2. SIZE-SPECIFIC DOSE ESTIMATES (SSDE)

Size specific dose estimate (SSDE) analysis for CT studies can be found in the study details by navigating to the *Series information* tab and selecting the desired acquisition protocol.



Series information tab with subcategories

The SSDE analysis also takes patient positioning into account. The real dimensions and positioning of the patient are measured and compared to the dimensions of a reference phantom to obtain a dose that depends on the morphology and positioning of the patient. In the image the white lines are the geometric dimensions of the patient that are used to calculate the effective diameter. The red line represents the mA modulation and the blue line represents the pixel density. If they are following each other, the tube current modulation is working well.



AP and LAT topogram


DOSE uses scout images for segmentation and slice information to define the actual scan range. The scan range is shown on the scout. The center of the scan range is defined as the reference (marked by the white line). The software segments the patient on the scout images and measures the segmented diameter of the reference slice. This

procedure is done both for AP and LAT (if available) topograms. The effective diameter is calculated either using both AP and LAT or one of them, if both are not available.

The calculated effective diameter is matched to the corresponding conversion factor according to AAPM Report No 204. The SSDE is calculated using the CTDI_{vol} of the examination corrected by the specific conversion factor.

The SSDE_{WED} is then calculated taking into account the patient attenuation (Water Equivalent Diameter, AAPM Rpt 220).

The calculation method of effective diameter and water equivalent diameter can be found in the information table in the SSDE tab (*SSDE calculation method*). The *calculation mode* indicates whether SSDE values are calculated using axial images, truncated axial images, or localizer (AP, LAT or both). The truncation percentage of the received axial images is also calculated, and when it exceeds 20%, the user is notified of the possibility of inaccurate calculations.




SSDE and SSDE_{WED}

5.1.2.3. PATIENT POSITIONING

In the graphical representation of the patient positioning tab, the blue circle shows the expected centering of the CTDI phantom and the orange ellipse represents the patient, based on the calculated effective diameter. The large white circle represents the gantry, the central point of the two axes is the isocenter, x-axis is the left-right direction, y-axis is the top-bottom direction.


The ellipse representing the patient dimensions (AP and lateral) is generated using the SSDE calculation according to the effective diameter. The horizontal, vertical and radial offsets are calculated and displayed in the blue box.

The source to isocentre distance and table height are extracted from the header, and the lateral positioning is calculated using image data. The difference between the gantry isocentre and the patient's isocentre is used to calculate the vertical, horizontal and radial patient offsets.



Patient positioning

As the primary option for the calculation is the axial images, it does not always correspond to a calculation using the table height. The gap between the estimated patient position and the tabletop is covered by an assumed pad.




Assumed pad

5.1.2.4. TUBE CURRENT MODULATION

Tube current modulation should adapt the tube current to the attenuation of the body region using the water equivalent diameter as reference, as follows:

- o increase tube current for higher attenuating area
- o decrease tube current for lower attenuating area

The overall goal is to reduce dose while maintaining image quality. All plotted data (tube current and three effective diameters) can be clicked to see the exact slice location on the localizer, represented by the red line. The data can also be hidden by clicking on the plotted value name(s) below the chart.




mA modulation

The chart can be exported using the options of the gear button in the upper right corner of the chart. The exported file includes all chart information (slice location, tube current, effective diameter from localizer, effective diameter from axial images, WED from axial images).

5.1.2.5. SERIES OVERLAP

When a study is composed of more than one series, overlap between scanned areas is calculated, and the overlap length is displayed in mm.



Series overlap

5.1.3. Peak Skin Dose (PSD)

This tool is designed to help users estimate the patient's skin dose after an interventional procedure. As the skin dose is directly linked to deterministic effects / tissue reactions, this tool allows for faster and improved follow-up strategy. The peak skin dose analysis can be found under the **Peak Skin Dose** tab, located within the **Study details** interventional procedures. Peak skin dose calculation is based on the parameters extracted from the Radiation Dose Structured Report (RDSR) of the particular interventional procedure. Correctness of calculated values depends on the correctness of the collected parameters. Additionally, several assumptions are taken into consideration in the calculation that could lead to a relatively different result compared to other calculation methods (e.g. table and pad attenuation can differ from system to system). For other limitations and assumptions please refer to the "DOSE Technical Information" document that you can find in the "About" button.

Peak Skin Dose is calculated for head and trunk examinations of adult and pediatric patients. Extremities are currently not supported.

The RDSR provides geometry (tube angles, table movement) and dose data (dose area product, air kerma at reference point) for every exposure event. The collimator shape is not available in the RDSR, and thus is predefined for the calculations. If the collimated irradiated area for each event is available, it is taken into consideration. If the area changes during the irradiation, then only the value mentioned in the RDSR for this event is used. The backscatter radiation factors are based on the thickness, field size and beam quality as described in literature (Benmakhlof et al 2011 and 2013) for each irradiation event. A standard attenuation from the table and table pad is considered. The factor corresponds to a specific type of table and table pad (total thickness of 8 cm).

Besides the adult phantom, five pediatric phantoms are included in the calculations (15 year-old, 10 year-old, 5 year-old, 1 year-old and newborn). The dimensions of the standard phantoms are defined as the corresponding head and body sections of the MIRD phantoms.


The phantom is initially selected based on the age of the patient. If the dimensions of the patient are filled in (weight and height), the phantom is scaled to match the actual patient size. The calculation is performed on the corresponding head or body phantom, based on the bodypart of the first irradiation event. All phantoms are considered headfirst and supine due to the lack of exact position and orientation on the table. If the information for patient age, patient size or the target region ('bodypart') as indicated in the data source are not correct, then this information can be corrected directly within DOSE, and peak skin dose map recalculated. Patient age and size can be edited by selecting *Edit study* in **Actions** (top right corner of the Study overview), while the body part can be edited within the *Series information* tab, by right clicking and selecting *Set series bodypart* (for all series of this study). This may happen in protocols not assigned to the correct body part by default.



Changing patient size and/or age

	Series Description	Radiation mode	kVp	Exposure [mAs]	XRy Tube
12	Exam Protocol SR	Fluoroscopy	81	80.56	
22	Exam Protocol SR	Fluoroscopy	78	40.21	

[Edit Serie Data](#)
Set series bodypart



Changing the body part for the study


Spatial distribution of the skin dose is shown as a map of 1x1 cm² pixels on the phantom surface. Each pixel value corresponds to the calculated dose. The map is recalculated with every dose event to add the associated dose value. This allows to follow the evolution of the skin dose distribution during the intervention. Peak skin dose is determined as the maximum dose value on the map.

As there no information indicating the exact positioning of the patient on the table is available, the displayed dose distribution map is shown centrally located on the patient. The exact location can differ, e.g. it can be closer to the head or closer to the feet than displayed. The physician should be aware of the center of irradiation, depending on the type of procedure.

- **Peak Skin Dose**


Peak skin dose map. The surface of the reference patient is represented in pixels of 1 cm². The x-axis is the transverse perimeter of the patient, while the y-axis corresponds to the length of the trunk (head-feet direction). The map represents the dose distribution on the surface of the patient, as if the surface is unfolded at the right part of the patient. To distinguish between back and front side of the patient, the front side is represented by grey pixels. The user must keep in mind that for the PSD calculations, the patient is considered supine and head first. For other types of positioning, the map must be interpreted accordingly. An additional note for the user is that the exact location of the dose distribution on the skin can differ from what is displayed on the map, as there is no information concerning the exact positioning of the patient on the table. The physician should interpret that, depending on the type of procedure. The pixels that have been

irradiated are represented by a color corresponding to the dose delivered there. The color-dose correspondence is shown in the legend. A sentence above the skin dose map shows the phantom (head or body), the patient age and size considered for the calculations. If no patient size information is available, then the standard dimensions of the phantom are considered and the phantom is selected only by age.




Peak skin dose map

At the side of the colored legend, the user can see interesting information about the skin dose and the area irradiated. In particular besides the peak skin dose value and the total irradiated area, the area that received 10%, 50% and 90% of the peak skin dose are calculated based on the colored pixels in the map. The user can also follow the evolution of dose distribution throughout the intervention (animated gif passing through events). In order to understand the 2D map easier, a 3D representation is introduced. The user can see it next to the animated gif and can also enlarge it by clicking it. The enlarged 3D phantom can be manually rotated by holding and moving the left button of the mouse.




Information about peak skin dose and irradiated area, animated gif and 3D phantom

At the right part, there is a menu where the user can select to see the impact on the map from each of the planes in case of a biplane system, or a combination of them.



Impact on skin dose map by different planes can be visualized.

Near this menu, there is an exclamation mark, where the user can click on the **PSD calculation method** to get extra information about the methodology used, as well as about the dimensions of the standard phantoms. By clicking **Potential radiation-induced skin reactions & advice to patient**, the user is guided to a table that provides information about radiation-induced skin injuries (with the corresponding references). One of the table rows is highlighted (colors according to the severity, with green less severe) to show in which category the particular PSD belongs.



PSD calculation method and potential skin reactions.

Band	Single-Site Acute Skin Dose Range (Gy)	NCI Skin Reaction Grade	Approximate Time of Onset of Effects				Advice to Patient
			Prompt (< 2 wk)	Early (2-8 wk)	Midterm (6-52 wk)	Longterm (>40 wk)	
A1	0-2	NA	No observable effects expected	No observable effects expected	No observable effects expected	No observable effects expected	No need to inform patient, because there should be no visible effects; if patient reports skin changes, then treat in response to the signs and symptoms
A2	2-5	1	Transient erythema	Epilation	Recovery from hair loss	No observable effects expected	Advise patient that erythema may be observed but should fade with time; Advise patient to call you if skin changes cause physical discomfort
B	5-10	1-2	Transient erythema	Erythema, epilation	Recovery, at higher doses, prolonged erythema, permanent partial epilation	Recovery, at higher doses, dermal atrophy or induration	Advise patient to perform self-examination or ask a partner to examine for skin effects from about 2 to 10 weeks after the procedure; tell patient where skin effects would most likely occur; if skin erythema and itching occur patient should call radiologist's office; skin reactions are often treated conservatively, might advise patient to be examined by dermatologist or other treating physician and to inform treating physician that injury may be due to radiation; radiologist should also provide that physician with medical details of where the radiation related skin effects are likely to occur
C	10-15	2-3	Transient erythema	Erythema, epilation, possible dry or moist desquamation, recovery from desquamation	Prolonged erythema, permanent epilation	Telangiectasia, dermal atrophy or induration; skin likely to be weak	Medical follow-up is appropriate; advice is same as that for band B but also advise dermatologist or other treating physician that skin effects may be prolonged due to radiation dose and that prophylactic treatment for infection and monitoring of wound progression may be required; pain could become a concern if doses were in the higher range of this band
D	15-...	3-4	Transient erythema; after very high doses, edema and acute ulceration; long-term surgical intervention likely to be required	Erythema, epilation; moist desquamation	Dermal atrophy; secondary ulceration due to failure of moist desquamation to heal; surgical intervention likely to be required; at higher doses, dermal necrosis, surgical intervention likely to be required	Telangiectasia; dermal atrophy or induration; possible late skin breakdown/wound might be persistent and progress into a deeper lesion; surgical intervention likely to be required	Medical follow-up is essential; nature and frequency of which, depending on estimated radiation dose, advice is same as that for band C, but advise treating physician that the wound could progress to ulceration or necrosis

***Notes:**

- Applicable to normal range of patient radiosensitivities in absence of mitigating or aggravating physical or clinical factors. Data do not apply to the skin of the scalp. Dose and time bands are not rigid boundaries. Signs and symptoms are expected to appear earlier as skin dose increases. Prompt (< 2 weeks), early, 2-8 weeks; midterm, 6-52 weeks; long term, > 40 weeks.
- Advice is applicable to normal range of patient radiosensitivities in the absence of mitigating or aggravating physical or clinical factors.
- Skin dose refers to actual skin dose (including backscatter). This quantity is not the reference point as kerma described by Food and Drug Administration (21 CFR § 1020.32 [2007]) or International Electrotechnical Commission (57). Skin dosimetry is unlikely to be more accurate than 5-50%. NA – not applicable.
- NCI – National Cancer Institute.
- Refers to radiation-induced telangiectasia. Telangiectasia associated with area of initial moist desquamation or healing of ulceration may be present earlier.


Citations:

- "Guidelines for Patient Radiation Dose Management", Stecker et al, J Vasc Interv Radiol 2009
- "Therapeutic Radiation Interventional Procedures: A Review of Radiation Effects on Patients' Skin and Hair", Beller et al, Radiology 2010


Table for potential skin reactions for each peak skin dose range.

- Angular dependency analysis

The graph indicates the primary and secondary angles that were used in the procedure. The colors show the frequency that each angle combination is used, with dark red being the most frequent. The signs of the angles are defined as in DICOM standards. The grey-colored pixels indicate the primary angles for which the tube is above the table. The x-axis corresponds to tube primary angle (left-right movement) and y-axis corresponds to tube secondary angle (craniocaudal movement).



Usage of primary and secondary angle during the procedure




Select a plane (or All planes) then select a parameter and Analyse

For the angular analysis, the user can select to see the angles used for a particular plane or for all planes. Next to that, there is a menu where the user can select the parameter for which the graph will provide the angular dependency and then analyse. For example, by selecting DAP, the graph will color the angles according to the distribution of DAP as coming from each angle combination (primary and secondary).


- **DAP trend analysis**

The graph differentiates the origin of the delivered DAP. The user can see the percentage of DAP that is delivered by fluoroscopy, stationary acquisition, rotational acquisition, etc. throughout the procedure. This is a relative chart, which means that the total DAP given by each moment in the procedure is considered 100% and the percentage of each irradiation type can be seen in the y-axis for the corresponding color.



Contribution to DAP by different irradiation modes

By passing the mouse above each line, the user sees the exact value of DAP delivered for each of the irradiation type until that moment. For example, the point shown in the figure below, corresponds to stationary acquisition DAP of 1382 dGy*cm² which is the 40% of the total DAP until that moment.



DAP source trend analysis

- Charts for dose analysis and angle / position

In all four charts, the x-axis is the real-time of the intervention. By passing the mouse above each line, the user can see the date and time of the procedure and the exact value of the represented parameter at this moment.




Changes in dosimetric and geometric parameters throughout the procedure

- Cumulative dose analysis

The graph follows the increase of maximum skin dose (peak skin dose) and cumulative dose at Reference Point (RP) during the intervention.

If the operators change the irradiation parameters (tube angle, table position, etc), this will be depicted in the graph, as the maximum skin dose will stop increasing, even though irradiation continues.




Cumulative dose analysis

- Max skin dose / Cumulative dose at RP analysis

This graph shows the ratio of Maximum skin dose and Cumulative dose at RP throughout the procedure and it is, as expected, a unitless number. From this graph, the user can see the changes in the technique used.

For example, when the ratio is shown to decrease, it means that even though the cumulative dose at RP increases, the maximum skin dose has stopped increasing or increases with lower rate. This can be due to changes in the tube angles, irradiation field




or table positioning.

Analysis of *max skin dose to cumulative dose at RP*

In order to track the changes in the technique, the user is also able to see the tube angles and table position analysis.

- Tube angle analysis

This graph enables the user to follow the changes of primary and secondary angles throughout the procedure. By passing the mouse above each line, the exact angle value at this moment is shown.



Tube angle analysis

- Table position analysis

The user is able to see the table position in all axes (longitudinal, lateral and height) throughout the procedure.





Table position analysis

For a general overview, refer to the video *How to use Peak Skin Dose* in our online training center.

5.1.4. Organ doses

This tab allows users to check the dose received by every single organ, replacing the effective dose tool. From an optimization/justification/follow-up point of view, this empowers users and can help raise their awareness.



Organ dose comparison

The user can compare the organ dose of the current examination with the 25th/75th percentile and a median organ dose. The comparison can be made with all the studies from the same device/modality and with the same BodyPart/Study description.

The calculation is done for each series and then added to show the total organ doses of the study. If one series has a body part not included in the following list, this acquisition will not be considered in the calculated organ doses.

ABDOMEN
ABDOMENPELVIS
CHEST
CHESTABDOMEN
CHESTABDPELVIS
HEAD
HEADNECK
MAXILLA
NECK
PELVIS
WHOLEBODY


List of suitable body parts for organ dose calculation

! Please note that the organ dose values are indicative, they consider standard settings and standard scanlength for all irradiation events based on the body part and they should be assessed by an expert.

5.1.5. History

History tab shows the dosimetry history of patient. Here the cumulative dose trend for one single patient is shown. The bars represent the annual cumulative effective dose, while the yellow dots represent the number of studies. By clicking on a year, the monthly contribution opens. In order to go back, click the gear and reset the chart. There is a possibility to export the patient dose history in XLS, PDF or CSV (click on the gear in the right corner).


Note: There is a possible security check if you belong to another Team (see [user management](#)), the user is not able to click further on the study details that is performed in another site.



Historical dose overview

5.1.6. Activity stream

In the **Activity stream** users can manually enter the **cause/action/justification** to explain any outliers and share best practices.



Study Details

Overview Series information Contrast details Organ doses History Activity stream *Milestone* Linked documents

Activities

Cause

- 20/11/2015 15:27 & ingria commented on this study
Ectrevasation
- 26/11/2015 15:30 & tom commented on this study
Allergic reaction (headache)

Action

No activity was found for this study. Be the first to comment on this study

[0 comments](#) [1 comment](#)


Justification

No activity was found for this study. Be the first to comment on this study

[0 comments](#) [1 comment](#)

Activity stream

The activity stream for studies within the same device can be evaluated using the **Activity Overview** chart on **Device** level.



Dashboard editor

On **Notifications** level, the *cause/action/justification* are displayed for studies corresponding to the search criteria.

When a comment is present for a study, it will be visible in several places inside the application within the following icon:




Study comment icon

Hovering over the icon will display comment text.

Cause	Action	Justification
Update protocol		
Obese patient		
Wrong protocol		

Tooltip showing the study comment

Comments can also be exported on **Device** level in *Actions/Export data*.




Export study comments from Device Level

For more information please refer to the section [Export](#).

5.1.7. Actions

Various options are available within the **Actions** button in Study Details, situated at the top right corner:

- Edit study
- Open in clinical viewer (*if configured*)
- Deidentify/upload to dose register (*only superadmin*)
- Create and share a link to the study
- Enabling (to ignore studies)
- Export study and / or series information as a csv, xls or pdf file



Action button

- Load the PACS integration service (when available)
- Add clinical quality evaluation
- Send HL7 dose message (*under licence*)
- Create notification
- Recalculate notifications


The user will see only aspects that he/she is allowed see (depending on [user management](#)).

6. Patient level

On Patient level users have the possibility to see the patient radiation dose passport and the study overview / patient history. It is possible to search for patients by Patient ID and/or Patient Name. It is also possible to search by patient name, however, due to the possible breach of privacy, the activation of the patient name search is not done automatically and must be requested.


6.1. Patient radiation dose passport

The cumulated effective dose is shown, together with a timeline of all the patient's studies across different modalities, including the interval time between studies. The bars for each year/month are clickable and exportable (via gear button).



The screenshot shows a search interface for patient radiation dose passports. At the top, there are search fields for 'Search Patient Id' and 'Search Patient Name'. Below these, a search bar contains the number '321' and a search button. To the right of the search bar are links for 'Load possibilities' and a checkbox for 'Restrict possibilities to date range'. Below the search bar, the title 'Patient radiation dose passport' is displayed. The main content area features four key statistics: Gender (F), total study count (7), Cumulated D_{eff} (79.49 mSv), and average study interval (123.17 days).


Patient Radiation Dose Passport - general information



Patient Radiation Dose Passport – more details on study level

6.2. Study overview


The Study overview shows the studies that apply to the chosen filter value. Each row represents a study: study details are shown when the user double clicks on the row.



Alerts	Study Date	Patient Gender	Patient Age	Modality	Device Name	Study Description	BodyPart Examined	Effective Dose	Parameter Type
▲▲▲▲▲	2013-03-20 08:54:38	F	68 Y	CT	Room 01 - SIEMENS SCON	RAD ct th abd[11-2-2-147]	CHEST_ABDOMEN	42.69 mSv	2.24E3 mGy.cm
▲▲▲▲▲	2015-03-10 15:27:41	F	68 Y	CT	Room 01 - SIEMENS SCON	RAD ct schedel[9-7-2-129]	HEAD	3.22 mSv	1.11E3 mGy.cm

Study overview


The user can activate and deactivate the visible parameters in the study overview as in the study overview (series analysis tab) by clicking on the little arrow on the right side.



Alerts	Study Date	Patient Gender	Patient Age	Modality	Device Name	Study Description	BodyPart Examined	Effective Dose	Parameter Type
▲▲▲▲▲	2017-2-173]						HEAD	626.14	1.82 mSv
▲▲▲▲▲	2017-2-173]						HEAD	626.14	1.82 mSv
▲▲▲▲▲	2017-2-173]						HEAD	407.89	1.18 mSv

Choose your filters

There is possibility to save the choice of visible columns for the device of the same modality, as a user setting via an “Action” button in the right corner below.



Save column layout

6.2.1. Color coding of studies

Color coding of the text is used to mark studies with different conditions.

<p>Study description RAD ct thorax 00,[10-11-2-146]</p>	<p>A study is colored white when all parameters are filled in, effective dose is calculated and the body part is available.</p>
<p>Study description RAD ct thorax 00,[10-11-2-146]</p>	<p>A study is colored orange if study effective dose is 0, no dose value is present and/or if the body part examined is missing/not recognized.</p>
<p>Study description RAD ct thorax 00,[10-11-2-146]</p>	<p>A study is colored red if the study is ignored (ignored studies are visible when <i>Show ignored studies</i> option is enabled).</p>


7. Notifications

The purpose of notifications within DOSE is to allow users to manage the vast amounts of data collected by the software, facilitating alert management and handling. Notifications can be generated automatically based on the **study group**, **Notification center** and **dynamic limit** configurations. Users can also create manual notifications, which go with their daily practices. It is possible to trigger a manual notification from Live Dashboard and from Study Details. This notification can be found and handled on Site Level, Device Level or on the Notifications Tab.

7.1. Notifications from Study Details

To create a notification from a specific study, the user must search a study, open the Study Details and then click on **Actions → Create Notification**.

NOTE: Automatic notifications can also be recalculated from here, by clicking on "Recalculate notifications". See [Compliance Configuration](#) for more details.



The screenshot shows the 'Study Details' page. On the left, there's a patient profile with a scan image and a series of warning icons. The main area contains a table with the following data:


Study date	Patient information	DLP [mGy.cm]
09/11/2015 11:12	M 62 Y	1.42E3
Study description	Accession Number	Effective Dose ICRP
CT schedel		4.12
Protocol name	Body part examined	Event #
N/A	HEAD	3

Below the table are links: 'Show study details', 'Show clinical quality evaluations', and 'Show study comparisons'. On the right, a sidebar menu lists various actions, with 'Create notification' highlighted in a red box.

Study Details


A window will open to enter the details of the notification, the user should choose a notification key and write down the Notification Message, for example, "update protocol".

Create notification ✕

 Add the parameters for a new notification

Identifier key <input type="text" value="Study notification"/>	Study UID <input type="text"/>
Device <input type="text" value="KLCT-AZ Groeninge"/>	Modality <input type="text" value="CT"/>
Notification message <input type="text" value="Update protocol"/>	Study overview <input type="text" value="CT schedel - 62 Y M"/>
Origin <input type="text" value="Study based notifications"/>	Supervision <input type="text" value="Supervision advised"/>
Notification status <input type="text" value="Unresolved"/>	Severity status <input type="text" value="ORANGE"/>
Creation date <input type="text" value="3/11/19 04:17:55.982 PM"/>	


Cancel

 Load study
Add notification

Create Notification

Then, after Adding the notification and refreshing, an alert will appear, and the details can be consulted in the Notification tab inside Study Details.

Overview
Series information
History
Activity stream
Notifications



⚠⚠⚠⚠⚠

Study date
09/11/2015 11:12

Study description
CT schedel

Protocol name
N/A

Show study details

Hide notifications

USER CREATED ⚠ Study notification

Show clinical quality evaluations | Show study comparisons



Orange sign for the created notification

Study Date	Origin	Limit type	Type	Study	Message	Value	Difference	Range	
11/10/2013 14:06:45	Study	STATIC	ORANGE	RAD ct comb 02.[18-2-	Outside achievable	658.41	20 %	[0.0 ; 550.0]	Handle Details

Notifications tab in Study Details

The notification can be handled directly from this Notification tab or, if the user has the proper rights (Notification Administrator) and it is a manually created notification, deleted.

Handle notification ✕

 Handle the selected notification 

Creation date : 3/19/19 10:41 AM

Study : RAD ct abd 22 - 29 Y F

Parameter : COMMENT


Description : NM: Patient motion

Handle notification with delete icon

For more information, refer to the video **How to create notifications from a study** in our online training center.

7.2. Notifications Module

This is where users can manage notifications at institution level. It is an extra tab next to the Location Level.




Notifications Tab

The user can select the devices they want to see, the parameter, the type of notification and its status. The numbers above change according to the filter they select.

The procedure would be as follows: select Study based notifications → select a site (If you want you can select them all, starting from the first one and clicking the mouse until the last one or selecting one and pressing CTRL+A) → select modalities → select devices → select a key → select a status → select severity type → click on REFRESH.

All the selected notifications will appear below. You can also see statistics of these notifications in the upper bar.



Notifications tab

The user doesn't have to do all the steps of the selection, but the more they select, the better the findings will be.

After the user has reviewed a notification by opening the study clicking on "Details", they can resolve them by clicking on "Handle", where they can add comments, declare if the physician has been contacted, update its status (resolved/under investigation, etc) and establish needed actions. The notification can also be deleted from here by the Notification Administrator.

Handle notification
✕

Handle the selected notification

Creation date : 7/14/19 1:04 AM
 Study : ALG.ANGIO - 68 Y M
 Parameter : OTHER
 Description : Sentinel event

Canned comments

Repeat examination

▼

Use comment

+

Physician contact

▼

Status

Unresolved ▼

Action

▼

✕ Cancel

↻ Load study

↻ Update

Handle notification

If the study that generated this notification had a comment, it will appear as a comment icon, which the user can hover over to see its content.


Handle notification
✕

Handle the selected notification

Creation date : 07/04/2017 19:05
 Study : N/A - 2 M O
 Parameter : OTHER
 Description :

1

Cause	Action	Justification
Update protocol		




NOTE: If a comment is entered, the status must be changed from "Unresolved" or "Unknown" to another status, e.g. "Under investigation", "Pending", etc. A pop-up message will remind the user to change the notification status.

If a comment is added to this notification, the following icon will be displayed for the study:



Notification comment icon

The user can filter by "Unresolved" notifications to see studies with notifications that have not yet been handled. In the same way, the user can filter the notifications that contain (or not) a comment in the study details, by using the "study comment" filter. These comments include the ones coming from the Live Dashboard, the Pop-up Agent and the dosimetry webservice.




Study comment filter

These comments will appear in different columns in the notifications list.


Study date	Identifier	Study	Description	Comment	Study comment: Cause	Study comment: Action	Study comment: Justification
07/04/2017 19:01	test comment	N/A - 2 M O	Obese patient		Update protocol		

Study comment columns in the notifications overview

The columns may not be visible by default: the user can select the desired columns by clicking on the small arrow at the top right of the table.




The columns can be made visible by default by clicking on **Actions** (top right) and then **Save column layout**.



The content of these columns can be exported via **Actions** → **Export data**.


Notifications can be searched for and handled on site level by accessing **Notifications** when clicking on the site name on **Device** level. In the same way, on Device level the user will have only the notification from the selected device.



For more information, refer to the video **How to use Notifications Level** in our online training center.


7.3. Notification Center

Any user whose role includes the “Notification administrator” functionality has access to the **Notification Center** with **Settings**, where various notification rules can be configured.




Active	Name	Type	Rule count
<input type="checkbox"/>	StriSchV Anlage 14 I+II (1)	Study based notifications	4
<input type="checkbox"/>	StriSchV Anlage 14 I (2) a	Study based notifications	2
<input type="checkbox"/>	StriSchV Anlage 14 II (2+3) a	Study based notifications	1
<input checked="" type="checkbox"/>	German Group Patients	Study based notifications	1
<input type="checkbox"/>	Test German Law Series Level	Series based notifications	2

Customized groups with specific rules can be configured here. These rules can be based on study, series or patient parameters. The user has also the possibility to make a notification to overrule the compliance alerts. Compliance monitoring will still apply, but the specified alert type will be overruled.



To configure a new notification center group, follow the steps below:

1. Access **Settings** → **Notification Center**
2. Click on **Add** in the upper-right corner to start a new group
3. Indicate the notification group description and select a **Group Type**
4. Mark it as “Active” for the rules to apply to incoming/existing studies
5. Select a **Scope** by indicating the devices it will affect



Active	Name	Type
<input type="checkbox"/>	StriSchV Anlage 14 I+II (1)	Study based notifications
<input type="checkbox"/>	StriSchV Anlage 14 I (2) a	Study based notifications
<input type="checkbox"/>	StriSchV Anlage 14 II (2+3) a	Study based notifications
<input checked="" type="checkbox"/>	German Group Patients	Study based notifications
<input type="checkbox"/>	Test German Law Series Level	Series based notifications

Description

Name:

Description:

Group type: Group is active


Scope

Apply rules to all devices
 Apply rules to selected devices

Device name:
 contrast test
 TOSHIBA_MEC_toyo
 Toshiba Platinum43

Notification Center

- if rules should be applied to all the devices, select **Apply rules to all devices** or **Apply rules to selected devices** and click on **Edit filter**
- A window will pop-up allowing the selection of specific devices
- The modality/modalities can be selected, then the >> button to select the desired devices from the resulting sites
- Finally, click on **Save**



Device filter configuration


6. The rules for the group can now be defined. To add a rule, click on the green button **+Add rule** and fill in the necessary parameters and values. Once a group is created (don't forget to click on the green disk icon to save), the rule will apply to new studies. Pay attention, because creating duplicate rules will cause duplicate notifications.



- Details of the notification can be configured by clicking on **Show notification configuration** (e.g. the severity, supervision status, message, notification key and status). By clicking on **Actions**, email alerts can be configured.

	<input checked="" type="checkbox"/> Acceptable	<input checked="" type="checkbox"/> Unacceptable
Notification severity	ORANGE	RED
Notification supervision	Supervision advised	Supervision obliged
Notification message	200% DRW Unterschreitung	200% DRW Überschreitung
Notification key	Anlage 14 I + II (1) DFP	Anlage 14 I + II (1) DFP
Notification status	Unresolved	Unresolved
Notification action	Actions	Actions

When notifications with advised/obliged supervision are created, an orange/red eye symbol will appear in the study list in a column called *Justification*.



The following table summarizes the different rules (comparison types) available in the Notification Center, in which levels (group types) they can be found and to which parameters they can be applied:

Rule in DOSE	Study	Series
% fixed value	CTDI, DAP, DLP, Dose at ref. Point, Effective dose, Entrance dose, Fluoro time, OrganDose [MG],	CTDI, Compression Force, DAP, DLP, Deviation Index, Effective dose,

	Organ doses, PSD, Radioactivity Administered, Radioactivity Administered per kg, Weighted CTDIvol	Exposure Index, OrganDose [MG], Organ doses
Equals	CTDI, DAP, DLP, Dose at ref. Point, Effective dose, Entrance dose, Fluoro time, OrganDose [MG], Organ doses, PSD, Radioactivity Administered, Radioactivity Administered per kg, Weighted CTDIvol	CTDI, Compression Force, DAP, DLP, Deviation Index, Effective dose, Exposure Index, OrganDose [MG], Organ doses
Range fixed value	CTDI, DAP, DLP, Dose at ref. Point, Effective dose, Entrance dose, Fluoro time, OrganDose [MG], Organ doses, PSD, Radioactivity Administered, Radioactivity Administered per kg, Weighted CTDIvol	CTDI, Compression Force, DAP, DLP, Deviation Index, Effective dose, Exposure Index, OrganDose [MG], Organ doses
German Regulation - Group of patients [Annex 14 I/II.1]	DAP, DLP, Organ Dose [MG], Radioactivity Administered, Radioactivity Administered per kg, Weighted CTDIvol	CTDI, DAP, DLP, Organ Dose [MG]
German Regulation - CT CTDI Single patient [Annex 14 I/II.2a]	CTDI	CTDI
German Regulation - DR DAP Single patient [Annex 14 I/II.2a]	DAP	DAP
German Regulation - XA DAP Single patient [Annex 14 II.2/3a]	DAP	DAP
German Regulation - Eff. Dose Single patient [Annex 14 I/II.2a]	Radioactivity administred	N/A

Highest organ dose	Organ doses	Organ doses
Group notifications	Notification	
Joint commission Fluoroscopy	DAP, Dose at ref. Point, PSD	
Available in list	Patient ID, Patient Name, Study description	Filter [MG], XRay target [MG]
Text contains	Patient ID, Patient Name, Study description	Filter [MG], XRay target [MG]
Text ends with	Patient ID, Patient Name, Study description	Filter [MG], XRay target [MG]
Text equals	Patient ID, Patient Name, Study description	Filter [MG], XRay target [MG]
Text starts with	Patient ID, Patient Name, Study description	Filter [MG], XRay target [MG]

Besides Study and Series, there are other 2 special group types: Patient and Compliance override:

	Patient	Compliance override
Cumulative value	Effective dose	
Compliance monitoring		Administered Radioactivity, CTDI, DAP, DLP, Dose at ref. Point, Effective dose, Entrance dose, Event count, Fluoro time, MGD, Organ Dose [MG], PSD, Prepared Radioactivity, SSDE, Scanlength


For further information, refer videos within the **Notification Center** section of our online training center.

7.3.1. Description of rules

7.3.1.1. % FIXED VALUE

In this rule, the user can enter a reference value and the different percentages that will be used for the different notifications.

For example, in the image below, the reference is set to 100dGy·cm², the achievable levels to 5% and the acceptable levels to 10%.



	Lower acceptable [%]	Lower achievable [%]	Reference (absolute value) [dGy·cm ²]	Upper achievable [%]	Upper acceptable [%]
	10.0	5.0	100.0	5.0	10.0
Rule 1					
	<input type="checkbox"/> Achievable	<input checked="" type="checkbox"/> Acceptable	<input checked="" type="checkbox"/> Unacceptable		
Notification severity	GREEN	ORANGE	RED		
Notification supervision	No supervision	Supervision advised	Supervision obliged		
Notification message	achievable DAP (% fixed value)	acceptable DAP (% fixed value)	unacceptable DAP (% fixed value)		
Notification key:	DAP_RANGE_PERC_FIXED_VALU	DAP_RANGE_PERC_FIXED_VALU	DAP_RANGE_PERC_FIXED_VALU		
Notification status	Resolved	Unresolved	Unresolved		
Notification action	Actions	Actions	Actions		

That means that the resulting levels will be as follows:

Lower Acceptable=90dGy·cm²
 Lower Achievable= 95dGy·cm²
 Upper Achievable= 105dGy·cm²
 Upper Acceptable=110dGy·cm²

This way, when a study arrives, the system will check the selected parameter (e.g. DAP) against all these levels, following this logic:

If the parameter is

1. above the upper acceptable (110dGy·cm²) or
 2. below the lower acceptable (90dGy·cm²) levels,
- it will create an "Unacceptable" notification, that by default is set to color RED.

If the parameter is


- c. between the upper acceptable (110dGy·cm²) and the upper achievable (105dGy·cm²), or
- d. between the lower acceptable (90dGy·cm²) and the lower achievable (95dGy·cm²)

It will create an "Acceptable" notification that by default is set to ORANGE.

If the parameter is between the lower achievable and the upper achievable, then the notification will be marked as "Achievable" with default color GREEN. Normally this notification is deactivated by default.

NOTE: If any of the above levels is set to zero, it will not generate notifications associated to it.

You have a summary of this rule in the following infographic:



7.3.1.2. EQUALS

If the received value

- a) equals the reference value, it will create a 'Matches' notification (e.g. GREEN)
- b) if it is not the case, it will generate a 'No match' notification (e.g. RED)


It is possible to just activate only one of the notifications. Like you can see in this example:

CTDI Equals 100.0 Reference [mGy]

Rule 1

	<input type="checkbox"/> Matches	<input checked="" type="checkbox"/> No match
Notification severity	GREEN	ORANGE
Notification supervision	No supervision	Supervision obliged
Notification message	acceptable CTDI	unacceptable CTDI
Notification key:	CTDI_SERIES_NUMBER_EQUALS	CTDI_SERIES_NUMBER_EQUALS
Notification status	Resolved	Unresolved
Notification action	Actions	Actions

Below you have an infographic that summarizes this rule:



7.3.1.3. RANGE FIXED VALUE

This rule is similar to the “% fixed value”, but instead of percentages, you need to enter absolute values. Here the “Reference” value can be any value between the lower ACH and the upper ACH, but no calculation will be done with this value.

DAP Range fixed value 50.0 75.0 100.0 125.0 150.0

Lower accep < Lower achievable Reference [dGy.cm²] Upper achievable Upper acceptable

Rule 1

	<input type="checkbox"/> Achievable	<input checked="" type="checkbox"/> Acceptable	<input checked="" type="checkbox"/> Unacceptable
Notification severity	GREEN	ORANGE	RED
Notification supervision	No supervision	Supervision advised	Supervision obliged
Notification message	achievable DAP	acceptable DAP	unacceptable DAP
Notification key:	DAP_RANGE_FIXED_VALUE	DAP_RANGE_FIXED_VALUE	DAP_RANGE_FIXED_VALUE
Notification status	Resolved	Unresolved	Unresolved
Notification action	Actions	Actions	Actions

So, in the example above, the values would be the following:

Lower Acceptable=50dGy.cm²

Lower Achievable= 75dGy·cm²
 Upper Achievable=125dGy·cm²
 Upper Acceptable=150dGy·cm²

And the value of 100 dGy·cm² is purely indicative.


As the previously described rule, there are 3 possible scenarios

1. If the parameter is
 - e. above the upper acceptable (150dGy·cm²) or
 - f. below the lower acceptable (50dGy·cm²) levels,
 it will create an "Unacceptable" notification, that by default is set to color RED.

2. If the parameter is
 - g. between the upper acceptable (150dGy·cm²) and the upper achievable (75dGy·cm²), or
 - h. between the lower acceptable (50dGy·cm²) and the lower achievable (75dGy·cm²)
 It will create an "Acceptable" notification that by default is set to ORANGE.

3. If the parameter is between the lower achievable and the upper achievable, then the notification will be marked as "Achievable" with default color GREEN. Normally this notification is deactivated by default.


Below the summary of this behaviour:



7.3.1.4. JOINT COMMISSION REQUIREMENTS FOR FLUOROSCOPY

The notification rules can be configured to look in a certain time period if a patient had an exam with a chosen threshold. In the image below you can see how we implanted the

Join Commission (JC) requirements for Fluoroscopy. In this example, every time a new study comes in the system, there is an automatic checking to see if there are other studies for the same patient and then the cumulative parameters (PSD, DAP, Dose at Reference Point) are checked to see if they fulfil the criteria to report to JC.



JC requirements as deployed by DOSE:

- Peak skin dose:
 - in 2 months, PSD > 3Gy → orange notification
 - in 6 months, PSD > 15Gy → red notification
- Dose at reference point: same rules as PSD (some interventional systems may be lacking RDSR, so there needs to be a second check)
- Dose area product:
 - in 2 months, DAP > 300Gy·cm² → orange notification
 - in 6 months, DAP > 1500Gy·cm² → red notification
- If one of the three above criteria produces an orange notification, then there is a orange group notification generated. If one of the three is red, then group notification is red. (To configure this, check the next rule explanation "Group notifications").

These rules are easily configurable by the user, just add the parameters DAP, Dose at RP and PSD with the comparison types "Joint Commission Fluoroscopy" and modify the applicable values, if needed.

Notification rules

Parameter	Comparison type	Configuration			
DAP	Joint commission Fluoroscopy	Interval 1 [months] 2	Interval 1 reference value [Gy.cm²] 0.4	Interval 2 [months] 6	Interval 2 reference value [Gy.cm²] 0.6
Dose at ref. point	Joint commission Fluoroscopy	Interval 1 [months] 2	Interval 1 reference value [Gy] 4.0	Interval 2 [months] 6	Interval 2 reference value [Gy] 6.0
Peak skin dose	Joint commission Fluoroscopy	Interval 1 [months] 2	Interval 1 reference value [Gy] 4.0	Interval 2 [months] 6	Interval 2 reference value [Gy] 6.0

Peak skin dose | Joint commission Fluoroscopy

Interval 1 [months] 5 | Interval 1 reference value [Gy] 3.0 | Interval 2 [months] 7 | Interval 2 reference value [Gy] 8.0

Interval 1

	<input type="checkbox"/> Complies	<input checked="" type="checkbox"/> Does not comply	<input type="checkbox"/> No data
Notification severity	GREEN	ORANGE	BLUE
Notification supervision	No supervision	Supervision advised	Supervision advised
Notification message	JC Peak skin dose Interval 1 compli	JC Peak skin dose Interval 1 does n	JC Peak skin dose no interval 1 date
Notification key	JC_FLUORO_PSD	JC_FLUORO_PSD	JC_FLUORO_PSD_NO_DATA
Notification status	Resolved	Unresolved	Unresolved
Notification action	Actions	Actions	Actions

7.3.1.5. GROUP NOTIFICATIONS

These rules work with an AND/OR logic and can be used to create notification escalations.

In the following example, a group notification has been configured to comply with the second part of the Joint Commission Requirement.



Second part of the Joint Commission Requirement

Notification: Group notifications

Notification key: JC_FLUORO_DAP | Severity: ORANGE


Notification key: JC_FLUORO_DOSE_AT_REF_POINT | Severity: ORANGE

Notification key: JC_FLUORO_PSD | Severity: ORANGE

Group of notifications matching

	<input type="checkbox"/> No match	<input checked="" type="checkbox"/> Matches
Notification severity	GREEN	ORANGE
Notification supervision	No supervision	Supervision advised
Notification message	Group notifications rule did not matc	Group notifications rule matched
Notification key	GROUP_NOTIFS	GROUP_NOTIFS
Notification status	Resolved	Unresolved
Notification action	Actions	Actions

Group notification ORANGE if any of the selected notifications (keys) is ORANGE



Notification key: JC_FLUORO_DAP Severity: RED

Notification key: JC_FLUORO_DOSE_AT_REF_POINT Severity: RED

Notification key: JC_FLUORO_PSD Severity: RED

Group of notifications matching

No match Matches

Notification severity: GREEN RED

Group notification RED if any of the selected notifications (keys) is RED

To create a group notification, select "Notification" as parameter and "Group notifications" as comparison type. Then choose an AND/OR logic.




Notification rules

Parameter: Notification Comparison type: Group notifications Logic: AND

+ Add notification key

Then click on **+Add notification key** and select a key from the drop-down menu. The notifications whose keys you can select are the ones inside the same notification group. That means, for example, that if you want to create the group notification (scalation) for the Joint Commission rules, you need to create the individual rules before, in the same group.



DAP Joint commission Fluoroscopy Interval 1 [months] 2 Interval 1 reference value [Gy.cm²]

Peak skin dose Joint commission Fluoroscopy Interval 1 [months] 2 Interval 1 reference value [Gy] 3.0

Dose at ref. point Joint commission Fluoroscopy Interval 1 [months] 2 Interval 1 reference value [Gy] 3.0

Notification Group notifications

Notification key: JC_FLUORO_DAP


Notification key: JC_FLUORO_PSD

Notification key: JC_FLUORO_DOSE_AT_REF_POINT

Once you have selected the keys, select the “Severity” and click on “Show notification configuration”.




There you will be able to choose the severity of the group notification.




7.3.1.6. TEXT-BASED NOTIFICATIONS


The parameters Patient ID, Patient Name and Study Description (at study level) and Filter [MG] and XrayTarget [MG] (at series level) can be tracked with text-based notifications. The possible rules are:

- Available in list: where you can click on Edit values and write down a list of possible names.





- Text contains/ends with/equals/starts with



You can ignore Upper-/Lowercase at any time by checking the corresponding box.

If you click on “Show notification configuration” you can choose what will happen if the study “Matches” or “Not matches” the entered text.

Rule 1

	<input type="checkbox"/> Matches	<input checked="" type="checkbox"/> No match
Notification severity	GREEN	ORANGE
Notification supervision	No supervision	Supervision obliged
Notification message	acceptable Patient name	unacceptable Patient name
Notification key:	PATIENT_NAME_TEXT_STARTS_V	PATIENT_NAME_TEXT_STARTS_V
Notification status	Resolved	Unresolved
Notification action	Actions	Actions

Notification when "No match" is produced

7.3.1.7. GERMAN REGULATION NOTIFICATIONS

Normally, if you are a German customer, all the notifications rules will be already configured for you, ordered by groups and with references to the text in the regulation.

- StriSchV Anlage 14 II (2+3) a Study based notifications
- StriSchV Anlage 14 I (2) a Study based notifications
- StriSchV Anlage 14 I + II (1) Study based notifications
- StriSchV Anlage 14 I (2) a Nuklearmedizin Study based notifications

Description

Name: StriSchV Anlage 14 I + II (1) Group type: Study based notifications

Group is active

Description: Bezogen auf eine Gruppe von Personen

Jede Überschreitung des Mittelwertes über die letzten 20 aufeinanderfolgenden Untersuchungen gleicher Untersuchungsart um mehr als 100 Prozent des jeweiligen diagnostischen Referenzwertes, sobald der diagnostische Referenzwert einer einzelnen Untersuchung um 200 Prozent überschritten wurde.

Scope

Apply rules to all devices

Apply rules to selected devices

Device name

DX-D100

DXD600ER


Hosp5 - SIEMENS SOMATOM Force

Notification rules

Parameter	Comparison type	Configuration						
DLP	German Regulation - Group of patk	Reference study count	20	Single deviation [%]	300.0	Group deviation [%]	200.0	?
TOTAL CTDI	German Regulation - Group of patk	Reference study count	20	Single deviation [%]	300.0	Group deviation [%]	200.0	?

German regulation notification groups

Nevertheless, these rules are useful to monitor outliers and they could be interesting for other users as well. If you wish to configure them by yourself, just add all the rules that contain « German Regulation » in their comparison type for each parameter.



7.3.1.7.1. GERMAN REGULATION - GROUP OF PATIENTS [ANNEX 14 I/II.1]

This rule works following the logic presented in this infographic:



Thus, 3 notifications are created:

- 1) an orange notification is created when the dose parameter is above 300% of the DRL (deviation of 200%)
- 2) A second orange notification is created if the average value of that dose parameter for the last 20 studies of the same type (study group) is above the 200% of the DRL (deviation of 100%).
- 3) If both cases happen, then a red notification will be created.

Note that for this rule no “Group notification” is needed, since the scalation is already included in this special rule. You can check this by clicking on “Show notification configuration”.

Combination


	<input type="checkbox"/> Not two warning	<input checked="" type="checkbox"/> Both rules give warning
Notification severity	GREEN	RED
Notification supervision	No supervision	Supervision obliged
Notification message	Kein bedeutsames Vorkommnis Anlage	Bedeutsames Vorkommnis Anlage 1
Notification key:	StriSchV Anlage 14 I (1) DLP	StriSchV Anlage 14 I (1) DLP
Notification status	Resolved	Unresolved
Notification action	⚙️ Actions	⚙️ Actions

Of course, if you wish to use this rule with other values, you can change them:


DLP	German Regulation - Group of patie	Reference study count 20	Single deviation [%] 300.0	Group deviation [%] 200.0	?
---	--	--	--	---	---

For more information about this rule and the criteria to select the last 20 examinations in both STUDY and SERIES level, please refer to the document information document ***GERMAN LAW NOTIFICATIONS - Group of patients [Annex 14 I/II. 1]***.

When handling this type of notification, a « Show references » button will appear and by clicking on it, the user will be able to see the last 20 studies that were used for the average calculation.




Show references for a German group notification




List of studies used for the German group notification

7.3.1.7.2. GERMAN REGULATION - CT CTDI SINGLE PATIENT [ANNEX 14 I/II.2A]

With this rule, the user can establish two different reference values for CTDIvol: one for the head and one for the rest of the body:



If any of the two values is exceeded for the corresponding phantom (16 or 32cm), a notification can be created. Here you can see the example with the values of the German Legislation:



This rule has a slightly different operation depending on whether it is configured in “Study-based notifications” or “Series-based notifications”.

a) Study-based notifications

The system will check the phantom of each series and will sum up the values for each phantom. So, the CTDI for

- i) head will be the sum of all the CTDIs of the series whose phantom is 16cm.
- ii) body will be the sum of all the CTDIs of the series whose phantom is 32cm.

Each one of these values will be checked against their corresponding reference one, and if they exceed it, a notification will be created.

b) Series based notifications


The system will compare individually the values from each series against the reference value, so

- i) each series with phantom 16cm will be checked against the reference value for head, if the value is higher than the reference, a notification will be created.
- ii) each series with 32cm phantom will be checked against the reference value for body, if the value is higher than the reference, a notification will be created.


Exclusion rules

It is possible to restrict this rule for CTDI on series level to specific series, based on the following parameters: Acquisition Protocol, Acquisition type, Series Description, Study Description-Acquisition Protocol, Study Description-Series Description.

There is a blue ‘+’ button next to the rule configuration:




After the blue '+' is clicked, the configuration appears, allowing the user to select rule restrictions for several parameters. These restrictions can be based on a list of values or on parts of text:



This way, constant angle acquisitions can be excluded from this rule, to not generate undesired notifications.


7.3.1.7.3. GERMAN REGULATION - DR DAP SINGLE PATIENT [ANNEX 14 I/II.2A]

This rule allows to set a limit for the DAP in general radiology devices:



The same way as in the previous section, **exclusion rules** can be added on study level (also possible for the XA DAP Single Patient Annex 14 I/II.2/3A described in the next section).

There is a blue '+' button next to the rule configuration, labelled as "1" below.




After the blue '+' is clicked, the configuration labelled above as "2" appears, allowing the user to select rule restrictions for event count, protocol name and study description.




In case of event count, the user can choose a comparison operator of less than, greater than, equal to, or not equal to. The desired value can be input as shown above.

In the case of protocol name and study description, the rule allows to enter a list of names or part of a name:




7.3.1.7.4. GERMAN REGULATION - XA DAP SINGLE PATIENT [ANNEX 14 II.2/3A]

This rule is summarized in the following infographic:



- If the DAP is higher than the first reference value (e.g. 20.000cGy·cm²), an ORANGE alert will be created.
- If the DAP is higher than the second reference value (e.g.50.000Gy·cm²), another ORANGE alert will be created.
- If the second warning is created and there are manual comments about deterministic effects within a certain number of days (e.g. 21 days), a RED notification will be created.

The values of the reference levels and the period can be adjusted.

DAP ref value [cGy*cm ²]	20000.0	DAP ref value Level 2 [cGy*cm ²]	50000.0	Period [day]	21	
--------------------------------------	---------	--	---------	--------------	----	---

Additional exclusion rules can be configured on study level, as described in [GERMAN REGULATION - DR DAP SINGLE PATIENT \[ANNEX 14 I/II.2A\]](#) above.

7.3.1.7.5. GERMAN REGULATION - EFF. DOSE SINGLE PATIENT [ANNEX 14 I/II.2A]

This rule refers only to the effective dose and organ doses in Nuclear Medicine.

Notifications for effective dose and organ doses are calculated based on study groups and on injected activity thresholds, as suggested by the Bundesamt für Strahlenschutz.


Tabelle 1 Aktionsschwellen für die zu Untersuchungszwecken verwendeten Radiopharmaka, für die DRW festgelegt worden sind

Radio-pharmakon	Untersuchung	Limitierende Dosis	Dosis-koeffizient	Aktions-schwelle	Aktions-schwelle
^{99m} Tc-Pertechnetat	Schilddrüse-Szintigraphie	Effektive Dosis	9,7E-03 ¹	2000 ²	
^{99m} Tc-MDP/DPD/HDP	Skelett-Szintigraphie	Organdosis (Harnblase)	2,8E-02 ¹		55 ³
^{99m} Tc-Sestamibi	Perfusion/Vitalität-Szintigraphie - Ruhe	Effektive Dosis	6,9E-03 ¹	3300 ²	
	Perfusion/Vitalität-Szintigraphie - Stress	Effektive Dosis	6,0E-03 ¹	3700 ²	
	Nebenschilddrüse-Szintigraphie	Effektive Dosis	6,9E-03 ¹	3400 ²	

Example of action levels for administered radioactivity by the Bundesamt für Strahlenschutz


The mapping of the study groups and the thresholds are configurable.

This rule can be found under group type "Study based notifications" and it's only available for the "Radioactivity administered" parameter.




German rule for NM effective and organ doses based on Radioactivity administered.

In "Configure", users can import the csv file with the action level values provided by the BfS for each study group by clicking on "Import".



Click on "Import".

This file can be provided by Qaelum.



Action levels based on pre-existing study groups.

Click on “Commit” to visualize the values:

Map study groups to Eff.Dose Single patient [Annex 14 I/II.2a] triggering values

Study group name	Eff.Dose Single patient [Annex 14 I/II.2a] triggering value [in MBq]	Version
Schilddrüse - Szintigraphie - [99mTc]Per technetat	2000.0	2022
Skelett - Knochenszintigraphie - [99mTc]MDP/DPD/HDP	55.0	2022
Herz - Perfusion/Vitalität - [99mTc]Sestamibi - Zweitagesprotokoll	3300.0	2022
Herz - Perfusion/Vitalität - [99mTc]Sestamibi - Eintagesprotokoll	3700.0	2022
Herz - Perfusion/Vitalität - [99mTc]Tetrofosmin - Zweitagesprotokoll	3800.0	2022
Herz - Perfusion/Vitalität - [99mTc]Tetrofosmin - Eintagesprotokoll	4100.0	2022
Nieren - Funktionsszintigraphie - [99mTc]MAG3	1900.0	2022
Lunge - Ventilation - [99mTc]DTPA	4200.0	2022
Lunge - Ventilation - [99mTc]Technegas	1400.0	2022
Lunge - Perfusion nach Ventilation - [99mTc]MAA	1900.0	2022
Lunge - Perfusion ohne Ventilation - [99mTc]MAA	1800.0	2022
Gehirn - DAT-SPECT - [123I]FP-CIT	1000.0	2022
Gehirn - PET - [18F]FDG	19.0	2022
Gehirn - PET - [18F]FET	21.0	2022
Nebenschilddrüse - Szintigraphie - [99mTc]Sestamibi	3400.0	2022

[Import](#)

Imported values

Close this window and save both the rule and the notification group.

7.3.1.7.6. HIGHEST ORGAN DOSE

This rule refers to Nuclear Medicine devices as well, it allows to set a limit for the highest organ dose.

Organdoses Highest organ dose OrganDose reference value [mGy] ?

- If the notification is created as STUDY-based, it will compare the maximum value of the organ doses at STUDY level (sum of doses for each organ at series level)
- If the notification is created as SERIES-based, it will compare the maximum value or the organ doses of each individual SERIES.

7.3.1.8 CUMULATIVE VALUE

This rule is only available in “Patient-based notifications” and for Effective dose (all modalities).

Effective dose Cumulative value Interval 1 [months] Interval 1 reference value [mSv] Interval 2 [months] Interval 2 reference value [mSv] ?

It works similarly to the one of Joint Commission but for effective dose.

7.3.1.9 COMPLIANCE MONITORING

This rule is found inside the group type “Compliance Override”. It can be used for any parameter that used in the study groups as a limit. The rule allows you to “override” the severity (color), supervision, message, key and status of the compliance notifications.

Rule 1

	<input type="checkbox"/> Achievable	<input checked="" type="checkbox"/> Acceptable	<input checked="" type="checkbox"/> Unacceptable
Notification severity	GREEN	ORANGE	RED
Notification supervision	No supervision	Supervision advised	Supervision obliged
Notification message	achievable DLP	acceptable DLP	unacceptable DLP
Notification key:	DLP_COMPLIANCE	DLP_COMPLIANCE	DLP_COMPLIANCE
Notification status	Resolved	Unresolved	Unresolved
Notification action			

Please note that compliance notifications are created based on the study group limits with the logic =>, that is, a notification is created if the parameter is equal to or greater than the achievable/acceptable level.

7.3.2. E-mail notifications

E-mail sending can be configured every time a notification is generated. To do so, click on "Actions" and select an e-mail sending job for the user(s).

Does not comply

ORANGE

Supervision advised

JC Peak skin dose Interval 1 does n

JC_FLUORO_PSD

Unresolved

Actions

Select actions from the table

	Name	Action Type	Users
<input type="checkbox"/>	Test	E-mail text	[bojan, Live]


For more information on the configuration of the e-mail sending job, please refer to [Notification configuration](#).

7.3.3. Import and export of notification rules

Notification rules can be imported and exported to save time on configurations or to keep a backup of the already configured rules.


7.3.3.1. EXPORT OF NOTIFICATION RULES

To export notification configurations from the **Notification center**, click on the *Export* button (blue button with a cloud downward icon).




Export notification groups

A popup will appear where the user can choose to export all groups or select individual groups.



Select groups to export

By clicking on *Export* (bottom right of the window), the download will begin.




Download window


By clicking on **Download**, the export will be downloaded to the predefined export folder according to browser settings (by default C:\Users*username*\Downloads) as a json file.

7.3.3.2. IMPORT OF NOTIFICATION RULES


To import previously exported rules in the **Notification center**, click on the **Import** button (blue button with a cloud upload arrow).




A popup will appear allowing to select a previously exported json file for import.



Select file to import



After clicking **Next**, the user can choose to import all groups or select individual groups.



Begin the download by clicking on **Import** (bottom right of the window). Additional configurations may be required if the export file originated from a different DOSE installation. Please find more details [below](#).

7.3.3.2.1. IMPORT FROM THE SAME DOSE INSTALLATION

When importing from the same DOSE installation, no additional configuration is required: all notification groups will be saved automatically, and the rules will be effective immediately.

7.3.3.2.2. IMPORT FROM A DIFFERENT DOSE INSTALLATION

When importing from another DOSE installation, additional configuration may be required, depending on the rule customization at the source:


A) ORIGINAL RULES WERE APPLIED TO ALL DEVICES AND NO EMAIL ACTION WAS SET UP

If the exported notification groups were applied to all devices (see [Notification Center](#)) and no email actions were configured for the rules (see [E-mail notifications](#)), no additional configuration is required. All notification groups will be saved automatically, and the rules will be effective immediately.


B) ORIGINAL RULES WERE APPLIED TO SELECTED DEVICES AND/OR EMAIL ACTIONS WERE SET UP

If the exported notification groups were applied to specific devices (see [Notification Center](#)) and/or email actions were configured for the rules (see [E-mail notifications](#)), then the following steps should be followed for notification group import:


1. Import one group at a time.
2. Once the group is imported, change the scope to *Apply rules to all devices*.




3. Click on *Show notification configuration* and activate the editing mode.




4. For each rule, verify that no email action is activated.



5. Save and repeat for each rule.



6. Once all rules have been imported, save the notification group.




7. Repeat for any other groups that need to be imported.

8. Device level

The device dashboard is represented as a tree with two levels: Site (e.g. hospital, department, group of devices) and Device.

Each device can have one or more data panels: users with edit permissions can choose the available panels in the Edit device section.




Tree structure on Device level


8.1. Site section

8.1.1. Add site

Users with the *Site management* functionality within their role can add a site by clicking on the **Add location** button. This button is located on the left side of the view and on the top of the site tree. The new site needs to consist of a name and a description. The site name needs to be unique in the application.




Configuration and Add location buttons



Register new site

A 3-character “tag” (see image above) can be added to easily identify the hospital in **Modality** comparisons.

By using the filter or grouping tool in the **Configuration** button, the tree structure overview will be clearer.



Filter and grouping screen

8.1.2. Edit site


A user with editing permissions can modify a site by clicking “Edit Site” following the “actions” button, which is located on the right-upper side of the view. Site ID is not editable.

8.1.3. Remove site

The remove button is located in the “Edit Site” window on the upper-left side. Removing a site will lead to removing all the site’s information and all the devices that apply to that site but will not erase dose data. User needs to have modification rights.

8.1.4. Dashboard

In this dashboard, the general overview of the site is shown: the user can find information about the amount of analyzed studies, the total effective dose for the devices, the 10 last analyzed studies and the list of the 10 highest dose studies.



Dashboard | Device list | Warning | Logbook | Dynamic reports

17526 # studies were analysed

93069.92 D_{eff} total dose given

10 most recent studies overview

STUDY DATE	STUDY	PATIENT	DOSE [DEFF]
2015-10-07 08:08	Coronary^Diagnostic Coronary Catheterization	F - 75 Y	152.51
2015-10-05 08:13	Electrophysiology^Electrophysiology study	M - 71 Y	4.34
2015-10-01 07:37	Coronary^Diagnostic Coronary Catheterization	M - 64 Y	74.08
2015-09-30 07:21	Coronary^Diagnostic Coronary Catheterization	M - 84 Y	87.26
2015-09-29 07:30	Coronary^Diagnostic Coronary Catheterization	M - 72 Y	69.74
2015-09-28 10:47	Electrophysiology^Electrophysiology study	F - 51 Y	0.86
2015-09-28 08:02	Coronary^Diagnostic Coronary Catheterization	M - 62 Y	41.53
2015-09-25 08:19	Coronary^Diagnostic Coronary Catheterization	F - 65 Y	41.58
2015-09-24 11:50	Coronary^Diagnostic Coronary Catheterization	M - 83 Y	76.68
2015-09-24 07:22	Electrophysiology^Electrophysiology study	M - 48 Y	0.29


10 highest dose studies overview

STUDY DATE	STUDY	PATIENT	DOSE [DEFF]
2013-12-01 11:28	RAD ct angio 11,[50-3-2-174]	F - 63 Y	223.54
2013-12-01 12:33	RAD ct angio 11,[50-3-2-174]	F - 47 Y	197.89
2013-11-07 22:51	RAD ct comb 03,[18-2-2-194]	F - 39 Y	152.25
2013-12-15 10:32	RAD ct abd 22,[1-2-2-187]	F - 72 Y	130.99
2013-09-25 12:23	RAD ct angio 06,[50-3-2-196]	M - 55 Y	115.59
2013-12-17 21:47	RAD ct angio 11,[50-3-2-174]	M - 85 Y	115.34
2013-11-28 15:35	RAD ct angio 11,[50-3-2-174]	M - 52 Y	113.49
2013-09-25 07:11	RAD ct thorax 24,[10-11-2-97]	M - 39 Y	111.58
2013-10-19 20:28	RAD ct angio 11,[50-3-2-174]	M - 66 Y	109.41
2013-11-05 12:02	RAD ct angio 05,[50-3-2-196]	M - 64 Y	109.17

General overview of studies performed on a specific device

8.1.5. Device list

List with all the devices available at the chosen site.



Dashboard | Device list | Warning | Logbook | Dynamic reports

NAME	MODALITY	VENDOR	MODEL	COLLECTION ID	#STUDIES
Room 03 - FujiFilm DEVO	DR	FujiFilm	DEVO	DEVO	622
Room 04 - Toshiba Aquilion 16	CT	TOSHIBA	Aquilion	ID_STATION_RADIODIAGNOSTICA 3 AOUPISA	406
Room 01 - Inspiration	MG			INSPIRATION01	192
Room 02 - Siemens Somatom Sensation 64	CT	SIEMENS	Sensation 64	siemens_sensation_64	116
Room 06 - Siemens Axiom Artis	XA	Siemens	AXIOM-Artis	AXIS05070	116
Room 05 - Siemens Somatom Definition Flash	CT	SIEMENS	SOMATOM Definition Flash	Def. Flash Z13	0
Room 07 - Siemens Fluorospot Compact FD	RF	SIEMENS	Fluorospot Compact FD	siemens_Fluorospot_Compact_FD	0
Room 08 - NM Hotlab	NM	Veenstra	Hotlab	nm_device	0

9 rows

Toggle Filter Bar visibility | Reset filter | Actions

Device list

8.1.6. Dynamic reports



Dynamic reports provide users with periodic reports of any dashboard in DOSE. This report can be scheduled to be created on a regular interval (days/weeks/months/...) and

contain data that is up to date for the specified time period. This allows users to monitor dashboard data without having to log in to DOSE and find or re-configure the required dashboards.

Dynamic reports provide two report options: a URL link to the dynamic report, and a static PDF report. The dynamic report is sent as a URL that will open the interactive report data within DOSE. The generated PDF of the report is sent as an email attachment and is static.

The remainder of this document explains how dynamic reports can be created. For certain actions, only specific roles have permissions to perform the required changes.

The dynamic report creation tab can be found under the upper-right action button on **Device** level.




Please follow the steps below to create a customized dynamic report:

1. The first step is to create a new dashboard (or use existing dashboard on device level), with the graphs presenting the data that you would like to receive in your report. The following two aspects are important before creating the dynamic report":
 - a. The date range. As the Dynamic Reports are often requested to be send on a periodic time interval, it is advised to put the DOSE daterange equal to this time interval. This way the report will always contain complete up-to-

date information for this time range (e.g. for a weekly report it is advised to put the time range on “last week”). If overlap of data is preferred, the data range can be chosen otherwise.

- b. The graphic filter. Each graphic allows you to put a specific graphic filter. Configure this filter correctly so that the data in the chart corresponds to the data that you want to be visible. The specific chart filter can be configured by clicking on the padlock symbol on the chart. Check that the general panel filter is temporarily configured to show all data. For more information about the graphic and general filters, check the user's manual. In the example, a filter has been configured to show only the effective dose of studies with a combination of specific series, since this type of studies can be the scope of a practical case.
2. After all the graphs in the newly created dashboard are presenting you the requested data, a dynamic report can be created from the current status. This report will save the current filter and date range settings. The filter and the date range can additionally be adjusted/double checked during the creating of the report.



Dynamic report creation [X]

Report name:

Date Filter: [\(Edit\)](#)
 📅 Date range (last 10 years) : 2008/07/18 - 2018/07/19

Data Filter: [\(Edit\)](#)
 Parameter type: DLP | mGy.cm

Schedule report:

Minutes Hourly Daily Weekly Monthly Yearly

Every minute(s)


Recipient(s)
 [Remove user](#)

[Add other user](#)

[Preview report](#)

[Save](#) [Cancel](#)


3. The last step is creating a scheduler job to create and send this report that you have just created. For this, activate the check box “schedule report”.



4. The creation of the job, consist of the following steps:
 - a. Provide a job name
 - b. Check/adjust the date and the filter
 - c. Now select the time interval for which you would like to receive an email with the Dynamic Report
 - d. Add the recipient(s)
 - e. Click on "save"

For more information, refer to the video *How to create Dynamic Reports* in our online training center.

5. Delete the dynamic report/scheduler job
 - a. The created report can be deleted via Settings → Scheduler → search for the right dynamic report → Actions (right below) → Edit job → Remove job



Qaelum

Expert mode | Welcome ECR_EN | Profile | Logout | Switch lights

Start Portal | Study | Patient | Device | Modality | Location | Quality Management System | Contrast | **Settings**

WARNINGS

Warning overview

Scheduler

ENGINE MONITOR

Engine Monitor

Study Analysis

USERS

User Management

DATA TABLES

Conversion Factors

Link Study Description - Bodypart exami

Study Groups

TOOLS

Scheduler (Running)

CREATION DATE	JOB NAME	ACTION NAME	USERNAME	FREQUENCY
2018-06-25 17:16:30	bxco	MailReportJob	qaelum	ADVANCED
2018-06-25 15:46:09	QMS Notifications silent	QmsAlertCalculator	qaelum	ADVANCED
2018-05-29 14:39:20	ssde	SSDECalculation	qaelum	YEARLY
2018-05-29 12:37:17	CTSeriesHandler	CTSeriesHandlerJob	qaelum	YEARLY
2018-04-27 14:21:33	MailReportJob Arnaud's report	MailReportJob	ECR_EN	WEEKLY
2018-04-04 16:15:31	test last studies	AnalysedStudiesDevicesJob	qaelum	YEARLY
2017-11-13 16:10:57	glenn test mail	MailReportJob	qaelum	YEARLY
2017-11-09 11:41:14	duplicate Studie	DuplicateStudiesJob	qaelum	YEARLY
2017-11-08 12:13:52	QMS Notifications no Mail	QmsAlertCalculator	ECR_EN	YEARLY
2017-10-24 11:27:42	Mail Report	MailReportJob	ECR_EN	WEEKLY
2017-09-04 13:59:57	Logbook follow-up reminder	LogbookFollowUpMail	qaelum	Advanced
2017-08-31 15:51:15	Marta	MailReportJob	qaelum	WEEKLY
2017-07-11 15:00:50	mac_hietnev_test	macCalculationHietnevInh	qaelum	YEARLY

22 rows


Toggle Filter Bar v

Add job

Edit job

Trigger now


© 2012-2018 QAELUM




8.2. Device section

8.2.1. Add Devices

A user with editing rights can add a device by clicking the **Add device** button, found in the **Actions** menu in the upper right corner of the site Dashboard. To add a new device the user will need to fill in the device name, collection ID, station name and the DICOM modality. The user can also decide which data panels are shown for the selected device.



Add device in Actions



Add device

Name * Siemens SOMATOM Definition 64G

Description Highly accurate device

Collection ID * 687

Vendor Siemens

Model SOMATOM

Serial number 5642158421

Internal reference

Modality * CT (Computed Tomography)

Station Name * Operation Room

Institution Name

Data panels Dosimetry overview (CT)
 Protocol Monitoring (CT)
 Live dashboard (CT)

Background Dashboard:

Demo Hospital ECR

Site Info

QAE LUM Long term data hospital

Dashboard Device list Warning Log

48645 # studies were analysed

10 last analysed studies overview

PATIENT	STUDY
M - 60 Y	Coronary*Diagnostic Coronary Cathete
M - 57 Y	Coronary*Diagnostic Coronary Cathete
M - 73 Y	Coronary*Diagnostic Coronary Cathete
M - 61 Y	Coronary*Diagnostic Coronary Cathete
F - 20 Y	Electrophysiology*Electrophysiology st
F - 18 Y	Electrophysiology*Electrophysiology st
F - 70 Y	Electrophysiology*Electrophysiology st
F - 69 Y	Coronary*Diagnostic Coronary Cathete
F - 79 Y	Coronary*Diagnostic Coronary Cathete
F - 70 Y	Coronary*Diagnostic Coronary Cathete

studies overview

gio 11,[50-3-2-174]
gio 11,[50-3-2-174]
kz,[12-7-2-51]
mb 02,[18-2-2-194]
mb 02,[18-2-2-194]
mb 02,[18-2-2-194]
d 22,[1-2-2-187]
d 22,[1-2-2-187]
brax 100,[10-11-2-146]
gio 11,[50-3-2-174]

Add new device configuration


8.2.2. Add/Edit station name

During configuration, it is necessary to fill in a primary station name for the device. It is possible to change, to delete or to add other station name for the same device by clicking on the "Edit Station Name" button. This button is located in the "Actions" menu on the right upper corner of the Device main dashboard. A device needs to have at least one station name. User needs to have modification rights to access to this section.

Remark: Please inform us if the station name of a device was changed after maintenance. If the station name is modified and it is not communicated to us, there might be missing data.

8.2.3. Edit device

A user with editing permissions can change device information, specifications, main parameter/units and data panels by clicking on "Edit Device", located within "Actions" in the upper right corner of the Device information page.




Room 03 - FujiFilm DEVO

1868 # studies were analysed | 0 # studies last 7 days | 0 # studies today | 97073.74 dGy cm² total dose given | 2382 # notifications

Study date	Study	Patient	Dose [dGy.cm ²]
27/02/2014	RX HANDEN	M - 65 Y	1.32
27/02/2014	RX SCHOULDERS	M - 65 Y	24.29
27/02/2014	RX THORAX FACE/PROFIEL	F - 57 Y	4.24
27/02/2014	RX RIBBENDOOSTER RECHTS	F - 30 Y	8.12
27/02/2014	RX THORAX FACE	F - 30 Y	1.35
27/02/2014	RX FEMUR RECHTS	M - 46 Y	73.85
27/02/2014	RX KNIJE RECHTS	M - 46 Y	12.25
27/02/2014	RX THORAX FACE/PROFIEL	F - 85 Y	5.77
27/02/2014	RX SCHOULDERS	F - 84 Y	11.35
27/02/2014	RX THORAX FACE/PROFIEL	F - 53 Y	3.06

Device Level



Edit device

Remove device

Name: Room 01 - Inspiration

Description:

Collection ID: INSPIRATION01

Vendor:

Model:

Serial number:

Internal reference:

Modality: MG (Digital Mammography)


Data panels:

- Live dashboard (MG)
- Dosimetry overview (Digital Mammography)

Commit Discard

Edit device

The user may be interested in changing the main parameter and/or units shown for the device for either functional or legal purposes. To do so, click on the “Parameter” dropdown menu to see the possible options, and select the appropriate one.



Change parameter/units

8.2.3.1. SPECIAL CASE: GENERAL (DX, DR, CR) AND INTERVENTIONAL (XA) RADIOLOGY


For these modalities, the units of the total DAP can be changed (however not for acquisition or fluoroscopy DAP). The unit for other parameters is predefined.

Parameter	Units
DAP	dGy·cm ²
	mGy·cm ²
	Gy·cm ²
	cGy·cm ²
	µGy·cm ²
DAP from acquisitions	dGy·cm ²
DAP from fluoroscopy	dGy·cm ²

For more information, refer to the video *How to change parameter units on Device Level* in our online training center.

8.2.4. Remove device

The “Remove device” button can be found in **Actions**, as well as within the “Edit Device” window. Removing a device will lead to the deletion of all device information, but the dose data will not be erased. Modification rights are needed to be able to delete the dose data.



Remove device

8.3. Configuration


The **Configuration** tab can be also found when clicking on the device name.

8.3.1. DICOM configuration

A user with the *Device Integration Manager* role can edit the DICOM configuration for the device.

In this tab, it is possible to change the data source. Also, the absence or the presence of a thumbnail can be selected here.

Note that changing the data source can lead to loss of data. A data source can be a combination of more than one source (i.e. DICOM header and MPPS).



Dashboard Configuration Notifications Logbook

Dicom configuration Correction factors Bodypart mapping


Device configuration (Don't forget to save your changes!)

Thumbnail Action: Thumbnail
Thumbnail Size: 65
Data Source: Radiation Dose Structured report (RDSR)
Dose multiplier: 1.0
Dose information: Dose available

Check Study Protocol name.
 Fill with series acquisition concat if no protocol name is present.
 Fill with N/A if no protocol name is present.
 Keep topograms

Dicom Configuration tab

It is possible to add or replace a DICOM tag. In **Actions**, by clicking on **Add tag**, the user is able to configure the new tag.



Name

Old tag:
Group, Element

New tag:
Group, Element


Level 1 DICOM
Level 2 DICOM
Level 3 DICOM

Data type

Save Close

Add tag

DOSE allows using subtags for this mapping, for example, if we want to replace dicom tag (0008, 1030) with the tag (0008,0100), which is a subtag of 0032,1064 and 0040,0275 respectively, the configuration should look as the following:




The dialog box contains the following fields and options:

- Name:** Study Code in Study Description field
- Old tag:** Group, Element: 0008,1030
- New tag:** Group, Element:
 - 0040,0275 Level 1 DICOM
 - 0032,1064 Level 2 DICOM
 - 0008,0100 Level 3 DICOM
- Data type:** Text
- Buttons:** Save (blue), Close (red)

If only 2 levels are present, level 3 should be empty. If there are no sublevels, only Level 1 should be used.

8.3.2. Correction factors/parameters



The screenshot shows the Qaelum interface with the following elements:

- Header:** Qaelum logo, Expert mode dropdown, Welcome ECR_EN |
- Navigation:** Start Portal, Study, Patient, Device, Modality, Location, Quality Management System, Contrast, Settings
- Left Panel:** Filter, Add Site, tree view including Default_Site, Demo Hospital 1, GE MEDICAL SYSTEMS BrightSpeed, Philips iCT 256, Room 01 - Inspiration Test, and Room 02 - Siemens Somatom Sensation 64 (highlighted).
- Main Content:** Room 02 - Siemens Somatom Sensation 64
 - Device info: Siemens Somatom Sensation 64
 - Tabs: Dashboard, Configuration (highlighted), Logbook, Calendar Overview, Logbook
 - Dicom configuration: Parameters (highlighted), Correction factors (highlighted)
 - Parameter correction, Range correction
 - Dose Correction Factor Configuration table:

LEVEL	PARAMETER	START DATE	END DATE

Correction factor configuration

- The user can add a correction (calibration) factor for dosimetric indices and the corresponding dose values arriving on the system from this point on will be multiplied by this factor.

This configuration can be performed on Device level by clicking on the device in the device tree > Configuration > Correction factors. This dashboard provides 2 types of correction factors: parameter and range correction factors.

- The **parameter correction** will apply a correction factor to a single parameter (e.g. global correction of DAP as DICOM data (0018,115E) is reported in cGy and officially DICOM expects dGy).
 - Using the **range correction**, the user can select parameters that this calibration factor will depend on. The parameter that will change is called a "target" and the parameter that the correction is based on is called a "trigger". For example, a calibration factor can be set for DAP (i.e. the target, called *Target* in the Dose Correction Factor Range Configuration) if tube voltage (i.e. the trigger, called *Parameter*) is in the range of 50-80 kV, and another one if it is in the range of 80-125kV.
- The user can set the correction date range starting from the date of the configuration (no historical data are affected, only new data arriving following the correction configuration). The end date can be set or remain empty.

Warning: this influences the dose calculation and we strongly recommend to contact Qaelum for any change.

Correction factors will only be applied to analyses of new studies, applying correction factors for historic studies is thus not possible through the application. Correction factors on SERIES will result in a recalculation of the total value on STUDY level. For that reason, we advise never to apply a correction factor on both levels. Take also into account that for RDSR, the study values are by default collected from the study values in the RDSR. Having correction factors on series level will then result in the sum of corrected series values being taken as study values.

Please consider that Correction Factors can be Added or Edited, but can't be deleted as this is the only way in which there is track of the applied correction factors to the studies that fall within the start and end date ranges.


8.3.2.1. PARAMETER CORRECTION

As explained above, this functionality offers a correction for a collected parameter on STUDY or SERIES level. To configure a correction factor, simply click on the Actions button in the lower right corner, Add row and make the correct configurations.

- Select STUDY or SERIES level in order to determine on which level the correction should occur. Consider that corrections on SERIES level will automatically result in a sum of the corrected SERIES values to the according STUDY level value.

- After having selected a level, DOSE will provide you with the possible parameters to which a correction factor can be applied. This is modality dependent.
- Optionally you can select a start or end date at which these corrections should be applied.

Last you can fill in the correction factor (multiplication value) that should be used.



8.3.2.2. RANGE CORRECTION

This correction is similar to the above but will only be applied if the specified trigger is valid. For example, the SERIES have a trigger parameter "kV" with a "Range start" 80 and a "Range End" 90 with a target parameter "DAP". In this case the correction factor will be applied for any series of this device with a value in the range of 80-90kV.

Please follow the steps below for range correction configuration:


- Select STUDY or SERIES level to specify on which level the correction should occur. Trigger and Target parameters will always be set up for the same level.

NOTE: Consider that corrections on SERIES level will automatically result in the sum of the corrected SERIES values for the corresponding STUDY level value when the data source is DICOM Header.

- After having selected a level, select a trigger parameter from the possible parameters offered by DOSE.
- Then select a target parameter from the dropdown menu. The target is the parameter to which the correction factor should be applied when the trigger falls within the configured range.
- Assign a start (lower limit) and end range (upper limit).
- Optionally a start and/or end date can be selected to limit the date range for which these corrections should be applied.

NOTE: Remember that correction factors will only be applied to analyses of new studies, thus applying correction factors for historic studies is not possible through the application.

- For RF and XA devices, a specific detector can be selected in the 'Detector' field to limit the correction to Single Plane, Plane A or Plane B detector, when applicable.
- Lastly the correction factor (multiplication value) can be filled in.



8.3.3. Bodypart mapping

The DICOM standard bodypart is used in DOSE for effective dose calculation. In some cases, the device reports a body part that is not a DICOM standard, e.g. 'Abd' instead of 'Abdomen'. Non-DICOM standard bodyparts can be mapped to standard bodyparts in DOSE. In cases when bodyparts are not reported by the device, bodyparts can be mapped to the acquisition protocol, series or study description, and combinations of study description-series description and study description-acquisition protocol.

Certain cases were determined where series data is not available:

- No series: HL7, MPPS, Dose Report, Exam protocol
- No data at series level: HL7 + DICOM, MPPS + DICOM

In these cases, an option to calculate effective dose based on the Study Description or Protocol Name on study level is possible.


Series level is the predefined and preferred method, as different anatomical regions may be present in each series, and other derived parameters (e.g. organ doses) can only be calculated at series level.

8.3.3.1. MAPPING AT SERIES LEVEL

The following options are available for the bodypart mapping at series level:

- **Standardized bodypart:** if the reported DICOM Header body part is not an official DICOM body part then it can be converted to an official bodypart based on the mapping between the DICOM Header reported bodypart and the official body part.
- **Acquisition Protocol:** the body part will be based on the reported Acquisition Protocol of the series.
- **Series Description:** the body part will be based on the reported Series Description of the series.
- **Study Description – Series Description:** the body part will be based on the reported combination of Study Description-Series Description of the series.
- **Study Description – Acquisition Protocol:** the body part will be based on the reported combination of Study Description-Acquisition protocol of the series.


The user can map the correct body part for every device in *Configuration/Bodypart mapping and effective dose/series* to the original bodypart, acquisition protocol, series description, combination of Study description and Series Description or the combination of Study Description and Acquisition Protocol.




Bodypart mapping at series level


The configuration layout has a filter at the top that allows the user to filter on the original value (the value being mapped to a bodypart), the rules mapping to a certain bodypart and the link state (entries that are or are not mapped).




-The  button applies the filter.

-The  button resets and removes the applied filter.

The table below the filter shows the current mappings and for most mappings also the (auto-generated) not yet mapped value options. By pressing the edit (pencil) button at the top right, you can enter edit mode to alter the mappings.




Now mappings can be changed by altering or removing the dicom bodypart description from the dropdown. After editing, changes can be saved or cancelled from the menu on the top right. Saving or cancelling exists the editing mode.



Add, mass map, delete, save, cancel changes


The "add new mapping" button opens a new window where a new mapping for an original value not currently already in the table.




After filling in the fields, a check for duplicate values is performed and depending on the result, it informs the user or adds the new entry to the table.

In order to simplify and speed up the user's procedure, it is possible to perform a multiple selection, similar to the mapping of study groups in the option "Study group combinations".

To access the body part mass mapping, click on **Edit (pencil)** and then on the **"Window"** button.




Mass map button



Mass map window

Here you can select several entries at the same time, with help from the Filter Bar that can be found on the lower-right corner. After the selection has been done, you can select a "Dicom bodypart" from the drop-down menu on the right side of the window and click on the "Map selected" button to save it.

There's also the possibility to add this mapping to the already existing studies in DOSE and recalculate the effective dose. To do so, the user must click on the "Actions" button and select "Apply settings to existing studies".




Apply settings to existing studies

8.3.3.2. MAPPING AT STUDY LEVEL

The same process as described above can be followed to map bodyparts at study level, based on either the:

- Study description
- Protocol name

In **Device Overview** → **Configuration** → **Bodypart mapping and effective dose** → **Study**, the user can perform the mapping at study level.



Bodypart mapping at study level

NOTE: We recommend that this mapping is done only in the cases where series data information is missing, since effective dose calculation will be more accurate at series level.


8.3.3.3. DEVICE CONFIGURATION

The default methodology for determining the bodypart in case the configuration has not been altered is as follows:

1. Attempt to map the bodypart to dicom bodypart at series level, as described above .
2. If no matching bodypart is found and the bodypart is actually filled in dicom header or RDSR, it's taken from there.

3. If the bodypart is empty or does not exist, it's filled in "UNKNOWN".

This methodology can be altered by creating a different configuration in **Configuration/Bodypart mapping/Device Configuration**.




For each device, the user can select the *level* for effective dose calculation via a dropdown menu:

- Series (predefined)
- Study



Based on the choice above , it is possible to choose the mapping in case the bodypart is not available in the data source.




Here a primary option can be selected and this method will be attempted first. Only if through mapping no dicom bodypart can be determined, the program will attempt the secondary option, if one was selected. If, in the end, no dicom bodypart could be found, steps 2 and 3 from above will be used.

8.4. Device dashboards

Data panels consist of different dashboards, themselves consisting of different components showing specific data for the current device in graphs, charts and tables in a very specific date range.

In the *Dosimetry Overview*, the user with “Dashboard management” rights can add other dashboards by selecting the voice “Add dashboard” located in the “Actions” menu on the right upper corner of the Data panel, or by clicking the plus-sign on the tab sheet.



Add dashboard (2 ways)

When adding a new dashboard, the user will have to select the number and the type of the components. The dashboard order determines the place the dashboard will have in the data panel.

In the first step (general dashboard info), the user is asked to fill in the dashboard’s name, the description and the dashboard order.

Dashboard editor

!WARNING! This dashboard is shared with different devices. Any changes will also occur on the other dashboards.

1. General dashboard info 2. Permission for access 3. Dashboard layout 4. Data components

Dashboard name
Statistics

Dashboard description

Dashboard order
1

The order of a dashboard determines the **index** of the current dashboard in the horizontal list of all dashboards.

Dashboard editor (step 1)

In the second page of the dialog, the dashboard can be set to be:

- Private: only the user can see it
- Group of users: the user can define which users can see it
- Group of teams: the user can define which users can see it, based on the teams they belong to
- Public (predefined): all users can see it.

Dashboard editor

!WARNING! This dashboard is shared with different devices. Any changes will also occur on the other dashboards.

1. General dashboard info 2. Permission for access 3. Dashboard layout 4. Data components

Devices and sites: [Filter](#)

Users:

PUBLIC

PUBLIC


PRIVATE

USERS


TEAMS


Dashboard editor (step 2)

To add some users/teams, the user just needs to add them to the list.




Adding a dashboard to specific users

It is also possible to select the devices and/or sites where the dashboard will be shown by clicking on .




Filter with the devices/sites where the dashboard will be shown

In the next step (dashboard layout) the user needs to select the layout of the dashboard. The dashboard layouts are grouped for each amount of data components per dashboard. With a maximum of 6 items.



Dashboard editor (step 3)


The fourth and last step (data components) of editor wizard contains the data components that can be added to the dashboard. The user can choose the position of the component by selecting the desired place on the dashboard preview at the right of the window, and hereafter selecting the component.



Dashboard editor (step 4)


After clicking on *Finish*, a new dashboard is created and will be available for all the devices or/and sites and/or users/teams (when chosen in step 2).

To edit any of the above mentioned aspects of an existing dashboard, the user must select the dashboard and click on *Actions/Edit current dashboard*.




Edit a dashboard

When editing a dashboard, the user can choose also to remove it.



Remove a dashboard

In that case the system will ask if the user would like to remove the dashboard only for that device or for all devices among which the dashboard is shared.




Confirmation for removing a dashboard

For more information, refer to the videos *How to create dashboards*, and *How to delete dashboards* in our online training center.


8.4.1. Import of dashboards

The user has the possibility to import dashboards already created in another device. First, the dashboard must be exported, in **Actions**, in the upper right corner on **Device** level, by clicking on **Export current dashboard**. A download link will appear, and the user can save it by right clicking on **Download** and **Save link as**.



Save a dashboard

In the device chosen to import the dashboard, the user can import it by clicking on **Import dashboard** in **Actions** button. Users can upload the dashboard and choose between 5 import options. The 2 last options enable the user to import the dashboard to one modality type.




Dashboards import option

Warning: imported dashboards are public by default, and visible to all DOSE users with 'Device level' access. Access to the dashboard(s) can be restricted after the import is complete by clicking on **Actions/Edit current dashboard**.

8.4.2. Export


The export button can be found in **Actions** in the upper right corner on **Device** level. The user has the possibility to export data from a certain device on study/series level by clicking on **Export data**. The parameters that will be exported on study and/or series level can be selected by the user. The export scheme can be saved, and preloaded for the next export to save time.

For more information, refer to the video *How to export data on Device Level* in our online training center.



Export data

The user can select the study descriptions or study groups (created in Settings/Study Groups e.g. for compliance purposes) in order to export the specific fields for these studies.




Export data based on study descriptions

In the field "Configuration" the user can select the parameters for export based on:


- Device fields
- Study fields
- Series fields

Please note that when exporting PATIENT_AGE, PATIENT_AGE_TYPE should also be selected to indicate if the age is in years, months, etc.



Field configuration of the data to be exported

Any filter applied to the dashboards also applies to the exported studies.



Filter configuration

The main parameter units will be taken into account for exporting (e.g. DAP in dGy·cm²).

The data is exportable in CSV or XLSX format, and easy to generate. The export schemes can be personalized and saved under a specific name for quicker exporting in the future.




Saving an export scheme for future use








For more information, refer to the video *How to export data on Device Level* in our online training center.




8.4.3. Advanced filter

The advanced filter (found in the upper-right corner) allows the users to go more in depth. Filter configurations can also be saved for future use, using the save and load functions, and can be accessed from other devices. The filters are user-specific.



Filter button

User action	Application	Print screen
User selects a device and clicks on <i>Dosimetry overview</i>	The system loads all the dashboards created for the device.	
User clicks on the Filter button	The system shows Edit filter and Reset filter dropdown box; when "Edit filter" is selected the user can load a <i>Saved filter</i> , change the <i>Parameter</i> , or add any number of filters using the <i>Filter On</i> dropdown menu.	
User clicks on the Load button	The system loads the saved filter configuration, showing a dropdown box and minus button for each filter (Age, Gender, etc.) in the selected filter configuration.	
User clicks on "minus" button	The system removes filter.	
User clicks on Apply filter button	The system applies the filter only for the current session.	
User clicks on the Save button	The system saves the filter configuration.	
User clicks on the Save button next to <i>Saved filter</i>	The <i>Save filter configuration</i> window opens, allowing the user to assign a name (using letters, numbers, underscore) and select a <i>Scope</i> (Global or Modality) for the filter configuration. If the file name already exists, an error message will appear. The Action will always be <i>Save as a new filter</i> for	

	<p>new filter configurations, or <i>Overwrite</i> for existing filters, allowing the user to edit the configuration.</p>	
<p>User clicks on Action/Set filter as my predefined</p>	<p>The system shows all the previously saved filter configurations.</p>	
<p>User clicks on one of the available filters</p>	<p>The system sets the chosen filter as <i>predefined</i> for current user.</p>	
<p>User clicks on Action/Load from History 10 most used filters</p>	<p>The system shows the 10 most used filters.</p>	
<p>User clicks on one of the 10 most used filters</p>	<p>The system loads it into the filter configuration.</p>	
<p>User clicks on Action/Load from History 10 last used filters</p>	<p>The system shows the 10 last used filters.</p>	
<p>User click on one of the 10 last used filters</p>	<p>The system loads it into the filter configuration.</p>	


The following filter parameters are available for all modalities:

Type of variable	Name
Events	Event count.
Patient characteristics	Body Mass Index, Body Mass Index Group, Patient Age (D, W, M, Y), Age Group (Adults, various Pediatric groups), Patient Gender (M,F,O), Patient ID, Patient Weight range (kg), Patient Weight Group (kg).
Study	Accession Number, BodyPart Examined, Hide All Zero Dose Studies, Protocol Name, Show Ignored Studies, Study Composition, Study Description, Study Group, Study ID, StudyUID.
Responsible	Contributing Physician Name, Operator Name, Performing Physician Name, Reading Physician Name, Requesting Physician, Requesting Service, Responsible Physician Name, Ward.
Software	Software version.

For more information, refer to the videos *How to filter data and change parameters* and *How to omit studies from the statistics* in our online training center.

8.4.4. Manual entry of an XA study

Users are able to manually add an interventional procedure for an XA device. A new study can be added by clicking on *Actions* on *Dosimetry Overview*, and selecting *Add study*.



Add a new XA study manually

The following fields are mandatory:

- Study date and time
- Study UID
- Patient ID
- Patient date of birth
- Bodypart
- Protocol name or study description
- At least one dose indicator

The feature is available for all users whose role contains the *Data management* functionality.

8.4.5. Dosimetry overview

This data panel shows to the user many analyses based on different parameters.



Main parameter in Dosimetry Overview

The user can choose the main variable by changing the parameter in the advanced filter mentioned before. The available parameters include the following list:

Type of variable	Name	Modality
Dose	DLP, max CTDIvol	CT
	Weighted CTDIvol	CT, CBCT
	DAP	CBCT, CR, DR, DX, RF, XA
	DAP (acquisition or fluoroscopy), dose at Reference point (total, acquisition or fluoroscopy), Entrance dose, Cumulated REX	CR, DR, DX, RF, XA

	Peak Skin Dose (PSD)	XA
	Organ dose, Glandular dose, MGD	MG
	Effective dose ICRP 103	CT, CBCT, CR, DR, DX, RF, XA, MG
	Effective dose ICRP 60	CBCT*
	Effective dose ICRP 128	NM
Radioactivity	Administered, Administered per kg, Prepared	NM
Absorption Rate	SAR	MR
Time	Fluoro time	RF, XA
	Total Exposure Time	CBCT
Scan range	Scanlength	CT
Events	Event count, , Number of sequences	CR, DR, MG, XA
	Number of series	CT
	Number of sequences	MR

**Effective Dose ICRP60 is also calculated in DOSE for all modalities and can be exported, but it's not included as main parameter in the filter.*

The following table shows the units available for each parameter.

Parameters	Units
DLP	mGy.cm
max CTDIvol, Weighted CTDIvol, Entrance dose, Organ dose, Glandular dose, MGD	mGy
DAP	dGy.cm ² , cGy.cm ² , mGy.cm ² , Gy.cm ²
DAP from acquisitions, DAP from fluoroscopy	dGy.cm ²
Dose at Reference point (total, acquisition or fluoroscopy), Peak Skin Dose	Gy
Effective dose ICRP 103, ICRP 60, ICRP 128	mSv
Radioactivity (administered and prepared)	MBq
Radioactivity per kg (administered)	MBq/kg
SAR	W/kg
Fluoro time	s
Total Exposure time	ms
Scanlength	mm


The data panel consists in many different dashboards, which users can edit or remove as they want. They can even create a new dashboard, as explained above in the first part of this chapter.

These are some examples out of the possible dashboards:

Name	Description	Modality
Trends	Daily, Monthly, Weekday and hourly trend charts	All modalities
Workflow	Workload analysis and Moment evaluation: Mean dose evaluation per weekday and hour and the usage of a device	All modalities
Statistical overview	Minimum, maximum, mean and percentiles per study description, age groups	All modalities
Workload	Hourly average volume chart and hourly dose spider chart: Charts showing the device actual usage	All modalities
Protocol distribution	BodyPart/Protocol distribution and BodyPart/Protocol dose distribution	CR, DR, DX, RF, XA
Analysis	Average dose per thickness and view, Summary of most used settings and dose/thickness scatter chart: Analysis based on specific variables for MG (e.g. thickness, type of anode)	MG
Settings	Breast thickness histogram, compression force histogram and image view distribution: Charts based on specific variables for MG (e.g. thickness, compression)	MG
Reference	Dose-compression chart, with reference to EU and BEDRL	MG
Volumetric breast density	Charts and analysis based on breast density values (applicable only when a third party software like Volpara® Density™ is available)	MG
Patient information	Patient age/gender distribution	CR, DR, DX, RF, XA
Responsible	Mean parameter per requesting physician & operator, volume per requesting physician & operator: Charts based on human variables	All modalities
Tools	See next subparagraph	All modalities

Note that statistics are shown only when the asked kind of data is available.

8.4.6. Tools dashboard




This dashboard contains many features which can be used by users:

Name	Description
Study overview	This is the list of all the examinations of the selected device. By clicking the Actions button, situated on the lower right part of the dashboard, a user with permissions can edit, remove or disable examinations. The same menu allows to set body part or to recalculate effective doses, too.
Alerts overview and Configuration	This is the list of all the alerts for the device
Study comparison	This feature helps the customer comparing different patient's populations
Dynamic reports	The reporting feature enables a user to create reports that contain the data components as they are shown in the device its dashboards
Configuration	Configuration of dynamic limits and linking of exciting study groups.

8.4.6.1. ALERTS

It is possible to create alerts for examinations which dose is outside acceptable range. A user can create dynamic acceptable and achievable ranges in the Tools → Configuration dashboard: here the user can choose the lower and higher percentiles for each group and the date ranges. By clicking on the Action button (upper right part of the page) and then Calculate, the system will automatically mark as alerts all the outliers.

The user will find a clickable alert in the Study details page. There is a differentiation between the static limit on certain parameter (study groups) and the dynamic limits (based on you own configuration). See example:




Dynamic and static notifications

The red bar under the thumbnail corresponds to the static limit (with the study group mentioned under “Category”) and the notifications below correspond to both dynamic and static limits.

If there is at least one acquisition protocol/series description red/orange, the color on the study details will be red/orange.


The list of all the alerts can be found in the Tools → Alerts overview dashboard.



STUDY DATE	TYPE	STUDY	MESSAGE	VALUE	DIFFERENCE	RANGE
2014-01-29 09:44:54	DAPTOTAL	RX SCHOUDEUR RECHTS - 60 Y M	Outside acceptable range	60.65	-88.82 %	[0.00 ; 32.12]
2013-12-11 09:26:31	DAPTOTAL	RX THORAX FACEPROFEL - 41 Y F	Outside acceptable range	38.85	-64.97 %	[3.09 ; 23.55]
2014-01-30 11:58:35	DAPTOTAL	RX THORAX FACEPROFEL - 46 Y M	Outside acceptable range	32.62	-38.51 %	[3.09 ; 23.55]
2014-01-23 10:32:49	DAPTOTAL	RX THORAX FACEPROFEL - 64 Y F	Outside acceptable range	20.45	-25.05 %	[3.09 ; 23.55]
2013-12-19 16:09:10	DAPTOTAL	RX THORAX FACEPROFEL - 55 Y F	Outside acceptable range	28.36	-20.42 %	[3.09 ; 23.55]
2014-02-25 17:29:08	DAPTOTAL	RX THORAX FACEPROFEL - 59 Y M	Outside acceptable range	28.13	-19.45 %	[3.09 ; 23.55]
2014-02-12 17:13:40	DAPTOTAL	RX SCHOUDEUR LINKS - 60 Y F	Outside acceptable range	55.19	-11.77 %	[3.10 ; 49.38]
2014-02-11 17:19:08	DAPTOTAL	RX THORAX FACEPROFEL - 48 Y F	Outside acceptable range	2.73	11.65 %	[3.09 ; 23.55]
2014-01-03 10:42:30	DAPTOTAL	RX THORAX FACEPROFEL - 48 Y M	Outside acceptable range	24.88	-5.65 %	[3.09 ; 23.55]
2014-02-20 09:15:08	DAPTOTAL	RX THORAX FACEPROFEL - 63 Y M	Outside acceptable range	2.92	5.50 %	[3.09 ; 23.55]
2013-12-24 12:15:12	DAPTOTAL	RX THORAX FACEPROFEL - 59 Y F	Outside acceptable range	41.83	3.56 %	[3.09 ; 23.55]
2013-12-07 09:06:43	DAPTOTAL	RX KNIJE LINKS - 57 Y F	Outside acceptable range	41.83	-2.73 %	[1.57 ; 40.72]
2014-01-03 15:20:27	DAPTOTAL	RX THORAX FACEPROFEL - 23 Y M	Outside acceptable range	3.01	2.59 %	[3.09 ; 23.55]
2014-02-05 12:21:04	DAPTOTAL	RX LUMBALE WERVELKOLOM EN BEKKEN - 26 Y F	Outside acceptable range	36.96	2.40 %	[37.87 ; 309.80]
2014-01-08 16:37:24	DAPTOTAL	RX BEKKEN - 74 Y M	Outside acceptable range	169.57	-2.27 %	[1.77 ; 165.80]
2014-01-23 14:43:07	DAPTOTAL	RX THORAX FACEPROFEL - 28 Y F	Outside acceptable range	3.04	1.62 %	[3.09 ; 23.55]
2014-02-19 10:22:15	DAPTOTAL	RX CERVICALE WERVELKOLOM - 58 Y M	Outside acceptable range	33.09	-1.41 %	[0.00 ; 32.63]
2013-12-26 11:43:30	DAPTOTAL	RX THORAX FACEPROFEL - 55 Y F	Outside acceptable range	23.84	-1.23 %	[3.09 ; 23.55]

Alerts overview


The same alerts list can be also found by clicking on the orange button with the number of alerts. If the user is only interested in Static or Dynamic alerts, there is a possibility to filter on this parameter via Action button on Device level.



8773 of 8773 notifications visible (All alerts) Hide filter Filter Actions

Parameter type: DLP | mGy.cm

Orange button (all alerts)



5 of 5202 studies visible **8773** of 8773 notifications visible (All alerts)


Legal Compliance Tools +

- Show Static Alerts
- Show Dynamic Alerts
- Show alerts with cause
- Show alerts with action
- Show alerts with justification

- Add dashboard
- Edit current dashboard
- Create draft from dashboards
- Export current dashboard
- Reset filter
- Reset to Predefined Filter
- Alerts Filter
- Save Alert Filter
- Export data
- Refresh study data
- Calculate

Choose static or dynamic alerts

You can also see the alerts trend. This graph is clickable but note that there all 2 types of alert mixed (static and dynamic).



Qaelum Expert mode Welcome qaelum Profile Logout Switch lights

Start Portal Study Patient Device Modality Location Quality Management System Contrast Settings

Filter Add Site

DEMO 90
Demo Hospital 1
Room 01 - Inspiration
Room 02 - Siemens Somatom Sensatio
Room 03 - FujiFilm DEVO
Room 04 - Toshiba Aquilion 16
Room 05 - Siemens Somatom Definitic
Room 06 - Siemens Axiom Artis
Room 07 - Siemens Fluorospot Compa
Room 08 - NM Hotlab
Room 09 - NIM NTSG
Demo Hospital 2
Room 01 - SIEMENS SOMATOM Defir

Room 01 - Inspiration | MG **325** of 326 studies visible **21** of 21 notifications visible (All alerts) Hide filter Filter Actions

Patient radiation dose monitoring Parameter type: OrganDose | mGy

Compliance per thickness x entrance x Scatter x operators x Glandularity x Operators x Informe SMS x t x compliance series x Alerts x Tools +

ALERT TREND PER DAY

ALERT TREND PER WEEK

Alerts

Error Warning

Date range (last 10 years) : 2008/09/11 - 2018/09/12

© 2012-2018 QAELUM NV. All rights reserved | About


You can also configure the graphs including the alerts trends into a dynamic report (see *Dynamic Reports*).

8.4.6.2. STUDY COMPARISON

This feature enables the user to compare groups of studies defined by some kind of population (e.g. age range, sex, body part, examination date...). This comparison can be shown with different tables and charts.


8.4.6.2.1. ADD A POPULATION

The user can select a population by using a preset or by creating a new one: to add a population the user must click on the Action button.




Study comparison

A population is defined with a name and a description.




Select population


The user can create his own population by setting up specified filters.



For example, you can select the "Study Date" filter, setting it to "Before or equal to" the desired date for the first population, you can choose "After" for the second population so you can compare 2 populations before and after an optimization.




You can also add additional filters; all the possibilities are listed on the screenshot below. For example, if you select "Body Part Examined" and choose ABDOMEN you will have a population that includes all the protocols linked to that body part.



Choose filter

Once you have all the filters you want for the population, you can Add Population.




Population filter

You can follow the same steps to create as many populations as you need.

NOTE: Be careful to click always on "Add population" every time you Edit a population, if you do not do this, the population will be removed.

8.4.6.2.2. ADD A COMPARISON

In the same way, the user can choose a tab or a chart to show the comparison. A new comparison can be added by clicking on the Actions button on the Comparisons side.



You must choose from a list of predefined comparisons.



8.4.6.2.3. COMPARE POPULATIONS

After having chosen at least two populations and a comparison, graphs and tables are shown when the user clicks the Compare button.

Room 05 - Siemens Somatom Definition Flash | CT
Patient radiation dose monitoring

3162 of 3162 studies visible | 598 of 598 notifications visible (All alerts) | Hide filter | Filter | Actions

Parameter type: DLP | mGy.cm

workflow tes Patient information Testing sajxc cldi compliance study Trends An_test moment evaluation Tools +

Study Overview Alerts overview Study comparison Dynamic reports Configuration

Populations


Populations Actions Comparisons Actions

CT Abdomen before
CT Abdomen after


Dose histogram comparison

Load combination Save combination Clear lists Compare

Clink on "Compare" button




Dose histogram




Dose vs Time scatter plot

8.4.6.2.4. SAVE AND LOAD COMBINATIONS

A combination of populations and comparisons can be saved by clicking on *Save combination* and indicating a name.



The populations can be preloaded using *Load combination* when needed.




For more information, refer to the video *How to use study comparison on Device Level* in our online training center.

8.4.6.3. DYNAMIC LIMITS

Dynamic limits can be configured based on the lower/upper acceptable (ACC) and lower/upper achievable (ACH) limits using device data collected by DOSE. Percentiles for the calculation of the limits are configurable by users with the *Alerts Management* functionality. Limits can be based on the study description, series combination (acq. protocols) or protocol name. When limit configurations are modified click on *Calculate* for the changes to go into effect.

Dynamic limits are automatically calculated using the data of the selected studies of the local device. They can be recalculated when necessary to allow for comparison with the current local practice. For configuration, follow the steps below:

- Select a device → **Dosimetry overview** → **Tools** → **Configuration** → **Dynamic limits**
- The following must be entered:
 - o Lower and upper ACH and ACC limits (in percentages).
 - o Date range.
 - o Select whether the limits should be based on the: Study description, series combination or protocol name. In the case of MG, limits can only be recalculated for a **series combination**, please verify this option has been selected. In the future, it will be possible to also choose 'series description' and 'acquisition protocol.'
 - o Minimum study count: minimum amount of a series combination in order to calculate the limit.



The screenshot shows the 'Dynamic Limits' configuration interface. Key elements include:


- Room 02 - Siemens Somatom Sensation 64 | CT
- 5194 of 5202 studies visible
- 3427 of 3427 notifications visible (All alerts)
- Parameter type: DLP | mGy.cm
- Navigation tabs: Statistics, Patient information, Legal compliance, Responsibles, Age analysis, Acquisitions, BoxPlots, test, test, glg, Comments, Tools, +
- Sub-navigation: Reference dashboard overview, Study Overview, Alerts overview, Study comparison, Dynamic reports, Static reports, Batch actions, Configuration
- Dynamic Limits section:
 - Lower ACC: 0% | Lower ACH: 0%
 - Upper ACH: 95% | Upper ACC: 98%
 - Start date: 11-09-2008 | End date: 31-12-2016
 - Minimum study count: 5
 - Dose types: All
 - Limits based on: Study description
 - Buttons: Calculate, Expand / Collapse
- Table of calculated limits:

STUDY DESCRIPTION	STUDY COUNT	LOWER ACC	LOWER ACH	UPPER ACH	UPPER ACC
▼ RAD ct abd 00.[1-2-2-17]					
17 - 150	20	0.0	0.0	1330.0	1340.0
RAD ct abd 01.11-2-2-1881					

Set dynamic limits

- Click on **Calculate**
- In the table below, the calculated dynamic limits per study description/series combination/protocol name will be displayed.

From now on, when a new study arrives, the system will check if its dose is above the dynamic levels and create a notification accordingly. If notifications of past studies need to be recalculated, the user must click on **Recalculate notifications** inside the **Actions** button of the “Study groups” tab (located next to the ‘Dynamic Limits’ tab).



Recalculate notifications


To remove the dynamic limits, all the lower/upper ACH/ACC fields must be to zero and recalculated by clicking on **Calculate**. To remove the associated notifications, click on **Recalculate notifications**.

For more information, refer to the video *How to configure dynamic limits & alerts* in our online training center.

8.4.7. Export/Save data

The system allows users to export or save dashboards’ data and charts in different

Export data as	Save a chart as
XLS	PNG
PDF	JPEG
CSV	SVG




formats.

Export and save dashboard

The user can save a graph or export data by clicking the machinery button which is located in the upper right part of every data component.

8.4.8. Favorite dashboard

To mark a chart or a table as favorite, the user can use the star button in the upper part of every element.



Set graph as favorite

The starred data components will be copied automatically to the personal Start portal: the user can decide to show the same graph or table in the Overview page by right clicking on it.



Add favorite graph to the Overview page (Start Portal)

8.4.9. Compliance monitoring

As many countries have different regulations about Diagnostic Reference Limits (DRLs) and the way to report each device's performances, it is important to choose the legal limits properly. In DOSE, compliance is done with the help of "Study Groups" that have some "Limits" (DRLs) associated to them and to which several kinds of studies can be mapped.

8.4.9.1. COMPLIANCE CONFIGURATION

Compliance configuration consists of 4 steps:

1. First the exam type categories ("Study groups") must be defined, where different studies are combined based on different parameters.
2. Then the DRL version is defined, indicating which limit to apply depending on the date range.
3. After that, "Limits" are defined, taking into account the version name.
4. Finally, the studies are mapped to the Study Groups on Device level.


8.4.9.1.1. STUDY GROUPS

The study groups are defined in Settings → Study Groups.

There are 4 sections: Study Groups CT, Study Groups General (RX, IR, CBCT, NM), Study Groups MG, CT SSDE Study Groups.

Each category has the following properties:

- Group name
- Modality (CT / MG / RX / IR / NM / CBCT)
- Start age (min age of the category)
- End age (max age of the category)
- Age type (D / Y) (Day/Year)
- Gender (M / F / MF / MFO / O)
- Start weight (CT, General) / Min Breast thickness (MG)
- End weight (CT, General)/ Max Breast thickness (MG)
- Min event count (min number of series)
- Max event count (max number of series)
- Minimum BMI (min Body Mass Index)
- Maximum BMI (max Body Mass Index)
- Minimum Scan Length (CT)
- Maximum Scan Length (CT)
- Minimum spiral count (CT)
- Maximum spiral count (CT)
- Minimum axial count (CT)
- Maximum axial count (CT)
- Start Eff. Diameter (SSDE)
- End Eff. Diameter (SSDE)
- Start WED (SSDE)
- End WED (SSDE)
- Start LAT Thickness (SSDE)
- End LAT Thickness (SSDE)
- Purpose (Legal/Customized)



Study Groups CT

Name	Modality	Start age	End age	Age type	Gender	Start Weight	End Weight	Min Event Count	Max Event Count	Minimum BMI	Maximum BMI	Minimum Scan Length	Maximum Scan Length	Purpose	Min Spiral Count	Max Spiral Count	Min Axial Count	Max Axial Count
Schelde (hermen) (5-10 jaar)	CT	5	10	Y	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1
Schelde (hermen) (<1 jaar)	CT	0	365	D	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1
Sinissen	CT	15	150	Y	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1
Sinissen (1-5 jaar)	CT	1	4	Y	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1
Sinissen (10-15 jaar)	CT	10	14	Y	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1
Sinissen (5-10 jaar)	CT	5	9	Y	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	0	1	0	1
Test nls series	CT	0	150	Y	MFO	0	250	0	999	-1	-1	-1	-1	CUSTOMIZED	-1	-1	0	2
Thorax	CT	15	150	Y	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1
Thorax (1-5 jaar)	CT	1	4	Y	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1
Thorax (10-15 jaar)	CT	10	14	Y	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1
Thorax (5-10 jaar)	CT	5	9	Y	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1
Thorax (<1 jaar)	CT	0	365	D	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1
Thorax-abdomen - volledig	CT	15	150	Y	MFO	0	99	0	999	-1	-1	-1	-1	LEGAL	-1	-1	-1	-1

Study Groups General


Name	Modality	Start age	End age	Age type	Gender	Start Weight	End Weight	Min Event Count	Max Event Count	Minimum BMI	Maximum BMI	Purpose
Bedden a.p.p.a. (1 bis < 5 Jahre)	RX	1	4	Y	MFO	0	999	0	999	-1	-1	LEGAL
Bedden a.p.p.a. (10 bis < 15 Jahre)	RX	10	14	Y	MFO	0	999	0	999	-1	-1	LEGAL
Bedden a.p.p.a. (5 bis < 10 Jahre)	RX	5	9	Y	MFO	0	999	0	999	-1	-1	LEGAL
Breusteekbuis a.p.p.a.	RX	15	150	Y	MFO	0	999	0	999	-1	-1	LEGAL
Breusteekbuis lat	RX	15	150	Y	MFO	0	999	0	999	-1	-1	LEGAL
Class (APPA) by (CSCS)	RX	0	5	Y	MFO	0	99	0	999	-1	-1	LEGAL
Class (APPA)	RX	0	150	Y	MFO	-1	-1	-1	-1	-1	-1	CUSTOMIZED
Dental (DVT FOW-Scan2)	RX	15	150	Y	MFO	0	999	0	999	-1	-1	LEGAL
Dental (DVT FOW-Scan2)	RX	15	150	Y	MFO	0	999	0	999	-1	-1	LEGAL
EPU - Onkariemysystem	IR	15	150	Y	MFO	0	999	0	999	-1	-1	LEGAL
EPU - Einlammsystem	IR	15	150	Y	MFO	0	999	0	999	-1	-1	LEGAL
EPU - Ziektemysystem	IR	15	150	Y	MFO	0	999	0	999	-1	-1	LEGAL
ERCP (intervention)	RX	15	150	Y	MFO	0	999	0	999	-1	-1	LEGAL

Configure study groups

Note: when a parameter is set to -1 it is not taken into account in defining the study group. By setting the value to an integer value above zero, the field is taken into account. For example for CT study groups, setting Min Spiral Count=2 and Max Spiral Count=4 means that only studies which contain either 2, 3 or 4 spiral series are included in this study group. Setting Min Spiral Count=-1 and Max Spiral Count=-1 means that studies are included in this study group independently from how many spiral series are present.

A new study group can be added in **Settings → Study Groups → Overview**.

- Select the modality tab and click on **Edit – Add**.
- An additional empty line is added below.
- Fill in all the required fields (name, modality, age etc.)
- The purpose can be LEGAL or CUSTOMIZED.
- Click on **Save**.



Study Groups Overview

NAME	MODALITY	START AGE	END AGE	AGE TYPE	GENDER	START WEIGHT	END WEIGHT	MIN EVENT COUNT	MAX EVENT COUNT
Thorax-abdomen Volledig onderzoek NO DRL (5-10	CT			Y	MFO	0	99	0	999
Thorax-abdomen Volledig onderzoek NO DRL (<1 j	CT			D	MFO	0	99	0	999
Thorax-abdomen Volledig onderzoek NO DRL (10-1	CT			Y	MFO	0	99	0	999
				Y	MFO				

Buttons: Remove, Add, Cancel, Save


A study group can be edited in **Settings → Study Groups → Overview**.

- Select the modality tab and click on **Edit**
- Edit the value which requires an edit

- c. Click on **Save to commit the changes**.
- d. If the MG modality is selected, it will appear in the **Study Groups MG** list. In that case, click on Edit and modify **Max. Breast Thickness** to indicate the desired range.

A study group can be removed in **Settings → Study Groups → Overview**.

- e. Select the modality tab and click on **Edit**
- f. Click on the Study Group line and this will be highlighted
- g. Then click on **Remove**
- h. Click on **Save**.



Name	Modality	Start age	End age	Age type	Gender	Start Weight	End Weight	Min Event Count	Max Event Count	Minimum BMI	Maximum BMI	Minimum
2018 Thorax und Oberbauch	CT	18	120	Y	MFC	0	150	0	999	-1	-1	
ANGIO-TC ABDOMEN-PELVIS	CT	16	150	Y	MFC	0	99	0	999	-1	-1	
ANGIO-TC CRANIAL (F. OF Willis)	CT	16	150	Y	MFC	0	99	0	999	-1	-1	
ANGIO-TC TORACICA PULMONAR	CT	16	150	Y	MFC	0	99	0	999	-1	-1	
ANGIO-TC TRONCOS SUPRAORTICOS	CT	16	150	Y	MFC	0	99	0	999	-1	-1	
Abdomen (10 bis <15 Jahre)	CT	10	15	Y	MFC	0	150	0	999	-1	-1	
Abdomen (5 bis <10 Jahre)	CT	5	10	Y	MFC	0	150	0	999	-1	-1	
Abdomen (z. B. Leber, Pankreas)	CT	15	115	Y	MFC	0	150	0	999	-1	-1	

Deleting a Study Group which is actively used may result in missing notifications. For this reason, the deletion of Study Groups is regulated in DOSE.

When a Study Group gets marked for deletion:

- If the Study Group is not mapped to any device, then it is considered as inactive and it gets deleted
- If the Study Group is mapped to one or more devices, then it is considered as active and the user is asked to confirm the deletion. For each Study Group, the user is provided with the list of mapped devices and with a description of the consequences of the deletion

Confirmation


Are You Sure?

Are you sure you want to remove this Study Group?
The study group with name **TEST ACQUISITON TYPE** is mapped to parameters on the following devices:

- An's hospital - GE - BrightSpeed
- An's hospital - GE - BrightSpeed

By removing this study group, the notifications based on the specific limits won't be calculated anymore. The removal of the study group won't affect the calculation of notifications for any non-involved study group. The already existing notification won't be deleted immediately, but may disappear if any user will trigger a notification recalculation. By clicking on Yes, the study group will be deleted. By clicking on Discard, the study group won't be deleted. Note that the other edits and deletions, except for any other notices appearing, are going to be carried on after this message will be closed.

[YES](#) [Discard](#)



Qaelum Expert mode Welcome qaelum Profile Logout Switch lights

Start Portal Study Patient Device Modality Location Quality Management System Contrast Settings

Torax-abdomen CT 15 150 Y MFO 0 99 0 999 LEGAL U


Study Groups MG

NAME	MODALITY	START AGE	END AGE	AGE TYPE	GENDER	MIN BREAST THICKNESS	MAX BREAST THICKNESS	MIN EVENT COUNT	MAX EVENT COUNT	PURPOSE
Mammography	MG	18	115	Y	MFO	0.0	10.0	0	999	CUSTOMIZED

© 2012-2018 QAEELUM NV All rights reserved | About

8.4.9.1.2. VERSIONS


Diagnostic reference levels (DRLs) are supposed to change throughout time; therefore, studies should be evaluated based on the applicable and most up-to-date DRLs. For this reason we have in DOSE a way to have a track of the different versions of DRL's.



Name	Start date	End date
DEFAULT	01/01/2002	31/12/2040

DRL versions


It is possible to filter by version on the Limits table by using a dropdown menu.




Name	Scope	Lower ACC	Lower ACH	Upper ACH	Upper ACC	Purpose	Version
Steve	STUDY	0.0	0.0	5.0	10.0	LEGAL	DEFAULT
Angiography CT (Adult)	SERIES	0.0	0.0	5.0	8.0	LEGAL	DEFAULT
TestBMI	STUDY	0.0	0.0	5.0	8.0	LEGAL	DEFAULT
COMPLIANCE TEST	SERIES	0.0	0.0	0.0	200.0	LEGAL	DEFAULT
CT pediatric chest	SERIES	0.0	0.0	400.0	600.0	CUSTOMIZED	DEFAULT
CT pediatric chest	STUDY	0.0	0.0	0.0	0.0	CUSTOMIZED	DEFAULT
Test Niki benchmark series	SERIES	0.0	0.0	200.0	300.0	CUSTOMIZED	DEFAULT
Test Niki benchmark study	STUDY	0.0	0.0	5.0	8.0	CUSTOMIZED	DEFAULT
Ana	SERIES	100.0	2000.0	5000.0	10000.0	CUSTOMIZED	DEFAULT
Thorax Angio CT (Adult)	STUDY	0.0	0.0	220.0	460.0	LEGAL	DEFAULT
Abdomen (Adult)	STUDY	0.0	0.0	420.0	700.0	LEGAL	DEFAULT
Heart (CCTA) (Adult)	STUDY	0.0	0.0	250.0	530.0	LEGAL	DEFAULT
Colon (Adult)	STUDY	0.0	0.0	230.0	530.0	LEGAL	DEFAULT
Thorax-abdomen (Adult)	STUDY	0.0	0.0	550.0	950.0	LEGAL	DEFAULT

When calculating the notifications, the date range is therefore taken into account.

To add a new version, go to **Settings** → **Study Groups** → **Versions**. Here you can click on **“Add”** to write a version name and define the date range.



Note: date range of the versions cannot overlap; however, it is possible to have date ranges without any version applicable. When this is the case, a message is displayed in "Study Details" stating that the study is linked to a Study Group (Category) but it has not limits applicable:




8.4.9.1.3. LIMITS

After you have defined the different groups and versions, you can assign for each parameter type static limits to the groups. To do so:

- Go to **Limits**
- Choose the right modality (CT, RX, MG...) and parameter (CTDI, DAP, MGD...) for the limit configuration.
- Click on **Edit - Add**
- In the dropdown menu you can choose customized made study group and fill in the corresponding limits.
- It is possible to create a study group with the limits based on SERIES/STUDY. You can even specify SERIES16 (cm phantom) or SERIES32 (cm phantom) only for CT of course or in general all SERIES.


*NOTE: For Mammography, it only makes sense to create a study group with the limits scoped on **SERIES**. The other options appear (STUDY, SERIES16, SERIES32) but they are not suitable since DRLs for mammography are defined only on **SERIES** level.*



- f. For limits with scope SERIES, the Acquisition Type (Helical, Spiral...) can be defined. In order to use this feature, a list of acquisition types must be defined first. Refer to [Advanced](#) for more details.

A study group limit can be edited in **Settings → Study Groups → Overview**.


- Select the modality tab and click on **Edit**
- Edit the value which requires an edit
- Click on **Save to commit the changes**.



A study group limit can be removed in **Settings → Study Groups → Limits**.

- Select the modality tab and click on **Edit**

- b. Click on the Study Group Limit line and this will be highlighted
- c. Then click on Remove
- d. Click on **Save**.



Study Groups

Versions Overview Limits Advanced

CT MG RX IR NM CBCT MR US XA

DLP (DLP) Scan Length (SCANLENGTH) CTDI (CTDI) Weighted CTDIvol (WEIGHTED_CTDI_VOL) Effective Dose ICRP103 (CT_EFFECTIVE_DOSE_103) SSDE (SSDE)

DRL Versions DEFAULT2023

In order to avoid misconfigurations and potential data loss, it's not possible to switch between DRL versions while on Edit mode. To switch the version safely, please save your current edits first.


Name	Acquisition type	Scope	Lower ACC	Lower ACH	Upper ACH	Upper ACC	Purpose	Version
Abdomen - volledig	N/A	STUDY	0.0	0.0	420.0	570.0	LEGAL	DEFAULT2023
Angio CT thorax	N/A	STUDY	0.0	0.0	230.0	300.0	LEGAL	DEFAULT2023
Hart (CCTA) - calciumscore	N/A	STUDY	0.0	0.0	50.0	70.0	LEGAL	DEFAULT2023
Hart (CCTA) - angiografie	N/A	STUDY	0.0	0.0	230.0	300.0	LEGAL	DEFAULT2023
Abdomen (Adult)	N/A	STUDY	0.0	0.0	380.0	430.0	LEGAL	DEFAULT2023
Colon - acquisitie rug	N/A	STUDY	0.0	0.0	160.0	200.0	LEGAL	DEFAULT2023
Colon - acquisitie buik	N/A	STUDY	0.0	0.0	120.0	160.0	LEGAL	DEFAULT2023
Cervicale wervelzuil	N/A	STUDY	0.0	0.0	350.0	450.0	LEGAL	DEFAULT2023
Lumbale wervelzuil	N/A	STUDY	0.0	0.0	450.0	600.0	LEGAL	DEFAULT2023
Schedel (hersenen)	N/A	STUDY	0.0	0.0	720.0	800.0	LEGAL	DEFAULT2023
Sinussen	N/A	STUDY	0.0	0.0	50.0	65.0	LEGAL	DEFAULT2023
Thorax	N/A	STUDY	0.0	0.0	175.0	210.0	LEGAL	DEFAULT2023

Remove + Add Cancel Save

Deleting a Study Group Limit which is actively used may result in missing notifications. For this reason, the deletion of Study Group Limits is regulated in DOSE.

When a Study Group Limit gets marked for deletion:

- If the Study Group is not mapped to any device, then it is considered as inactive and it gets deleted
- If the Study Group is mapped to one or more devices, then it is considered as active and the user is asked to confirm the deletion. For each Study Group, the user is provided with the list of mapped devices and with a description of the consequences of the deletion




The following parameter types can be assigned for each modality:

CT	MG	RX	IR	NM	CBCT
DLP ScanLength CTDI _{max} wCTDI _{vol} EffDose ₁₀₃ SSDE	OrganDose MGD EffDose ₁₀₃	DAP EntranceDose FluoroscopyTime EffDose ₁₀₃	DAP EntranceDose FluoroscopyTime Event count EffDose 103 PSD Dose at ref point	Radioactivity (administered) Radioactivity (prepared) EffDose ICRP128 Radioactivity per kg (administered)	DAP EffDose ₁₀₃

If studies are outside the acceptable (ACC) limits this should be considered as critical for the specific category and parameter type. If the studies are outside the achievable (ACH) limits and still below the ACC limits, this should be considered as a warning. If only one upper limit is defined, we assign the same value to ACH and ACC and if no lower limit is

applicable, we assign the value 0. All the studies between the lower and upper ACH boundaries are considered to have a normal parameter value.



Lower/higher ACC/ACH

Study Groups							
Overview		Limits					
CT	MG	RX	IR	NM	CBCT	MR	US
DAP		Entrance Dose	Fluoroscopy time				
NAME	SCOPE	LOWER ACC	LOWER ACH	UPPER ACH	UPPER ACC	PURPOSE	
Abdomen Meenvoudig onderzoek (Adult)	STUDY	0.0	0.0	17.5	40.0	LEGAL	
Abdomen Enkelvoudig onderzoek (Adult)	STUDY	0.0	0.0	8.5	27.5	LEGAL	
Bekken face (AP) (Adult)	STUDY	0.0	0.0	11.0	35.0	LEGAL	
Thorax PA (Adult)	STUDY	0.0	0.0	1.0	3.0	LEGAL	
Thorax volledig (Adult)	STUDY	0.0	0.0	3.5	11.0	LEGAL	

Set limits for the study groups

8.4.9.1.4. ADVANCED

Constant angle acquisition notifications can be avoided by filtering by “Acquisition type”. This filter is present in the definition of limits since different acquisition types may have different DRLs.

Study Groups

Versions Overview **Limits** Advanced


CT MG RX IR NM CBCT MR US

DLP (DLP) Scan Length (SCANLENGTH) CTDI (CTDI) Weighted CTDIvol (WEIGHTED_CTDI_VOL)

DRL Versions **DEFAULT**

Name	Acquisition type	Scope
Test Dimitri 2	N/A	SERIES
Sinuses (Adult)	N/A	STUDY
Head (Adult)	N/A	STUDY
Lumbar Spine (Adult)	N/A	STUDY
Cervical Spine (Adult)	N/A	STUDY
Thorax-abdomen (Adult)	N/A	STUDY
Colon (Adult)	N/A	STUDY
Heart (CCTA) (Adult)	N/A	STUDY
Abdomen (Adult)	N/A	STUDY

In order to use this feature, a list of acquisition types must be defined first. Qaelum provides a list with the most typical cases, but the user may change this at will. To do so, go to **Settings/Study Groups/Advanced** and click on **Import**.



Indicate the list of acquisition types that should be included in some rules:



Import acquisition types

Template:Name;


```
Spiral Acquisition;
Constant Angle Acquisition;
Helical;
Localizer;
```

Discard Commit

Then, a list of certain values must be created by in **Acquisition types lists**, by clicking on **Actions/Add List**:



If the Limits are edited, the user will see the following dropdown menu:



If the "Import" option is used, an existing acquisition type list must be indicated.

Import Study Group Limits by CSV

TEMPLATE: name; **acquisition type list**; modality; parameter type; scope; lower ACC; lower ACH; upper ACH; upper ACC;

Lines with an incorrect number of arguments will be ignored!
Any previous study group limits will be deleted!

```

PSD;N/A;IR;PSD;STUDY;0.0;0.0;0.0;0.0;3.0;LEGAL;Version 2019
PSD;N/A;IR;DAP;STUDY;0.0;0.0;0.0;0.0;20.0;LEGAL;Version 2019
CT abdomen and pelvis without contrast;N/A;CT;DLP;SERIES;0.0;0.0;0.0;0.0;100.0;LEGAL;Version 2019
CT chest + C standard size;Spiral/Helical;CT;DLP;SERIES;0.0;0.0;0.0;0.0;10.0;LEGAL;Version 2019
Mammography- craniocaudal (CC) 45 mm breast thickness;N/A;MG;ORGANDOSE;SERIES;0.0;0.0;0.01;0.03;LEGAL;Version 2014-2019
Thorax (Adult);N/A;CT;DLP;STUDY;0.0;0.0;200.0;320.0;LEGAL;Version 2014-2019
COMPLIANCE TEST;N/A;CT;DLP;STUDY;0.0;0.0;0.0;200.0;LEGAL;Version 2014-2019
Becken AP/PA (1 bis < 5 Jahre);N/A;RX;DAP;STUDY;0.0;0.0;0.0;1.2;LEGAL;DEFAULT
Becken AP/PA (5 bis < 10 Jahre);N/A;RX;DAP;STUDY;0.0;0.0;0.0;2.5;LEGAL;DEFAULT
AP Abdomen (ACR);N/A;RX;DAP;STUDY;0.0;0.0;0.0;3.4;LEGAL;DEFAULT
AP Abdomen (HCSC35);N/A;RX;DAP;STUDY;0.0;0.0;7.0;15.0;LEGAL;DEFAULT
AP L-spine (ACR);N/A;RX;DAP;STUDY;0.0;0.0;0.0;4.2;LEGAL;DEFAULT
AP L-spine (HCSC35);N/A;RX;DAP;STUDY;0.0;0.0;7.0;10.0;LEGAL;DEFAULT
Chest (AP/PA) 5y (HCSC35);N/A;RX;DAP;STUDY;0.0;0.0;0.05;0.15;LEGAL;DEFAULT
PA Chest (ACR);N/A;RX;DAP;STUDY;0.0;0.0;0.0;0.15;LEGAL;DEFAULT
PA Chest (HCSC35);N/A;RX;DAP;STUDY;0.0;0.0;0.2;0.3;LEGAL;DEFAULT
Pediatric Chest (PA) (w/wo grid) (ACR);N/A;RX;DAP;STUDY;0.0;0.0;0.06;0.12;LEGAL;DEFAULT
Steve;N/A;CT;WEIGHTED_CTDI_VOL;STUDY;0.0;0.0;20.0;30.0;LEGAL;DEFAULT
AP Abdomen (ACR);N/A;RX;ENTRANCEDOSE;STUDY;0.0;0.0;0.0;3.4;LEGAL;DEFAULT
AP Abdomen (HCSC35);N/A;RX;ENTRANCEDOSE;STUDY;0.0;0.0;7.0;15.0;LEGAL;DEFAULT
AP L-spine (ACR);N/A;RX;ENTRANCEDOSE;STUDY;0.0;0.0;0.0;4.2;LEGAL;DEFAULT
AP L-spine (HCSC35);N/A;RX;ENTRANCEDOSE;STUDY;0.0;0.0;7.0;10.0;LEGAL;DEFAULT
Chest (AP/PA) 5y (HCSC35);N/A;RX;ENTRANCEDOSE;STUDY;0.0;0.0;0.05;0.15;LEGAL;DEFAULT
PA Chest (ACR);N/A;RX;ENTRANCEDOSE;STUDY;0.0;0.0;0.0;0.15;LEGAL;DEFAULT

```

Discard Commit

This way, if a limit is defined for “helical” acquisitions only, and there are two series with the same name, one helical and one topogram, only the helical acquisition will be taken into account.

8.4.9.1.5. COMPLIANCE MAPPING

Once the study group and its corresponding limits are created, studies from the device can be linked to the user defined study groups. Several connection options are possible:


- At **study** level (groups with limit scope STUDY or SERIES can be selected):
 - Study description
 - Series combination
 - Protocol name
 - Radiopharmaceutical (only for NM)
- At **series** level (groups with limit scope SERIES can be selected):
 - Acquisition protocol
 - Series description
 - Study Description – Series Description (combination)
 - Study Description – Acquisition Protocol (combination)

The following table presents a summary and explanation of these rules:



DRL SCOPE	STUDY	SERIES (all) SERIES32 (CT) SERIES 16 (CT)
Linked parameter		
STUDY DESCRIPTION PROTOCOL NAME SERIES COMBINATION RADIOPHARMACEUTICAL (NM)	✓	✓ All series of the study will be considered as belonging to the study group and checked against the series DRL
SERIES DESCRIPTION ACQUISITION PROTOCOL STUDY-SERIES DESCRIPTION STUDY DESC. – ACQ. PROTOCOL	X Not possible as it would lead to multiple repeated notifications	✓

To perform the mapping, access **Dosimetry Overview** (Device level) ⇒ **Tools** ⇒ **Configuration** ⇒ **Study Groups**.

The desired connection option from the list above can be selected from the dropdown menu.




For example, if “Acquisition protocol” is selected (only possible if the scope of the study group(s) is SERIES), click on **Edit** (5), select the desired study group, click on the **>>** (7) button to select the study group(s) and click on **Save** (8). Each acquisition protocol name can be linked to as many study groups as needed (e.g. adult and pediatric, encompassing different ranges of BMI, event count etc.), as long as the study group population is not exactly the same (as defined in Settings/Study Groups/Overview). The “Study group filter” (6) can be used to filter out specific study groups. Repeat the steps to map as many device acquisition protocols as needed.



NOTE: Only one connection option can be mapped for each study group, i.e. if a study group is mapped using the “acquisition protocol” connection option, the same group cannot also be mapped using the “series description” connection.


The data displayed to the user during study group mapping corresponds to the selected device and date range.

Alternatively the *Make study group combinations* tool in *Actions* can be used:




Actions menu

Filters can be used (*Toggle filter bar visibility*) to find the desired studies and study groups and make several mappings at the same time.




Make Study description/Study Group combinations

After the study groups have been mapped, alerts can be recalculated by clicking on **Actions** ⇒ **Recalculate notifications**.



Recalculate notifications

NOTE: Every time the user recalculates notifications, they will be prompted to select the notification types to be recalculated based on notification resolution type (e.g. “unknown”, “unresolved”, etc.). This way when notifications are recalculated the resolution will not be set back to “unresolved”.




*The default recalculation settings (deletion and data migration) can be configured in Settings. See [Notification Recalculation Settings](#) for more information. **Please be aware that recalculation of notifications may lead to lose some preexisting notifications if the recalculation configuration has not been carefully done.***

8.4.9.1.6. COMPLIANCE MONITORING TABLE

If the studies are linked to a category, we can add the **compliance monitoring** widget (a dashboard on device level) in order to evaluate the defined categories with the assigned studies. We can choose between the overview on study level or series level (acquisition protocol, series description, combinations between study description and series information). The compliance table is sorted automatically from maximum to minimum amount of studies per study group. The comparison column at the right side shows if the category is between the earlier defined limits.



Each record in the compliance monitoring has now two "traffic light" statuses: one which is conform with the date range and another one conforms with date range and last (default 100) studies ready for export to a template. The user can adjust it via changing the "Study count comparison" in the gear button. Here you can also save the column layout.



A Sparkline give the customer a sense in what direction (trend) the dose is going to. This Sparkline is showing the trend of the full date range with a resolution of 30 days, combined inside the traffic lights.

Change study count comparison

- RED – the studies are outside the acceptable limits
- ORANGE – the studies are outside the achievable but between the acceptable limits
- GREEN – the studies are between the achievable limits


COMPLIANCE MONITORING - STUDY LEVEL (DLP)

Parameter type: DLP | mGy.cm Dri Version: Version 2014-2018


Study Group	Count (#)	Min	Perc. (25)	Mean	Median	Perc. (75)	Max	Comparison
Head	943	174.56	626.14	661.83	626.14	626.14	2.59E3	
Abdomen	310	95.74	432.42	664.58	538.86	742.91	3.56E3	
Chest	194	25.2	269.4	369.75	356.77	417.29	1.53E3	

Compliance monitoring based on study level

The parameter type (DLP, CTDI, DAP, etc) can be changed in the dropdown menu and the default/favorite one can be saved by clicking on the engine button on the upper-right part and selecting *Layout Settings/Save column layout*.



Select parameter type




Save favorite parameter type

The DRL version against which the studies are being compared can also be selected from the dropdown menu for *Drl version*.



Select the DRL version

If a category is double clicked in the compliance table, studies assigned to the category can be viewed.



Studies assigned to compliance category

NOTE: the compliance bar color displayed in this list corresponds to the compliance status of the selected parameter and NOT the color of the highest notification, as default in the rest of DOSE.



Details of studies not assigned to a study group

At top, an indication is made of the amount of studies which are not assigned to a study group yet and which studies are not applicable to a certain limit. Clicking on the indications gives an overview of the studies which are not assigned yet or don't correspond to a certain limit. There is the possibility to exclude the free, stationary and/or series with small scan length (<20mm) out of the compliance statistics table, to do so, click on the blue "i" button and uncheck the boxes.



When the compliance is configured to SERIES level, you must use the Compliance monitoring table at series level and when clicking on a study, color bars will appear also at series level. In this case, the color bar at study level will be the one with the highest severity (e.g. 2 series, one orange and one red, color at study level will be red).

Study Details

Patient Name: Unknown

Overview | **Series information** | Organ doses | History | Activity stream | Notifications

DLP (Head16 phantom) : N/A | CTDI(Head16 phantom) : N/A | DLP(Body 32 phantom) : 547.76 | CTDI(Bod


Compliance	Date time	Acq. Protocol	CTDI [mGy]	DLP [mGy*cm]	Peak vo
	21/10/2013 22:27:17	Topogram AP	0.14	6.1	120
	21/10/2013 22:27:17	Abd venus	11.67	541.66	100

Color bars at series information

NOTE: In general the compliance bar will correspond to the most severe color of all the parameters mapped to that series. For example, if there is an orange alert for DLP and a red alert for CTDI, then the color bar of the series will be red. The only exception to this is opening studies from the compliance monitoring table, in this case the color of the parameter selected in **Parameter type** will prevail.


There are different user rights concerning Compliance. A user only with **Compliance** rights is able to visualize the compliance table, but not edit its configuration (mapping). **Compliance Management** rights are needed to edit the configuration.

The **Export to template** option within the gear button can be used to export the compliance table to one of the predefined legal export templates.




Export to template

All the available templates can be found here (*e.g. Belgian, Polish, Danish, Swedish*) and the export filename can be changed. For CT, localizers (topogram), free and stationary acquisitions and short length acquisitions can be excluded from the export. Desired adult and pediatric group templates can be selected.




Export in legal template

When exporting the mammography compliance to the Belgian template, tomosynthesis can be included. Tomosynthesis series are currently identified as containing “tomo” in their series description.



Include tomosynthesis series in export to legal template



Generating file

The result will be an excel file where each category is defined in a worksheet. For every category, the studies are sorted ascending by study date.


For more information, refer to the videos in the *DOSE Advanced Training/Legal compliance* section of our online training center.

8.4.10. Optimization Gain Prediction widget

Among many other charts, tables and widgets, there is a particular one whose working needs some explanation in this manual. It is the "Optimization gain predictor", that can be added in any CT dashboard.

This widget evaluates the impact on DLP when CTDI or scanlength are reduced for a particular type of exams (based on bodypart/study description/etc) in an optimization procedure.

For example: by selecting CTDI 5% (relative) and scanlength 10mm absolute you can see how the total doses and the dose of each study description or body part would change (in DLP) if you lowered the CTDI of each case by 5% and/or shortened the scanlength of an examination.



You can fill one of the fields with "0" to simulate variations only on CTDI or scanlength. The widget takes into account only helical/axial acquisitions (to avoid having localizers and stationary acquisitions affecting the results).

The change in DLP is calculated for all axial/helical acquisitions for the selected category (e.g. bodypart, study description, etc).

$$Total\ DLP(mGy \cdot cm) = \frac{CTDIvol(mGy) * Scanlength(mm)}{10}$$

Analysis results

	Total dose DLP (mGy*cm)	Total scanlength (mm)
Before optimization	3692538	2925934
After optimization	3333322	2823392
Gain (absolute)	359216.779	102542.000
Gain (relative)	9.728 %	3.505 %

Category	Count (#)	Relative count (%)	Total DLP (mGy*cm)	Total Scanlength (mm)	Gain Total (in DLP %)
RAD ct hersen 01,[31-7-2-173]	945	18 %	5.27E5	1.58E5	16 %
RAD ct comb 02,[18-2-2-194]	848	16 %	7.02E5	7.39E5	7 %
RAD ct abd 22,[1-2-2-187]	309	6 %	1.88E5	1.92E5	8 %
RAD ct comb 05,[18-2-2-194]	229	4 %	1.88E5	1.95E5	7 %
RAD ct thorax 24,[10-11-2-87]	193	4 %	6.24E4	7.05E4	8 %
RAD ct thorax 23,[10-11-2-87]	190	4 %	8.27E4	1.41E5	8 %

All those values are summed up and displayed above under "Total Dose DLP (mGy*cm)" for before and after the optimization. The Gain (absolute) is calculated by the difference of "Before-After", the Gain (relative) is calculated as follows:

$$DLP\ Gain\ (relative) = \frac{Total\ dose\ DLP\ before - Total\ dose\ DLP\ after}{Total\ dose\ before} * 100(\%)$$

The same logic applies to the gain in scanlength.

8.5. Optional features


These tools are enabled depending on customer's licenses agreements, and they can require extra connections and / or data provided by the customer.

8.5.1. Reject analysis

This tool is designed to help users reviewing the errors that lead to rejected images in general radiology and can be added as data panel to this kind of devices. This feature requires direct connection to the device, so that the system receives all the images performed by it.

Rejected analysis data panel consists of 4 preset dashboards, and others can be added as for the Dosimetry overview data panel.

The general overview contains one chart for daily rejection trend and the lists for reject reasons and series descriptions.



Reject analysis general overview


The KPI overview tabs show trends in terms of Study Related Rejection Rate (SRR), Instance related Rejection Rate (IRR), Dose related Rejection Rate (DRR) and Effective dose related Rejection Rate (ERR).



Daily SSR

Daily IRR

Validation overview is the only tab which cannot be edited. In this dashboard, the rejected image and its histogram are shown, so that the user can fill in the reason for it to be



rejected.

Validation overview of the rejected analysis

8.5.2. Contrast Administration Monitoring


When detailed contrast media data is available, it is displayed in the tab of interventional radiology or CT Study details page. If the data is available, patients' length/weight, the injected volume (ml), peak flow (ml/s) and peak pressure (Pa) can be analyzed.




Contrast administration monitoring

Additionally, direct integration with the contrast injectors from Ulrich Medical and Nemoto is possible. These injectors will appear as different "devices" of the CONTRAST modality. Integration with other vendors may be developed upon request.

Among the available charts are the injected volume per operator (can be study specific), daily/weekly/monthly trends for the used contrast volume, scatter plots of contrast volume versus patient age, weight or BMI and pie charts for contrast/operator name.




Average contrast volume and number of studies per operator




BMI vs contrast volume

Furthermore, for Ulrich injectors, Ulrich reports are integrated into the study details:



Pressure and flow rate graph



Injection sequence graph



8.5.3. Virtual Dosimetry

8.5.3.1. VIRTUAL DOSE FOR CT

In the 'Virtual Dose' tab, within 'Dose analysis' that appears when the user clicks on a CT series, phantoms and scan parameters can be changed to recalculate organ doses. By default, the system selects the appropriate phantom for the particular patient, the scanner used for the examination and parameters like kVp, mAs, collimation and pitch. The scan length is also shown as a red area on the phantom, and is automatically taken from the series information, together with the body part. If the bowtie filter type (head or body) is omitted, the head filter is assumed for scans of the head and neck region, and the body filter is assumed for all other regions.

The user can select one of the phantoms from the library (including adult/paediatric, normal weight/overweight/obese I and II/morbidly obese, male/female and pregnant woman in three gestational stages 3m-6m-9m). Then, select an automatic body region or manually modify the scanned area on the phantom. The software includes all vendors and many scanners. When all the fields are filled in, the user can click on the "Calculate Dose" button, and resulting organ doses will be displayed in the Dose Analysis tab.

For example, for a pregnant woman, fetal doses can be calculated and exported as a report to provide to gynaecologists.



Simulation results of organ doses

- The scan length of the series is taken into account together with the body part, which is prerequisite. The user must check if the boundaries that are shown on the phantoms are as scanned on the patient, as this is performed automatically, and may have differences from the actual scanning of the patient. If there is tube current modulation, then the reconstructed volume is sent as scan length (and not the DICOM value) to match the tube current modulation entries. Patient height is included in the parameters, to better match the scan length on the patient. If this is not available, the default height of each phantom is shown.
- Tube current modulation (TCM) is taken into account in Virtual Dose calculations if it is available for the patient. In this case, option yes/no for TCM is open in the Virtual Dose input parameters and the user can check the doses in both cases. As the dosimetric calculation will be affected by the TCM, the user should ensure that the display on the phantom matches the scanned area on the patient. The user must be careful when sending scan length information, because the DICOM value for scan length includes the overscan, while the reconstructed volume does not. The comparison of doses should be performed on the same scan length.

Please note that the dose values are indicative and should be evaluated by an expert (medical physicist or radiation protection expert).

For more information, refer to the videos *All about Virtual Dose CT by Qaelum in 30 seconds!* and *How to access Virtual Dose for computer tomography* in our online training center.

8.5.3.2. VIRTUAL DOSE FOR IR

As interventional procedures usually consist of multiple irradiation events, it is not easy to perform the calculation manually. This task is facilitated in DOSE, as the integration of Virtual Dose for Interventional Radiology (VDIR) allows the automatic calculation without the need for manual input by the user.

The calculation of organ doses in VDIR requires fields that are either general for the whole study or they are specific to each irradiation event. Some of the necessary fields, like DAP, tube voltage, primary and secondary angle etc are extracted by the RDSR. There are however some fields that are not included in the RDSR, like the source to skin distance, and these fields are calculated by the in-house Peak Skin Dose module.

The phantom is selected based on age and gender for pediatric patients and age, gender and size (height/weight) for adult patients. Please note that patient height and weight are mandatory fields for the calculation of virtual dose IR for all age groups.

The main dosimetric parameter that is sent as input to VDIR is the DAP. The direction of the beam is defined for each irradiation event by the primary and secondary angle. The center of irradiation is defined by the target organ, which is a generic field applied to the whole study calculation.

8.5.3.2.1. VISUALIZATION

In a separate tab on study details from a XA study, the user can see the calculated effective and organ dose values, as well as the selected phantom and the target area.

Overview Series information Peak Skin Dose **Virtual Dose** History Activity stream Notifications

The selections below apply to all events


Patient Phantoms
 Obese-level-I Adult Male

Ray Direction
 Postero Anterior

Target Region
 Heart


Recalculate

Export report



Description of VirtualDoseIR in DOSE

The automatically selected phantom is the Obese-I Adult Male and the target region "Heart".




The dose values are displayed in a chart and in a table.

The user can perform simulations by changing the selected phantom or the target region. This can be used to evaluate the effect on a pregnant patient or if initially the size of the patient is not known.

8.5.3.2.2. MISSING FIELDS

In case some of the necessary fields are not included in the RDSR or cannot be calculated, then the calculation is performed without these irradiation events. In this case, the events


that were not considered for the calculation are displayed in a table. This way, the user is always informed about which dose percentage was not included in the calculation.



Example of a case where the DAP was zero in the RDSR for some irradiation events. In this case, these events were not included in the calculation. The user is notified by a window on the right side, where each event and the corresponding fields are reported

8.5.3.2.3. REPORTING

The user can export all the calculated values and the settings/parameters in excel. This can be performed by clicking the **Export Report** button.



A	B
Calculation result	
Phantom	Adult Male - RPI
R60	34.77
E60 (mSv)	47.93
E103 (mSv)	63.57
Organ doses (mGy)	
Bladder	0.09
Bone surface	55.89
Bone-marrow (red)	29.52
Brain	13.16
Breast	210.73
Colon	2.64
Gonads	0.04
Liver	24.28
Lung	180.81
Oesophagus	57.62
Salivary glands	2.46
Skin	23.6
Stomach	39.88
Thyroid	8.32
Peak skin dose	2552.33
Remainder Dose	34.02
Remainder Organs	
Adrenals	45.29
Extrathoracic Region	1.75
Gall bladder	25.39
Heart	127.74
Kidneys	12.46
Lens of the eye	0.51
Lymphatic nodes	43.8
Muscle	37.87
Oral mucosa	2.84
Pancreas	34.03
Small Intestine	1.09
Spleen	45.68
Thymus	64.24
Uterus/Cervix (female) or Prostate (male)	0.05
Peak skin dose	6873.52

Example of a generated report. The user can see the calculation but also the settings that were considered for the calculation

In the case of missing fields, an additional excel sheet includes the irradiation events that were not considered for the calculation.

	A	B	C	D	E	F	G	H
1	Excluded series							
2	Series Number	Series Datetime	Plane	Acquisition type	DAP	SID	SSD	Field size
3	10	41190.57014	Plane A	Fluoroscopy	0	103	55	10
4	12	41190.57029	Plane B	Fluoroscopy	0	106	55	10
5	36	41190.57846	Plane B	Fluoroscopy	0	106	55	10
6	53	41190.5808	Plane A	Fluoroscopy	0	103	55	10
7	57	41190.58119	Plane A	Fluoroscopy	0	103	55	10
8	82	41190.58655	Plane A	Fluoroscopy	0	107	55	10
9	83	41190.58655	Plane B	Fluoroscopy	0	106	55	10
10	93	41190.59095	Plane B	Fluoroscopy	0	101	55	10
11	119	41190.61179	Plane B	Fluoroscopy	0	101	55	10
12								
13								
14								
15								

Settings | Calculation result | **Excluded series** | +

An additional sheet in the exported report shows the irradiation events (series) that were not included in the calculation.

For more information, refer to the video *All about Virtual Dose IR by Qaelum in 30 seconds!* in our online training center.

8.5.4. Global Noise Level

As it is important for medical imaging to combine dose information together with a satisfactory image quality, the automatic calculation of an image quality index as described in the literature has been implemented in DOSE. In particular, the selected image quality parameter is a noise index called Global Noise Level (GNL), published in well-known journals. The purpose of this feature is to evaluate image quality in combination with other parameters, like dose, patient size and positioning, and provide all the important information for identification of outliers and optimization.


8.5.4.1. CALCULATION

The calculation methodology of GNL is described in the literature by Christianson et al. (AJR. 2015) and Ria et al. (Med Phys. 2017):

1. First the images are segmented into tissue types according to the Hounsfield Unit (HU) values:
 - a. Soft tissue
 - b. Bone
 - c. Air

Soft tissue and fat were merged by extending the considered HU range (Ria et al, Med Phys 2017).

2. For each pixel of a tissue type, the standard deviation is calculated in a central kernel. This results in a noise map with the standard deviation (in HU) of the surrounding kernel for each tissue pixel (Figure B).
3. A histogram of the tissue noise map (with the standard deviations) is generated. The global noise of that slice is determined as the mode of the histogram peak (Figure 1C).




Schematic overview of the noise detection algorithm from Christianson et al. (AJR 2015). A) Original Image, B) Noise map of soft tissue, C) Histogram of soft tissue noise map.

The calculation is performed for every convolution kernel and every slice thickness of each acquisition. These values are then normalized to a 5mm slice thickness, as described in Ria et al (AJR 2019).

8.5.4.2. VISUALIZATION


Global Noise Level can be used as an absolute number, according to recent studies (Ria et al. 2019). However, it is mainly meant to be used as a relative value to identify outliers and potential errors and initiate optimization projects (Dedulle et al. 2019).

For every CT acquisition (series), under the “Image Quality” tab, a table shows the values of GNL for every convolution kernel and every slice thickness. The GNL values are normalized to 5mm, although the original slice thickness is also reported. By clicking on each row, the user can see a colored chart correlating the CTDIvol to patient size (Water Equivalent Diameter-WED) and GNL for all series with the same acquisition protocol, the same convolution kernel and the same original slice thickness for this device. The chart follows the general date range selected on the system. The GNL is the parameter indicated by colors. The colored legend changes dynamically, which means that the outer colors correspond to the maximum and minimum values of the shown data points. There is the option to see the same chart for soft tissue, air and bone, by selecting a different tab. The chart is clickable and allows to open the corresponding studies.




Visualization CTDIvol versus WED and GNL to facilitate identification of outliers. The current acquisition is indicated by a larger point of different shape as now the color has a different meaning (color indicates GNL value based on the colored legend).

The GNL is also displayed as a function of CTDIvol and water equivalent diameter (WED) of the patient (attenuation-based metric for the patient size). The current acquisition is visualized in yellow in the cloud of data points.




Visualization of the GNL versus CTDIvol and WED for soft tissue, with the current acquisition indicated in yellow.

GNL is exportable in two ways; for the current acquisition, the table with GNL values can be exported from the actions button below the GNL table.



Series number	Convolution kernel	Slice thickness (mm)	GNL Soft tissue (HU)	GNL Air (HU)	GNL Bone (HU)
2	Br40d	3	12.1	31.4	113.9
3	Br40d	1	11.2	23.8	78.3

For a large-scale export of GNL values and combination with other parameters, GNL can be exported on Device level using the **Export data** functionality (See [Export](#)). Under Field configuration, GNL can be found as **Global Noise Level** in the column with series level fields.



By selecting the "Global Noise Level" and exporting, all the relevant fields are automatically exported to allow for further handling.


	A	B	C	D	E	F
1	SERIES NUMBER (Series)	SLICE THICKNESS (mm) (Series)	CONVOLUTION KERNEL (Series)	GNL SOFT TISSUE (HU) (Series)	GNL AIR (HU) (Series)	GNL BONE (HU) (Series)
2						

Example of exporting the Global Noise Level field in excel. Six columns are automatically exported.

For more information, refer to the video **All about Global Noise Level by Qaelum in 30 seconds!** in our online training center.

8.5.5. Nuclear Medicine


The Nuclear Medicine department can be monitored and managed by DOSE. The source of information can be the hotlab management system or the NIS (via HL7 connection), or the PACS (via DICOM header analysis). The source of information (e.g. the hotlab management system, the PET etc) will be shown as a separate “room” or “device” in the device level tree.



STUDY DATE	STUDY
2015-02-26 10:30	OSSELUX - Osted
2015-02-26 10:30	OSSELUX - Osted
2015-02-26 10:18	OSSELUX - Osted
2015-02-26 09:58	OSSELUX - Osted
2015-02-26 09:30	OSSELUX - Osted
2015-02-25 12:18	COEUR EFFORT
2015-02-25 10:30	OSSELUX - Osted
2015-02-25 10:18	OSSELUX - Osted
2015-02-25 09:58	OSSELUX - Osted
2015-02-25 09:30	OSSELUX - Osted

Different levels of analysis on the upper horizontal line (study, patient, device). Tree of devices shown vertically on the left side, where the Hotlab/PET etc is also included as a separate room/device.

Devices can be configured as “hybrid modalities” (e.g. PET/CT, SPECT/CT) so that DOSE can show linked studies of different modalities (e.g. NM and CT) in **Study details** (sharing the study UID but of a different modality). Users whose role contains the “Device management” functionality can edit the configuration in **Device Overview** → **Configuration** → **Hybrid modality**.



Configure device as hybrid modality

From the different levels of system, the user has the option to search for a specific study, a specific patient or a specific “device” (in this example: the hotlab) and get a full overview of the hotlab and the department in general.


The Nuclear Medicine module is only available to customers with the corresponding license.

8.5.5.1. DOSIMETRY

8.5.5.1.1. EFFECTIVE DOSE AND ORGAN DOSE CONVERSION FACTORS

Effective dose and organ doses are calculated for all radiopharmaceuticals available in the ICRP Publication 128 (Radiation Dose to Patients from Radiopharmaceuticals: a Compendium of Current Information Related to Frequently Used Substances).

The **Conversion factors** are displayed in **Settings**, under **Data tables**, in a separate tab, next to the one of Radiology. The unit of the conversion factors for effective dose is mSv/MBq and for organ doses mGy/MBq.



General Conversion factors		Nuclear Medicine Conversion factors			
Start Age	End Age	Study description	Radiopharmaceutical	Condition	Conv ICRP128
1	3	All	111In-labelled human immunoglobulin	N/A	0.99
3	7	All	111In-labelled human immunoglobulin	N/A	0.58
8	12	All	111In-labelled human immunoglobulin	N/A	0.41
13	17	All	111In-labelled human immunoglobulin	N/A	0.22
18	150	All	111In-labelled human immunoglobulin	N/A	0.17
1	3	All	111In-labelled monoclonal tumour-associated antibodies	F(ab') ₂ fragments	1.0
1	3	All	111In-labelled monoclonal tumour-associated antibodies	F(ab) ₂ fragments	1.8
1	3	All	111In-labelled monoclonal tumour-associated antibodies	Intact antibody	1.2
3	7	All	111In-labelled monoclonal tumour-associated antibodies	Intact antibody	0.64
3	7	All	111In-labelled monoclonal tumour-associated antibodies	F(ab) ₂ fragments	0.55
3	7	All	111In-labelled monoclonal tumour-associated antibodies	F(ab) ₂ fragments	0.54
8	12	All	111In-labelled monoclonal tumour-associated antibodies	F(ab) ₂ fragments	0.35
8	12	All	111In-labelled monoclonal tumour-associated antibodies	F(ab) ₂ fragments	0.36
8	12	All	111In-labelled monoclonal tumour-associated antibodies	Intact antibody	0.41
8	12	All	111In-labelled monoclonal tumour-associated antibodies	Intact antibody	0.33

Table of conversion factors for Nuclear Medicine under Settings tab.

In the displayed table, each row shows the effective dose conversion factor for a particular age range and radiopharmaceutical. In order to see or edit the corresponding organ dose conversion factors, **right click on the row** and **edit conversion factors**. A new window will show all the conversion factors of the particular case. In the case of an organ, for which a conversion factor is not provided by ICRP 128, the conversion factor is displayed as -1. This is often the case for salivary glands, for which only some radiopharmaceuticals have a conversion factor.

123I-iodide
123I-iodide
123I-iodide
123I-iodide
123I-iodide
123I-iodide
123I-iodide
123I-iodide
123I-iodide

Clear Medicine Conversion factors

Study description	Radionuclide
All	11C-labelled amino acids (generic model)
All	11C-labelled amino acids (generic model)

Edit conversion factors

Effective dose [mSv/MBq]	Kidneys [mGy/MBq]	Liver [mGy/MBq]
0.41	1.3	2
Adrenals [mGy/MBq]	Bone surfaces [mGy/MBq]	Brain [mGy/MBq]
0.53	0.57	0.11
Breasts [mGy/MBq]	Gall bladder wall [mGy/MBq]	Stomach wall [mGy/MBq]
0.14	0.6	0.31
Small intestine wall [mGy/MBq]	Colon wall [mGy/MBq]	Upper large intestine wall [mGy/MBq]
0.28	0.27	0.32
Lower large intestine wall [mGy/MBq]	Heart wall [mGy/MBq]	
0.21	0.29	
Lungs [mGy/MBq]	Muscle [mGy/MBq]	Oesophagus [mGy/MBq]
0.26	0.18	0.15
Ovaries [mGy/MBq]	Pancreas [mGy/MBq]	Red marrow [mGy/MBq]
0.22	0.53	0.69
Salivary glands [mGy/MBq]	Skin [mGy/MBq]	Spleen [mGy/MBq]
-1	0.1	2.2
Testes [mGy/MBq]	Thymus [mGy/MBq]	Thyroid [mGy/MBq]
0.095	0.15	0.12
Urinary bladder wall [mGy/MBq]	Uterus [mGy/MBq]	Remaining organs [mGy/MBq]
0.16	0.2	0.2

Save Discard

Window with the corresponding conversion factors that appears when right clicking on a row. Conversion factors can be edited. In the example, salivary glands have a factor of -1, meaning that there is no available conversion factor for this case.

Below the table and at the right, the **action** button includes “Insert data” option.

Insert Data
Refresh
Edit

Insert data by the actions button below the table allows for visualization and editing of the conversion factors or the corresponding parameters.

By selecting it, the user opens the window with all the entries of the table and above this, the valid template for new entries.

Template:

age_start,age_stop,study

description,radionuclide,radionuclideCondition,conv_icrp_128,adrenals,bone_surfaces,brain,breasts,gall_bladder_wall,stomach_wall,small_intestine_wall,colon_wall,upper_large_intestine_wall,lower_large_intes

heart_wall,kidneys,liver,lungs,muscle,oesophagus,ovaries,pancreas,red_marrow,salivary_glands,skin,spleen,testes,thymus,thyroid,urinary_bladder_wall,uterus,remaining_organ


Empty - or lines with a wrong format are ignored.

```
1 0.3 0. All, 111In-labelled human immunoglobulin, NA, 0.89, 1.0, 1.1, 0.58, 0.5, 0.88, 0.76, 0.74, 0.7, 0.74, 0.65, 1.4, 1.1, 1.0, 1.3, 0.58, 0.65, 0.69, 1.0, 0.67, 1.0, 0.39, 3.3, 1.8, 0.65, 0.76, 0.58, 0.69, 0.61,
3 0.7 0. All, 111In-labelled human immunoglobulin, NA, 0.58, 0.58, 0.55, 0.33, 0.27, 0.58, 0.44, 0.42, 0.41, 0.44, 0.36, 0.81, 0.64, 1.1, 0.69, 0.31, 0.37, 0.38, 0.58, 0.37, 1.0, 0.21, 1.9, 1.3, 0.37, 0.41, 0.33, 0.39, 0.34,
8 0.12 0. All, 111In-labelled human immunoglobulin, NA, 0.41, 0.38, 0.35, 0.2, 0.17, 0.39, 0.29, 0.27, 0.26, 0.27, 0.24, 0.54, 0.42, 0.75, 0.45, 0.2, 0.24, 0.25, 0.38, 0.25, 1.0, 0.13, 1.2, 1.0, 0.24, 0.25, 0.24, 0.26, 0.21,
13 0.17 0. All, 111In-labelled human immunoglobulin, NA, 0.22, 0.25, 0.23, 0.12, 0.11, 0.26, 0.19, 0.17, 0.17, 0.18, 0.15, 0.36, 0.28, 0.5, 0.29, 0.13, 0.17, 0.17, 0.25, 0.16, 1.0, 0.03, 0.81, 0.22, 0.17, 0.16, 0.18, 0.17, 0.14,
15 0.16 0. All, 111In-labelled human immunoglobulin, NA, 0.22, 0.25, 0.23, 0.12, 0.11, 0.26, 0.19, 0.17, 0.17, 0.18, 0.15, 0.36, 0.28, 0.5, 0.29, 0.13, 0.17, 0.17, 0.25, 0.16, 1.0, 0.03, 0.81, 0.22, 0.17, 0.16, 0.18, 0.17, 0.14,
```


Table C.41. Absorbed doses for ^{51}Cr -ethylenediaminetetraacetic acid.

Organ	Absorbed dose per unit activity administered (mGy MBq^{-1})				
	Adult	15 years	10 years	5 years	1 year
Intravenous administration, normal renal function					
Adrenals	$7.3\text{E}-4$	$9.1\text{E}-4$	$1.4\text{E}-3$	$2.1\text{E}-3$	$3.9\text{E}-3$
Bone surfaces	$8.2\text{E}-4$	$1.0\text{E}-3$	$1.5\text{E}-3$	$2.1\text{E}-3$	$3.8\text{E}-3$
Brain	$4.8\text{E}-4$	$6.0\text{E}-4$	$9.8\text{E}-4$	$1.6\text{E}-3$	$2.9\text{E}-3$
Breast	$4.3\text{E}-4$	$5.6\text{E}-4$	$8.2\text{E}-4$	$1.3\text{E}-3$	$2.5\text{E}-3$
Gallbladder wall	$7.9\text{E}-4$	$1.0\text{E}-3$	$1.7\text{E}-3$	$2.3\text{E}-3$	$3.4\text{E}-3$
Gastrointestinal tract					
Stomach wall	$6.9\text{E}-4$	$8.5\text{E}-4$	$1.3\text{E}-3$	$2.0\text{E}-3$	$3.4\text{E}-3$
Small intestine wall	$1.1\text{E}-3$	$1.4\text{E}-3$	$2.1\text{E}-3$	$3.1\text{E}-3$	$4.8\text{E}-3$

Example from the ICRP Report 128 of what is mentioned as Condition in DOSE.



General Conversion factors		Nuclear Medicine Conversion factors			
Start Age	End Age	Study description	Radiopharmaceutical	Condition	Conv ICRP128
18	150	All	^{51}Cr -ethylene	Show All	Show All
18	150	All	^{51}Cr -ethylenediaminetetraacetic acid	Intravenous administration - normal renal function	0.002
18	150	All	^{51}Cr -ethylenediaminetetraacetic acid	Intravenous administration - abnormal renal function	0.0047
18	150	All	^{51}Cr -ethylenediaminetetraacetic acid	Oral administration (H=0)	0.031

Example of a radiopharmaceutical with multiple "Conditions".

Thus, the same radiopharmaceutical and the same age group can have different conversion factors depending on the condition. This will be indicated in DOSE by separate rows in the conversion factor table.

Not all radiopharmaceuticals have different conditions in ICRP128. In the case of a radiopharmaceutical with only one set of conversion factors, the "Condition" field in DOSE will be "N/A".

"Condition" can also refer to the patient condition (e.g. normal - abnormal renal function) or to the type of examination (e.g. intravenous or oral administration, rest/exercise etc.).

The **Study description** field allows to differentiate conversion factors based on the study description. As the pharmacokinetics is not always the same, the conversion factors will not be the same either. This can be helpful in cases when the study description indicates the type of examination and the corresponding conversion factors can be matched. If there is no need to use different conversion factors for different study descriptions, then the "Study description" field in DOSE should indicate "All".

General Conversion factors		Nuclear Medicine Conversion factors			
Start Age	End Age	Study description	Radiopharmaceutical	Condition	Conv ICRP128
[x] = 18	Show All	Show All	99mTc-labelled phosphates	Show All	Show All
18	150	All	99mTc-labelled phosphates and phosphonates	Normal uptake and excretion	0.0049
18	150	Bone imaging_impaired kidney	99mTc-labelled phosphates and phosphonates	High bone uptake and/or severely impaired kidney function	0.0043

Example of mapping a study description to a radiopharmaceutical and to specific conversion factors.


The study description for a particular row can be edited either by clicking on **Insert data** or the **Edit** button (found in **Actions**). When using the latter, a dropdown list of study descriptions of nuclear medicine devices is shown.

Start Age	End Age	Study description
1	3	NM Lungventilation/lungperfusion, SPECT
3	7	THERAPIE THYRO:DE -131-I(capsule)
8	12	PET whole body
13	17	REFLUX GASTRIQUE 99mTc MAA SENTINELLES MAMAIREs 99mTc Nanocollo:des THERAPIE THYRO:DE -131-I(capsule)

Edit study description in NM conversion factors

8.5.5.1.3. MAPPING OF RADIOPHARMACEUTICALS


Nuclear medicine conversion factors for effective dose and organ doses use the standardized naming of radiopharmaceuticals, as reported in ICRP Publication 128. As different hospitals may use different names for the radiopharmaceuticals, a tool to map the original radiopharmaceutical, as reported in the data source, to the standardized radiopharmaceutical, as indicated in the conversion factors table, exists in **Settings**.



Radiopharmaceutical mapping under Settings.


Radiopharmaceutical mapping is crucial in order to perform dosimetric calculations. Otherwise the system cannot match the original radiopharmaceutical coming from the data source to the conversion factors for effective and organ dose calculations. This process needs to be performed only once unless there is a new radiopharmaceutical or a new radiopharmaceutical name.

The mapping can be performed as follows: click on the **Edit** (pencil) button in the bottom left corner.




Edit button in Radiopharmaceutical mapping.


Then select the + button to add a new row. First select the device in the “device” field on the right. Based on the device, the available original radiopharmaceuticals will appear in a dropdown list.



Perform the mapping by selecting the official name and condition using the dropdown menu for every original radiopharmaceutical, as shown in the example.



Select original radiopharmaceutical



Select condition

Then click the **Save** button (floppy disk).

If multiple original radiopharmaceuticals refer to the same official radiopharmaceutical, then this step must be performed for each one (e.g. Tc-99mMDP, HDP, 99mTc-Osteocis, Tc-99mMDP must correspond to 4 lines in the mapping, all mapped to the same official radiopharmaceutical).


For the condition, it is advised to use a default condition for a radiopharmaceutical (when multiple are available). This can be changed in an individual study depending on the condition of the patient or the particular examination. For example, when normal and abnormal renal function is available, it could be prudent to map the normal function and then change the individual conditions of patients with abnormal renal function, as otherwise dose values will be under- or overestimated for all patients. Thus, it is good to know that changing radiopharmaceutical and condition is also possible for an individual study, and particularly from series level ([Series information](#)).

Users whose role includes “Physical parameters management” or “Dose management” functionality can access the radiopharmaceutical mapping in Settings.

8.5.5.2. STUDY LEVEL

8.5.5.2.1. STUDY OVERVIEW

Details of the study are shown in the **Overview** tab. This includes patient weight, administered radioactivity, administered radioactivity per kg, prepared radioactivity and radiopharmaceutical administration type (e.g. injection, ingestion, etc.). The radionuclide and radiopharmaceutical are indicated in separate fields. Effective dose and organ doses calculated by DOSE are displayed. Effective dose is displayed in the **Overview** tab, while organ doses are displayed in a separate tab “Organ doses” (see [Organ Doses](#)).




Study Details [Cancel] [Actions]

Overview | History | Activity stream | * Validation * | Linked documents


<p>NO PREVIEW AVAILABLE</p> <p>⚠️ ⚠️ ⚠️ ⚠️ ⚠️</p> <p>Category Bone Imaging</p> <p>Radioactivity (Administered) Limits > ACH (upper): 640.0 MBq > ACC (upper): 740.0 MBq</p>	<p>Study date 17/02/2015 10:30</p> <p>Study description OSSEUX - Osteocis 99mTc</p> <p>Hide study details</p> <p>Prepared activity [MBq] 650.00</p> <p>Body part examined TSPINE</p> <p>Manufacturer VEENSTRA</p> <p>Study ID [REDACTED]</p> <p>Operator Operator 579</p>	<p>Patient information F [REDACTED] 2A M 68 Y 95 kg</p> <p>Accession Number [REDACTED]</p> <p>Radionuclide 99m - Tc</p> <p>Modality NM</p> <p>Modelname HOTLAB</p> <p>Study UID [REDACTED]</p> <p>Performing physician Perf_Physician 74</p>	<p>Radioactivity (Administered) [MBq] 650.00</p> <p>Effective Dose ICRP103 [mSv] 3.71</p> <p>Administration type INJECTION</p> <p>Event # 1</p> <p>Software version SV_100</p> <p>Device Name Room 08 - NM Hotlab</p> <p>Requesting physician Req_Physician 32134</p>
	<p>Hide notifications</p> <p>⚠️ Obese patient ⚠️ NM: Artefact due to contamination STATIC LIMIT ⚠️ MBQ_INJECTED</p> <p>Hide clinical quality evaluations Hide lot purity Show study comparators</p>		

Overview of the study: parameters, such as patient weight, administered activity and effective dose are shown in Study Details. In this example the study description has been added to the study group “Bone imaging” and an alert for administered radioactivity has been generated.




Example of an examination where the administered activity is lower than the prepared activity. The user can add a comment (e.g. the patient showed up late for the appointment). Effective dose is calculated based on the activity that was actually administered to the patient

Additional studies of other modalities (e.g. CT) are shown in the **Hybrid** tab, appearing when a study with the same UID of the linked hybrid modality is present.




Hybrid study

It is possible to manually add NM studies to DOSE. To do so, the user must go to Device Level/NM device/Dosimetry Overview and click on **+Add study**:



A window will pop up allowing the user to fill in all the details. As many injections as necessary can be added (1) and activity details can be edited by clicking on the pencil button (2):



Add study

Study details

Patient ID: Patient Name: Patient DOB:


Gender: Weight [kg]: Length [cm]:

Study date: Accession Nr: Operator name:

Study description:

Events 1 **+ Add**

Nr.	Series description
1	Injection 1
2	Injection 2

Event details 2 

Event description:

Administration type:

Radionuclide:

Radiopharmaceutical:

Radiopharmaceutical condition:

Lot number:

Bodypart:

Prepared activity (MBq):

Initial volume (ml):

Remaining volume (ml):

Radiopharmaceutical volume:

Usage rate (%):

Administered activity (MBq):

Save **Cancel**



The usage rate and the administered activity are calculated according to the formulas below:

$$\text{Usage rate} = \frac{\text{initial volume} - \text{remaining volume}}{\text{initial volume}} * 100$$

$$\text{Administred activity} = \text{Prepared activity} * \text{Usage rate}$$

8.5.5.2.2. LOT NUMBER PURITY

When users perform a quality control of the pharmaceutical, LOT number purity can be registered in the corresponding DOSE Logbook. Examinations with this LOT number can then be linked and the pharmaceutical purity will be shown in the study overview.




- a) LOT number linked to this study together with the measured purity, b) when purity is below 95% an orange light is shown


For more information, refer to the video ***How to view nuclear medicine study information on Study Level*** in our online training center.

8.5.5.3. ALERTS & JUSTIFICATIONS

Notifications or alerts can be manually or automatically created for Nuclear Medicine and displayed in the study overview. The automatic notifications that relate to study groups and compliance are explained in [Compliance configuration](#). The automatic notifications that are generated by the Notification center are explained in the [Notifications Chapter](#). For manual notifications, the user can add a comment for the study, e.g. to justify a high activity or an artefact either through the [live dashboard](#) or via the study details (Tab [Activity stream](#)). Comments do not always have to create a notification; the user can decide. There is already a list of predefined comments especially for Nuclear Medicine for easy and quick handling (e.g. impaired renal function) but additional ones can be also added.




Live dashboard for Nuclear Medicine. The user can add a comment in the comment field and select if a notification is needed or not. Under the options button, predefined comments, general and specific to Nuclear medicine, are included



Predefined comments for Nuclear Medicine. The user can select the one that fits or create a new one and Save.

Alerts / notifications can be generated from the comments and handled through the [Notifications Module](#).



Study Details

Overview | History | Activity stream | *Validation* | Linked documents

NO PREVIEW AVAILABLE

Study date: 17/02/2015 10:30

Study description: OSSEUX - Osteocis 99mTc

Category: Bone Imaging

Radioactivity (Administered) Limit(s):
 > ACH (upper): 940.0 MBq
 > ACC (upper): 940.0 MBq

Radioactivity (Administered) [MBq]: 650.00

Effective Dose ICRP103 [mSv]: 3.71


Alerts: **NM: Artefact due to contamination**, **Obese patient**

Example of alerts in a study. The contamination could cause an artefact that affects the image.

For more information, refer to the video *How to create nuclear medicine notifications using Live Dashboards* in our online training center.

8.5.5.4. SERIES INFORMATION


The evaluation of Nuclear Medicine examinations usually happens at study level. However, there are cases where the examination consists of multiple syringes with the same or different radiopharmaceutical. To avoid missing information, DOSE can show multiple syringes as different series. If the data source only includes information about one administered radiopharmaceutical, then there is one series with the same radioactivity and radiopharmaceutical as mentioned on study level. The administered radioactivity of study level is the sum of all the series. If there are different radiopharmaceuticals involved, then on study level the radiopharmaceutical is a concatenation in alphabetical order and the effective dose is the sum of effective dose from different series. Details of the series can be exported by the Actions button below the series table.



Content datetime	Effective dose	Injected activity	Prepared activity	Lot number	Radiopharmaceutical Start time	Radiopharmaceutical Stop time	Radiopharmaceutical volume	Injection time	Radiopharmaceutical	Radiopharmaceutical condition	Syring
2/07/2018 01:49 PM	3.45E-3	0.5	0.45	95000	3/09/2020 10:44 AM	3/09/2020 10:44 AM	0.78	3/09/2020 10:44 AM	99mTc-labelled tetrofosmin (Myoview)	Exercise	syring001
2/07/2018 01:49 PM	3.45E-3	0.5	0.45	95000	3/09/2020 10:44 AM	3/09/2020 10:44 AM	0.78	3/09/2020 10:44 AM	99mTc-labelled tetrofosmin (Myoview)	Exercise	syring001
2/07/2018 01:49 PM	3.45E-3	0.5	0.42	99000	3/09/2020 10:44 AM	3/09/2020 10:44 AM	0.3	3/09/2020 10:44 AM	99mTc-labelled tetrofosmin (Myoview)	Exercise	syring002


Series analysis for one examination. Information for every syringe, like radiopharmaceutical, administered activity, prepared activity, volume, administration time, effective dose etc. are displayed.

If a series needs to be edited, this can be done by opening the 'Series Information' tab and right-clicking on the particular series to select 'Edit Serie Data'.



Effective dose is calculated for every series based on the age of the patient, the radiopharmaceutical of the series, and the study description (in case a link with the specific study description has been made in the [Conversion factor table](#) (as described in [Dosimetry for different conditions and study descriptions](#)). For more flexibility in the dosimetry of each patient, there is the option to change the radiopharmaceutical and/or the condition for each series of the individual study. By right clicking the selected series and **edit radiopharmaceutical**, a pop-up window allows to change the radiopharmaceutical (if the data source was not correctly filled) or the condition (if the patient differs from the default condition as selected in the [mapping of the radiopharmaceuticals](#)). The user can select if the change will apply to the current series or to the whole study and click **Save**.


Content date/time	Effective dose	Injected activity	Lot number	Radiopharmaceutical Start time	Radiopharmaceutical volume	Injection time	Radiopharmaceutical	Radiopharmaceutical condition	Syringe id	Prepared activity ref date time	Radionuclide
07/02/2018 13:49:09	0.01	0.5	99000	09/03/2020 10:44:44		09/03/2020 10:44:33	131I-labelled monoclonal antibodies	F(ab)' fragments	syring001	09/03/2020 10:43:49	125I
07/02/2018 13:49:09	0.01	0.5	99000	09/03/2020 10:44:46	Edit radiopharmaceutical	09/03/2020 10:44:35	18F-fluoride	N/A	syring002	09/03/2020 10:43:51	18F



Right click on the series to edit radiopharmaceutical and/or condition for the selected series and/or study helps to make the calculation of effective and organ doses more patient-specific.


8.5.5.5. ORGAN DOSES

Organ doses are calculated on study level by summing the organ doses of individual series. They are displayed in a tab called “Organ doses” next to the Series information. The visualization is similar to the in-house Organ dose calculation in CT (see [Organ Doses](#)).



Organ dose calculation for a specific study (red bars). The results displayed are the sum of absorbed doses from individual series for each organ.


The organ doses of the study (red bars) can be compared against the same radiopharmaceutical, the same study description or the same radiopharmaceutical and condition to allow for more relevant comparisons. The selected filter applies on the calculation of the 25th, 50th and 75th percentile, which can be compared with the current study.



The organ doses of the current study can be compared with similar studies based on the criteria of the right drop-down menu.

8.5.5.6. MANAGEMENT LOGBOOKS

Multiple logbooks must be kept in Nuclear Medicine in order to manage all the aspects of daily routine. Management logbooks centralize all data, link it to the studies, and export when needed. An overview of issues is possible by checking the statistics above every logbook.



REFERENCE DATE	LOT NUMBER	PURITY (%)	EXPIRATION DATE	QC COMMENT	ACTIVE
2015-02-11	1258945	95.4	2017-02-08		NO
2015-02-16	1258999	96.4	2017-02-22		YES
2016-06-14	122879	97.4	2016-06-14		NO
2016-09-27	123456	99.0	2016-09-27		YES
2016-09-29	147895	94.4	2016-09-29	Contacted the vendor	YES
2016-09-29	147896	95.6	2016-09-29		YES

Quality Control of pharmaceutical. The purity is recorded. Double click on an entry will provide the studies that were performed under this LOT number. Purity below 95% is highlighted with orange



CONTAMINATION DATE	RADIONUCLIDE	HALFLIFE (H)	LOCATION	SIZE	CONTAMINATION LEVEL	CONTAMINATION AFTER DECONTAMINATION	DECONTAMINATED	DECONTAMINATION DATE	CLEARED BY	COMMENT
2016-10-21 12:42	I-131	192.48	Hollab floor	medium	500.0 Bq/cm ²	10.0 Bq/cm ²	NO			
2016-10-05 18:24	I-131	192.48	Hollab floor	medium	500.0 Bq/cm ²		YES	2016-10-06 18:24	Thomas	
2016-09-17 15:07	I-131	192.48	Hot wc floor	small	260.0 Bq/cm ²	50.0 Bq/cm ²	NO			
2016-09-16 12:24	Tl-201	73.0	Hot wc floor	large	450.0 Bq/cm ²	60.0 Bq/cm ²	NO			Inform cleaning lady
2016-09-13 14:37	Tl-201	73.0	Hollab floor	medium	700.0 Cps	120.0 Cps	YES	2016-09-15 22:02		
2016-09-10 09:18	F-18	1.83	Camera table	small	500.0 Bq/cm ²	80.0 Bq/cm ²	YES	2016-09-11 06:04		From clothes

Contamination logbook. Every contamination event is recorded with the radionuclide, the contaminated area and the person who handled the problem. This could give interesting information, like for which radionuclide the personnel could use extra training



REFERENCE DATE	ALERT LEVEL	LOCATION	EVENT TYPE	EFFECT	CONTACTED COMPANY	EVENT STATUS	RESOLVED	SOLUTION	SOLUTION DATE	COMMENT
2016-10-04 14:39	Error	Room 1	PC	System works but with problems	YES	Under Investigation	NO			Error 1789
2016-09-29 15:05	Malfunction	Room 1	Collimator	No effect	NO	Under Investigation	NO			Hardware
2016-09-28 11:15	Malfunction	Room 2	Table	System still works ok	YES	Resolved	YES	Technician visit	2016-09-29 11:15	
2016-09-27 11:15	Error	Room 1	Collimator exchange system	System works but with problems	YES	Reported	NO			Error 512

Functional logbook. The user records errors, malfunctions etc by the dropdown menus (for all the cameras). The user can check if the company was contacted and what is the status of the events. The time needed to solve malfunctions or damages is estimated



REFERENCE DATE	BIN NAME	RADIONUCLIDE	HALFLIFE (H)	CONTAINER	ACTIVITY ON REF DATE (MBQ)	INITIAL VOLUME (ML)	EXTRACTED VOLUME (ML)	REMAINING VOLUME (ML)	CURRENT ACTIVITY (MBQ)	REF ACTIVITY (MBQ)	ESTIMATED LIBERATION DATE	LIBERATED	LIBERATION DATE
2017-01-07 15:57	Cardio bin	Ga-67	78.24	Syringe	185.0	1.2	1.1	0.1	2	0.37	2017-01-25	NO	
2017-01-01 15:24	Iodine bin	I-131	192.48	Vial	200.0	1.0	0.9	0.1	5	0.37	2017-02-17	NO	
2016-12-23 15:24	Cardio bin	Ga-67	78.24	Vial	200.0	2.0	1.2	0.8	0	0.37	2017-01-18	NO	
2016-10-03 16:27	General	Sm-153	46.3	Vial	2220.0	1.4	1.3	0.1	0	0.37	2016-10-20	YES	2016-11-25

Nuclear Waste Management. The user fills in the waste from radionuclides with half-lives that exceed 1 day. The system calculates the date that the waste for each bin can be liberated, when the reference activity is below 0.37MBq (10µCi).

For more information, refer to the video *How to view nuclear medicine data on Device Level* in our online training center.


9. Live dashboard

The Live dashboard can be found within each device on Device level. The Live dashboard tab shows the latest studies performed using the selected device. For each study, users can see the study information and the dose histogram. The orange arrow in the dose histogram shows the position of the current study dose among other studies with the same study description.



Live dashboard


This allows immediate direct feedback to the radiographer following an examination.



Add comment to the examination


There is a possibility to look into the study details by clicking on the information button and add a comment to this study. You see immediately if this study contains alerts or not. In addition, a clinical quality evaluation can be done immediately from the Live Dashboard.

The users can configure the Live Dashboard view via Action Button in the right corner.



Configure dashboard

Another way to access the Live Dashboard is from the dropdown menu on the top-right corner. Users whose role contains “Live Dashboard” and/or “Live Dashboard Modality” functionalities can see this. This is recommended for radiographers.



In “Live Dashboard” select the desired device, and click “Load Dashboard”.











While in “Live Dashboard Modality” select a modality and click “Load”.

A list of the last studies performed by any device of the selected modality will appear. The toggle bar filter can be used to filter by device and see as many studies as necessary (in “Live Dashboard” a maximum of 10 studies can be displayed).


Live Dashboard Modality

Select a modality: Computed Tomography (CT) Load Overview last updated at 27/05/2021 16:50:39 (100 studies loaded) Actions

Thumbnail	Alerts	Justification	Study date	Patient gender	Patient age	Modality	Device Name	Study description	BodyPart Examined	Total dose	E
			06/06/2018 12:26	M	80 Y	CT	TOSHIBA Aquilion ONE	RAD ct thorax 23,[10-11-2-87]	ARM	1.32E3	1.0
			06/06/2018 14:02	M	47 Y	CT	TOSHIBA Aquilion ONE	RAD ct thorax 23,[10-11-2-87]	ARM	1.47E3	1.1
			06/06/2018 14:06	M	47 Y	CT	TOSHIBA Aquilion ONE	RAD ct abd 22,[1-2-2-187]	ABDOMENPELVIS	300	4.6
			06/06/2018 14:40	F	57 Y	CT	TOSHIBA Aquilion ONE	RAD ct abd 22,[1-2-2-187]	ABDOMENPELVIS	588.9	9.1

100 rows Toggle Filter Bar visibility Reset filter | Actions


The number of displayed studies can be configured in **Actions/Configure Study Overview**.




9.1. Notifications from the Live Dashboard

Notifications can be created in the Live Dashboard using predefined canned comments, or by typing free text. These canned comments can be configured to automatically generate notifications.

Using the **Add comment** window: enter a comment → tick the **Create notification from comment** box → indicate an 'Identifier key' (this can be used to categorize notifications and find them when using the search on the **Notifications** level).



If using a canned comment: choose one of the canned comments → create notification.




List of canned comments

For more information, refer to the video *How to use the Live Dashboards* in our online training center.

10. Modality level


This module contains two tools: modality comparison and modality export.




Modality level

10.1. Modality comparison

This dashboard can help users to compare the performance of different devices. Users can customize their own comparison, to do so, they need to add first a comparison or to load one of the predefined ones, by clicking on “+Add comparison”.




Add a new comparison




Add a predefined comparison

For a new comparison, choose the modality and the devices intended to compare from the left side of the screen.




*List of modalities and the corresponding sites/devices. Select by clicking the boxes or select all from the button **Select**.*

Then click on Configuration, here the user will be able to choose which kind of graph he/she wants to use for comparison.




Configuration window

In the Actions button, the user can save the comparison in order to make it predefined.



Actions button for handling the comparisons (save/remove)

For a more specific comparison, it is possible to apply filters: by filtering, it is possible to choose the type of variable the user wants to compare, or to compare studies with a specific study description, study group or involving a specific body part. The filter works with the same logic as on **Device** level (see [Advanced Filter](#)).



Filter for the selected configuration

If the selected devices or the configuration are changed, then the refresh button will update the comparison.




Refresh button for changes in the comparison

For more information, refer to the video *How to use the Modality Level* in our online training center.


10.1.1. List of configurations

In the configuration window, a list of predefined comparisons is available. Histograms, tables, statistical analysis, trend analysis, either general or based on study description, protocol name, body part are available.




Window to select the comparison type

Every time a new comparison is selected, the **Add** button will include it in the list. Multiple comparison types can be added at once. At the end of the selection, **Update the configuration**.




Add button to add the selected comparison type in the current comparison page. When the desired comparison types have been selected, the button update configuration shows the comparison page.

In some cases, a comparison type may need adjustments. For example for a bin size that is too big, a smaller one can be selected. This happens by the **configure** button of the specific comparison type.




The 'Configure comparisons' dialog box features a title bar with a close button. Below the title bar is a header with a pencil icon and the text 'Select and configure all different comparisons'. The main area contains a 'Comparison' dropdown menu set to 'Frequency histogram and Protocol name' with a '+ Add' button to its right. Below this is a 'Description' field with the same text. A horizontal separator line follows. The next section contains a 'Percent histogram comparison' icon, a 'Remove comparison' button with a 'Down' arrow, and a blue 'Configure' button highlighted with a red rectangle. Below that is a 'Statistical comparison' icon, a 'Remove comparison' button with an 'Up' arrow, and a blue 'Configure' button. At the bottom, there is a red 'Cancel' button on the left and a blue 'Update configuration' button on the right.



The 'Modality overview configuration' dialog box has a title bar with a close button. It contains several input fields: 'Panel title', 'Chart title', and 'Chart subtitle', each with a dark text input field. Below these is a 'Chart bin size' field with a numeric input set to '150'. The 'Chart data restriction' field is a dropdown menu currently set to 'No data restriction in charts'. At the bottom, there is a blue 'Change' button on the left and a red 'Cancel' button on the right.

To configure a specific comparison type, select the configure button. This can be useful for example to change the bin size of a histogram.



Example of a histogram comparison

To easily identify the site (hospital) to which each device belongs, the 3-character tag that was defined in the [Add Site](#) section is included before the name of the device.



Site tags before the name of the devices

10.1.1.1. COMPLIANCE BENCHMARK

One of the available comparison types is the **compliance benchmark**. This table is the same as on device level (see [Compliance Monitoring](#)) but includes and compares data for all study groups and all devices (linked to these study groups) for the specific modality.

Study groups	Count	Min	Perc. (25)	Mean	Median	Perc. (75)	Max	Comparison
06b_CT_Stokdarm_Abdomen	14	275	534.25	637.286	660	724.75	1078	ALL 14
03b_Abdomen_Totaal_2F	5	662	674	829.6	689	1055.5	1353	ALL 5
▼ Demo Hospital 1								
▼ Room 05 - Siemens Somatom Definition Flash (Study description)	246	69.69	459.258	897.137	600.125	984.05	7277.33	ALL 100
RAD ct abd 21,[1-2-2-187]	121	69.69	461.245	896.39	589.65	986.595	5955.22	ALL 100
RAD ct abd 22,[1-2-2-187]	125	199.62	454.05	897.86	602.69	1018.915	7277.33	ALL 100
▼ Colon (Adult)								
▼ Demo Hospital 1								
▼ Room 05 - Siemens Somatom Definition Flash (Study description)	127	88.11	1010.11	1446.768	1303.17	1748.18	4025.42	ALL 100
RAD ct uro 20,[48-19-2-177]	127	88.11	1010.11	1446.768	1303.17	1748.18	4025.42	ALL 100
▼ Room 02 - Siemens Somatom Sensation 64 (Study description)	123	65.5	103.87	138.516	123.3	136.14	1800.77	ALL 100
RAD ct uro 13,[48-19-2-179]	123	65.5	103.87	138.516	123.3	136.14	1800.77	ALL 100
▼ Head (Adult)								
▼ Demo Hospital 1								
▼ Room 05 - Siemens Somatom Definition Flash (Study description)	605	172.7	556.37	748.437	781.73	804.655	4161.02	ALL 100
RAD ct hersen 04,[31-7-2-173]	137	172.7	307.59	435.142	333.18	360.12	4161.02	ALL 100
RAD ct hersen 01,[31-7-2-173]	468	357.75	693.322	840.15	803.56	815.62	3616.96	ALL 100
▼ CT pediatric chest								

Compliance benchmark for CT.

The table shows, not only how the different scanners perform for each study group, but also how the mapping was done (based on protocol name, acquisition protocol etc) and how the individual protocol performs.

For each study group (in this example “Thorax adult”), the devices linked to this study group are displayed. In the following example, 1 CT scanner from Demo Hospital 2 and 2 CT scanners of Demo Hospital 1 are linked to “Thorax adult”. In the case of Room 1, the mapping was based on protocol name, while for the other two scanners, the mapping was based on the Study description.



Study groups	Count	Min	Perc. (25)	Mean	Median	Perc. (75)	Max	Comparison
▼ Thorax (Adult)								
▼ Demo Hospital 2								
▶ Room 01 - SIEMENS SOMATOM Definition Flash (Protocol name)	29	67	249.5	317.552	294	362	707	ALL 29
▼ Demo Hospital 1								
▶ Room 05 - Siemens Somatom Definition Flash (Study description)	302	111.57	179.667	250.844	205.1	247.625	4605.13	ALL 100
▶ Room 02 - Siemens Somatom Sensation 64 (Study description)	194	25.2	269.4	369.748	356.77	417.287	1530.26	ALL 100

As seen, for Room 1, there is 1 protocol linked to the study group, which can be seen below the name of the CT scanner.



Study groups	Count	Min	Perc. (25)	Mean	Median	Perc. (75)	Max	Comparison
▼ Thorax (Adult)								
▼ Demo Hospital 2								
▶ Room 01 - SIEMENS SOMATOM Definition Flash (Protocol name)	29	67	249.5	317.552	294	362	707	ALL 29
07b_Thorax_TPS_3mm	29	67	249.5	317.552	294	362	707	ALL 29
▼ Demo Hospital 1								
▶ Room 05 - Siemens Somatom Definition Flash (Study description)	302	111.57	179.667	250.844	205.1	247.625	4605.13	ALL 100
RAD ct thorax 28,[10-11-2-176]	302	111.57	179.667	250.844	205.1	247.625	4605.13	ALL 100
▶ Room 02 - Siemens Somatom Sensation 64 (Study description)	194	25.2	269.4	369.748	356.77	417.287	1530.26	ALL 100
RAD ct thorax 24,[10-11-2-87]	194	25.2	269.4	369.748	356.77	417.287	1530.26	ALL 100

In the following example of the study group “Abdomen”, the particular scanner (Room 5) has two study descriptions linked to this study group and thus the values shown for the specific scanner (upper row) take into account both study descriptions. (The number of studies is the sum of studies of both scanners, the minimum value is the minimum of both, the average is the average of both).



Study groups	Count	Min	Perc. (25)	Mean	Median	Perc. (75)	Max	Comparison
▼ Room 05 - Siemens Somatom Definition Flash (Study description)	246	69.69	459.258	897.137	600.125	984.05	7277.33	ALL 100
RAD ct abd 21,[1-2-2-187]	121	69.69	461.245	896.39	589.65	986.595	5955.22	ALL 100
RAD ct abd 22,[1-2-2-187]	125	199.62	454.05	897.86	602.69	1018.915	7277.33	ALL 100
▼ Demo Hospital 1								
▶ Room 05 - Siemens Somatom Definition Flash (Study description)	246	69.69	459.258	897.137	600.125	984.05	7277.33	ALL 100
RAD ct abd 21,[1-2-2-187]	121	69.69	461.245	896.39	589.65	986.595	5955.22	ALL 100
RAD ct abd 22,[1-2-2-187]	125	199.62	454.05	897.86	602.69	1018.915	7277.33	ALL 100


The table can be exported in the Actions button, right above the table.



Study groups	Count	Min	Perc. (25)	Mean	Median	Perc. (75)	Max	Comparison
Room 05 - Siemens Somatom Definition Flash (Study description)	246	69.69	459.258	897.137	600.125	984.05	7277.33	ALL 100
RAD ct abd 21,[1-2-2-187]	121	69.69	461.245	896.39	589.65	986.595	5955.22	ALL 100
RAD ct abd 22,[1-2-2-187]	125	199.62	454.05	897.86	602.69	1018.915	7277.33	ALL 100

10.1.1.2. CONFIGURABLE BENCHMARK FOR CT AND MAMMOGRAPHY

For CT and mammography, there is the option to add a configurable comparison. This feature is called “**Benchmark**” and can be found in the **Configuration**.



Configure comparisons

Select and configure all different comparisons

Comparison: **Benchmark** + Add

Description: **Benchmark**

Benchmark Remove comparison


Configure

Cancel Update configuration

Benchmark: configurable comparison.

This configurable comparison gives the option to select many variables; protocol or study group, series level or study level, chart type, parameter(s), statistical measure.

Type: depending on the device and possibilities, users can select between Study group (on study or series level), Study or series description, Protocol or acquisition protocol, depending on what they want to compare. Only devices that have the selected group/description/protocol will have values in the chart. Selecting between study or series level will affect the list of available parameters. For example, for CT devices, SSDE is only available on series level.



Category: depending on the selection in “Type”, the options will appear.


Chart: Column chart, Scatter plot or Box and whisker plot. Depending on the selection of chart, different options open for parameter selection. A scatter plot needs two parameters; for x and y-axis. A box and whisker plot needs one parameter. A column chart gives the option to select two parameters (left y-axis and right y-axis) and also the statistical measure to compare (mean, median, percentiles, total etc).



The **Generate chart** button will provide the configured comparison.



Example of CT benchmark: a column chart with two parameters, for the left one (effective dose) the median value is compared and for the right one (volume) the total number is compared.



Example of mammography benchmark: a scatter plot representing the organ doses as a function of breast thicknesses.


10.2. Modality Export

This tool enables users to export study data from several devices of the same modality at the same time. There are different options available depending on your country: Modality Export (general), SSM Export (Sweden) and Full Export (Norway).

10.2.1. Modality export

This feature is similar to the one described in section [Export](#).


To access it, click on "**+Add Export**", enter a name and select "**Modality Export**" from the dropdown menu.



Add export


Select the modality and the desired devices before clicking "**Export**", where a new window will open where the desired export fields can be selected.

Please note that when exporting PATIENT_AGE, PATIENT_AGE_TYPE should also be selected to indicate if the age is in years, months, etc.



Site	
▼ Demo Hospital 1	<input checked="" type="checkbox"/>
▼ Devices	
▶ Room 05 - Siemens Somatom Definition Flash	<input checked="" type="checkbox"/>
▶ Room 02 - Siemens Somatom Sensation 64	<input checked="" type="checkbox"/>
▶ Room 04 - Toshiba Aquillon 16	<input checked="" type="checkbox"/>
▶ SIEMENS SOMATOM Force Data An	<input checked="" type="checkbox"/>
▶ GE MEDICAL SYSTEMS BrightSpeed	<input checked="" type="checkbox"/>
▶ SIEMENS SOMATOM Definition Flash	<input checked="" type="checkbox"/>
▶ SIEMENS SOMATOM Definition Flash 2	<input checked="" type="checkbox"/>
▶ SIEMENS SOMATOM Definition AS+	<input checked="" type="checkbox"/>

Select a modality and the devices



Select the desired fields to export


Then click on **“Export”** and the data spreadsheet will be generated.

For more information, refer to the video *How to export on Modality Level* in our online training center.


10.2.2. SSM Export

The SSM file is meant to comply with the Swedish regulation. Each report consists of multiple Excel files (one per device) that are grouped together in a zip file for one download. SSM export can be accessed in Modality level following the steps below:

1. Click on **Add export** in the upper right corner.
2. Provide an export name, and select **‘SSM Export’** from the dropdown menu in **‘Export type’**.



3. Click on **Add export**.
4. Select the modality from the dropdown menu.
5. Select the necessary devices from the list and click on **'Export'**.



6. Once the file has finished generating, a **'Download'** button will appear.


10.2.3. Full export

The NPR report was developed to be exported on **Modality Level**, so information about several devices and modalities can be exported at once. The resulting export is in XML format and can be included in the required general XML report for the Norwegian authorities.


```
<?xml version="1.0" encoding="UTF-8" ?>
<MsgHead xmlns="http://www.kith.no/xmlstds/msghead/2006-05-24" xmlns:xs1="http://www.w3.org/2001/XMLSchema-instance" xmlns:xsd="http://www.w3.org/2001/XMLSchema">
  <MsgInfo>
    <Type V="NPR_RAD" DN="NPR Innrapportering av data fra bildediagnostikk, intervensjon og nukleærmedisin"/>
    <MIGVersion>v1.2 2006-05-24</MIGVersion>
    <GenDate>2022-04-29T14:52:49</GenDate>
    <MsgId>3c9e69cd-3f8f-418a-ae56-447548cfc3fe</MsgId>
    <Sender>
      <Organisation>
        <OrganisationName>DEMO HOSPITAL</OrganisationName>
        <Ident>
          <Id>123456789</Id><!-- Organisation ID -->
          <TypeId V="ENH" S="2.16.578.1.12.4.1.1.9051" DN="Organisasjonsnummeret i Enhetsregister"/>
        </Ident>
        <Ident>
          <Id>123456</Id><!-- HER ID MHH -->
          <TypeId V="HER" S="2.16.578.1.12.4.1.1.9051" DN="HER-id"/>
        </Ident>
      </Organisation>
      <OrganisationName>DEMO HOSPITAL</OrganisationName> <!-- 924445114 org Number -->
      <Ident>
        <Id>123456</Id> <!-- HER ID MHH Radiology -->
        <TypeId V="HER" S="2.16.578.1.12.4.1.1.9051" DN="HER-id"/>
      </Ident>
    </Organisation>
  </Sender>
  <Receiver>
    <Organisation>
      <OrganisationName>AUTHORITY</OrganisationName>
      <Ident>
        <Id>1234</Id>
        <TypeId V="HER" S="2.16.578.1.12.4.1.1.9051" DN="HER-id"/>
      </Ident>
      <OrganisationName>Test</OrganisationName>
      <Ident>
        <Id>123456</Id><!-- 121017 (Test), 94378 (Prod) -->
        <TypeId V="HER" S="2.16.578.1.12.4.1.1.9051" DN="HER-id"/>
      </Ident>
    </Organisation>
  </Receiver>
</MsgHead>
```

Example of NPR export


Any user with a Norwegian license whose role includes the “Modality level” and “Export” functionality can access this feature in **Modality Level** → **+Add Export** → **Full Export**.



The necessary devices must be selected, and the report is generated when clicking on **Export**.




It consists of two files:




- Dose data file (RAD) that contains an internal PID, study and dose information
- Patient data file (IDE) that contains a link between the internal PID and the Social Security Number-

11. Location level

This dashboard can help users comparing different location performances. The user has to choose the location, then the modality they want to compare: by clicking the **Compare** button, which is situated in the upper right part of the page, the system will show a bar graph and a table.



The 3-character tag will displayed next to each site name, as configured in the section [Add Site](#).




The user can analyze the data in a similar way as in modality comparison.

For more information, refer to the video *How to use the Location Level* in our online training center.

12. Logbook

The logbook is available at institution, hospital and device level. It is the place where the user can register events that happen unexpectedly or are scheduled. The user can filter and export the logbook entries, attach documents and track everything that happens. This can be also used to set tasks with reminders and receive an email e.g. Event: unscheduled maintenance → Indicate from when to when the system was down and who needs to be notified.




The screenshot shows the Qaelum Logbook interface. At the top, there are navigation tabs: Start Portal, Study, Patient, Device, Modality, Location, Notifications, Logbook, COACH, Protocol Management, and Settings. The Logbook section is active, showing a summary of 132 entries and 42 unresolved entries. Below this, there are filter options for Created date, Reminder date, Deadline date, Description, Comment, and Collapse tree. A table of event types and their statuses is visible, including Preventive Maintenance, Unscheduled Maintenance, Task, Protocol task, Note, Error, Incident, and Other. The main table below has columns for Name, Creation date, Event type, Event status, Event effect, Reminder date, Deadline date, Event solution, Actions, and Attachments. The table contains several entries, including Preventive Maintenance for GE MEDICAL SYSTEMS BrightSpeed and Siemens Somatom Sensation 64 (BS).

Logbook at institution level

The filter can be used to find events of interest. Besides the main filter, additional filter fields can be activated: Created date, Reminder date, Deadline date, Description and Comment.

For example, users can search for events containing “Quality Control” in their description by activating the “Description” filter field, typing “Quality” and clicking on “Refresh”.




Search by "Description" of the event

For more information, refer to the video *How to use Logbook on Device and Hospital Levels* in our online training center.

12.1. Device level

Device specific logbooks can be accessed on **Device Level**, by selecting the device and clicking on **Logbook**.




Logbook on *Device Level*

Users can select between the *Tree table* and the *Calendar* views.




View selection in the upper right corner



Calendar view


To register a new event, click on **+Add entry** in the upper-right corner.

A new window will open where specific event details can be customized. The **Start event date** will be the used to place the event in the calendar, and will automatically take the creation date when the field is omitted by the user.




Create logbook entry

Attachments can be added (documents, PDFs, pictures, etc.) by dragging and dropping or selecting the file after clicking on the “...” button.




Attaching a file

After the upload is complete, it will be visible in the attachments list, where it can be downloaded or removed.




List of attachments

When an event has an attached file, a paperclip sign will be displayed under **Attachments..**




Paperclip icon indicating an attachment

The column can be added or hidden using the dropdown arrow:



Add "Attachments" column

Reminder dates and deadlines can also be indicated; they can be used to e-mail the involved users when the event date is approaching.




Add reminder date

E-mail alerts must first be configured in **Settings**, see [Configuration of e-mail sending in Logbook](#) for further details.

12.2. Hospital level

A site level logbook can be accessed on **Device** Level by clicking on the site and selecting the **Logbook** tab. All the site-specific logbook entries will be displayed here.




Logbook at hospital level

New logbook entries can be created here and will be automatically linked to the selected site.


12.3. Configuration of e-mail sending in Logbook

For reminder e-mails to be sent, a specific job must first be configured in **Settings -> Scheduler**. To configure a new job, access **Actions/Add job** at the bottom right corner, and after indicating a job name and description, select **Logbook follow-up reminder**:



Create a new scheduler job


E-mails will be sent to all involved users according to the configured settings. Users can choose to receive emails to notify all assigned users upon event creation (Creation Date), as reminders on the start date (Start Date), reminder date (Reminder Date) or deadline (Deadline Date). The configuration can be customized by checking the required boxes, as shown below:



The screenshot shows the 'Scheduler job' configuration window. It has two tabs: '1. General job info' and '2. Scheduled job action'. The 'Action' field is set to 'Logbook Date Checker'. The 'Interval' is set to '24 hours'. Below this, there are four checkboxes: 'Creation Date' (unchecked), 'Start Date' (checked), 'Reminder Date' (checked), and 'Deadline Date' (checked).

Job settings

As all scheduler jobs, they can be programmed to run with a certain frequency. Upon each trigger of the job, based on the set frequency, DOSE will check for any logbook entries with selected datetime within the configured time interval i.e. the last 5 minutes to the last 24 hours before the triggered job time.



The screenshot shows the 'Scheduler job' configuration window, specifically the 'Scheduled job action' tab. It features a navigation bar with tabs: 'Minutes', 'Hourly', 'Daily', 'Weekly', 'Monthly', 'Yearly', and 'Advanced'. The 'Daily' tab is selected. Below the tabs, there are two radio button options: 'Every 1 day(s)' (selected) and 'Every Week Day'. At the bottom, there is a 'Start time' field with dropdown menus set to '12' and '00'.

Set the job frequency

This can be configured by users with the functionality "Scheduler management" in their role.

13. Webservice integration

Webservice is part of DOSE by Qaelum’s Expert license: however, setting up these services requires third party support in order to configure the URLs. The configuration of the webservices is done upon request.

All Qaelum’s webservices requires the port 10525 to be opened in order to show results to final users. On request, Qaelum may set another port.

13.1. PACS: Patient Dose passport

This webservice is designed to show the patient’s dosimetric passport, as well as dose for the specific examination. There is a possibility to have a webservice called by patient ID that includes:

- total dose
- number of studies
- phantom with most irradiated areas
- list of studies performed.

Individual studies can be opened from the study history by clicking on *Show details*.




Webservice by patient ID

NOTE: those installations which include the license DOSE_PATIENT_REPORT can export the Japanese patient dose report directly from the webviewer, by clicking on the "Download" link, available in the study table (bottom of the webservice).



Download the Japanese patient dose report



このレポートは、放射線科の診断に使用されたCT検査の結果を示しています。患者の総有効線量（mSv）は10.7 mSvと算出されました。この値は、患者の放射線曝露の歴史に基づいて計算されたものです。また、このレポートには、検査の技術的パラメータ（534.8 mGy * cm）も記載されています。患者の個人情報（氏名、生年月日、ID）も確認できます。このレポートは、放射線科の診断に使用するための重要な情報源です。

検査ID	CT
検査日時	2022-05-18 00:00:00
検査日	2022-05-18
生年月日	1970-01-01
氏名	田中 太郎
ID	123456789

放射線科の診断に使用するための重要な情報源です。

Japanese patient dose report

There is also the possibility of calling a study by its StudyUID or a combination of Accession Number and Patient ID. This will result in a webservice showing the details of that specific study, including the compliance colour.

Category: Thorax-abdomen (Adult) DLP: 2.24E3 mGy * cm

DLP Limit(s): Show details

Upper acceptable: 960 mGy * cm

Lower acceptable: 0 mGy * cm

321

Patient ID

PACS6207374

Accession number

1

1

1

studies found for this patient

2.24E3


mGy * cm

dose of the selected study

41.27

mSv

total effective dose for this study



Webservice by Study UID or accession number



A "Comments" icon is present, and comment content is displayed when the user hovers over the icon.

1

Cause	Action	Justification
extravasation		


Study comments

If the icon is clicked, comments can be entered to be added to *Cause/Justification/Action* in the *Activity Stream*.



Add a comment


When notification comments are present the icon shows a warning (loudspeaker) sign.



Webservice by Study UID or accession number with study comments (above) and notification (below) comments

13.2. Proactive warning

This webservice provides some information that is useful when ordering or performing the exam, including but not limited to possible pregnancy risk, patient radiological history, patient cumulated dose and contrast related issues.




Proactive dose warning webservice

13.3. Clinical quality evaluation

This interactive webservice allows the radiologist to evaluate the quality of the images. The evaluations are based on a -2 to +2 scale and refer to several modality-specific quality indicators.

The webservice can be integrated into the radiologist’s PACS viewer, allowing users to score image quality alongside the images.

After clicking a link previously enabled link in the PACS interface, a window (webservice) with the layout shown below appears, allowing the user to click on the color scale and score the three criteria of image quality. A free text box for comments is also available.



Visualization of the webservice before and after image quality scoring

After submitting the score, the user will see a summary of the evaluation. This evaluation cannot be edited but a new evaluation can be added. This way multiple observers can evaluate the same study and the average scoring is automatically calculated.

Clinical Image Quality

Patient ID	60502848
Accession number	PACS6759293
Username	qaelum
Rating date	Dec 6, 2017 1:52:57 PM
Image contrast	2
Image sharpness	1
Patient positioning	-2
Comments	Image quality OK. Possibility to reduce the acquisition dose. Due to patient movement, the positioning is off.


Window with the scoring results as shown to the user

After submitting the score, the evaluation will be automatically added to the study details in DOSE. The user must click on **Show Clinical quality evaluations** to view it.

Study Details

Patient Name: Unknown

Overview
Series information
Organ doses
History
Activity stream



⚠⚠⚠⚠⚠

Study date
09/11/2015 10:10

Study description
CT thorax + abdomen

Protocol name
N/A

[Show study details](#)

[Show clinical quality evaluations](#)

[Show study comparisons](#)

Show clinical quality evaluations

Four tags will appear with average values from all submitted evaluations. The colors red/blue/green, respectively, indicate if the average value is unacceptable/acceptable/optimal.

Study Details
Patient Name: Unknown

Overview | Series information | Organ doses | History | Activity stream

Study date: 09/11/2015 10:10
Study description: CT thorax + abdomen
Protocol name: N/A

Patient information: 1009279 | F | 64 Y
Accession Number: 1234
Body part examined: CHEST

DLP [mGy.cm]: 1.90E3
Effective Dose ICRP103 [mSv]: 38.06
Event #: 4

Hide clinical quality evaluations | Show study comparisons

1 # quality evaluations (red) | -2 (-2, -1, 0, 1, 2) average contrast score (red) | 1 (-2, -1, 0, 1, 2) average sharpness score (blue) | 2 (-2, -1, 0, 1, 2) average patient positioning score (green)

Study with one quality evaluation

Study Details
Patient Name: Unknown

Overview | Series information | Organ doses | History | Activity stream

Study date: 30/03/2016 17:12
Study description: RAD ct thorax 28,[10-11-2-176]
Protocol name: G14_DE_Longenbool

Patient information: 72782212 | F | 50 Y
Accession Number: PACS6787663
Body part examined: CHEST

DLP [mGy.cm]: 177.97
Effective Dose ICRP103 [mSv]: 3.56
Event #: 5

Hide clinical quality evaluations | Show study comparisons

3 # quality evaluations (blue) | -0.7 (-2, -1, 0, 1, 2) average contrast score (blue) | 0.7 (-2, -1, 0, 1, 2) average sharpness score (blue) | 0.7 (-2, -1, 0, 1, 2) average patient positioning score (blue)

Study with three quality evaluations showing the average values

In the dosimetry overview of a device, the user can see all the added evaluations under the tab *Clinical Image Quality* in *Tools*.

Room 01 - Philips Ingenuity Core | CT
Patient radiation dose monitoring

1235 of 1236 studies visible | 1 of 1 notifications visible (All alerts) | Parameter type: DLP | mGy.cm

Compliance monitoring - Serie level | Statistic Spec. | Trends | Dynamic report | Trends Overview | SSDE | workload | Protocols | **Tools** | +

Reference dashboard overview | Study Overview | **Clinical Image Quality** | Alerts overview | Study comparison | Dynamic reports | Static reports | Batch actions | Configuration

Rating date	Study UID	Accession No	Patient ID	Image contrast	Image sharpness	Patient positioning	Comments
26/11/2015 15:25:11	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-1	2	-2	
26/11/2015 15:25:28	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	2	2	2	
26/11/2015 15:25:43	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-1	2	1	
26/11/2015 15:24:32	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-2	2	-2	
26/11/2015 15:24:52	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-2	-1	2	
28/01/2020 11:32:05	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-2	2	0	Contrast not acceptable

Details of every performed evaluation in the device, i.e., rating date, study UID, Acc. number, patient ID, image contrast, image sharpness, patient positioning and comments can be viewed here.

This data can be exported by clicking on the **Actions** button in the lower right corner.



Export to XLS

1	Rating date	Study UID	Accession No	Patient ID	Image contrast	Image sharpness	Patient positioning	Comments
2	26/11/2015 15:25:11	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-1	2	-2	
3	26/11/2015 15:25:28	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	2	2	2	
4	26/11/2015 15:25:43	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-1	2	1	
5	26/11/2015 15:24:32	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-2	2	-2	
6	26/11/2015 15:24:52	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-2	-1	2	
7	29/01/2020 11:32:05	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-2	2	0	Contrast not acceptable
8	26/11/2015 15:25:11	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-1	2	-2	
9	26/11/2015 15:25:28	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	2	2	2	
10	26/11/2015 15:25:43	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-1	2	1	
11	26/11/2015 15:24:32	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-2	2	-2	
12	26/11/2015 15:24:52	1.2.124.113532.80.22201.5366.20151022.162217.331974146	123	79283578	-2	-1	2	

Details as exported to XLS

An **Image quality** table is also available in the **Dashboard editor**.




Image quality table in the Dashboard editor

This table presents a summary of all performed image quality evaluations on the current device.

Compliance monitoring - Serie level Statistic Spec. Trends Dynamic report Trends Overview SSDE workload Protocols Image quality Tools

IMAGE QUALITY

Category	# Evaluation	Contrast	Sharpness	Patient Positioning
Not Acceptable	10	6	2	2
Sub Average	10	5	2	3
Average	20	7	6	7
Achievable	38	10	14	14
Optimal	33	9	13	11

*Image quality table within a **Device** level dashboard*

14. Settings


All the actions described here are available depending on the role of the user. So not everyone has access to see all the information. For further information consult section [User Management](#).

14.1. Warnings

14.1.1. Unconfigured data

When the system receives a study from an unknown device, the images are not processed until the device is configured. Users can find the list of images from unknown devices in this overview.

A user can add a new device copying its information to a site's device list . To copy the device's information the user must right click the row and then choose Copy to internal clipboard from the menu.



Warning overview copy to internal clipboard

After the information is copied, it is possible to copy it by clicking "Insert from clipboard" in the Add device page.

Add device ✕

+ Configure the parameters of this new device

Visible name *

Description

Collection ID *

Vendor

Model

Serial number

Internal reference

Modality *

Parameter

Station Name *

Institution Name

Data panels

✖ Cancel
📄 Insert from clipboard
+ Create device

Add device by inserting from clipboard

After the warning is reviewed, the user can remove it by right clicking on it and then choosing Remove Warning from the menu. By clicking on Remove similar warnings the user will delete all the warning concerning the same device. It is also possible to remove

Start Portal Patients Devices Modalities Settings

WARNING ?

Warning overview

DATE	MESSAGE
▲ 2014-03-05 12:22:41	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:22:41	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:21:21	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64LUZ_Gasthuisberg)
▲ 2014-03-05 12:21:21	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64LUZ_Gasthuisberg)
▲ 2014-03-05 12:21:21	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64LUZ_Gasthuisberg)
▲ 2014-03-05 12:21:21	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64LUZ_Gasthuisberg)
▲ 2014-03-05 12:20:47	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64JDD)
▲ 2014-03-05 12:20:47	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64JDD)
▲ 2014-03-05 12:19:07	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:18:49	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:18:49	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:18:47	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:18:03	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:18:03	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:18:02	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:18:02	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:17:52	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:17:52	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:17:50	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:17:50	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)
▲ 2014-03-05 12:17:38	Unknown system found during analysis (UJZLGRADISCT11 SIEMENS Sensation 64KIL_Gasthuisberg)

438 rows

Toggle Filter Bar visibility | Reset filter | Actions

all warnings.

Warning overview

All these tasks can be performed also by clicking on the Action button, which is located in the lower right part of this page.

14.1.2. Scheduler


The scheduler is a stand-alone application which is running next to the main web application. By using the scheduler, the user is able to schedule reoccurring jobs (e.g. reports, alerts calculation, DICOM Q/R ect.), which are going to be done automatically.

In the left top corner, the user can check the status of the scheduler (Running / Not running). The table shows scheduled jobs.

CREATION DATE	JOB NAME	ACTION NAME	USERNAME	FREQUENCY
2014-09-02 13:40:28	Monthly Dose Report	MailReportJob	qaelum	MONTHLY


Scheduler

The user can add, edit or trigger jobs by clicking on the Action button, which is located in the lower right part of this page. When adding a job, the user needs to configure the frequency the job needs to be triggered and an action. The frequency of the job can be configured in very detail. The user can select a frequency category or write one`s own regular expression



The screenshot shows the 'Scheduler job' configuration window with three tabs: '1. General job info', '2. Scheduled job action', and '3. Job frequency'. Under '2. Scheduled job action', the 'Action' is set to 'Limits calculation'. There is a checkbox for 'frequency is similar to selected interval'. The 'Interval' is set to '3 months' and the 'Grace period' is set to '1 hour'. Below this, there is a table for 'LIMITS' with columns 'DEVICE NAME' and 'LIMITS'. The 'DEVICE NAME' is 'Philips CT (SSDE_CT)'. The 'LIMITS' section has four checkboxes: 'Study Description', 'Series Combination', 'Protocol Name', and 'Study Groups', all of which are currently unchecked.

Scheduled job action



The screenshot shows the 'Scheduler job' configuration window with three tabs: '1. General job info', '2. Scheduled job action', and '3. Job frequency'. Under '3. Job frequency', there are tabs for 'Minutes', 'Hourly', 'Daily', 'Weekly', 'Monthly', 'Yearly', and 'Advanced'. The 'Monthly' tab is selected. Below the tabs, there are two radio buttons. The first is 'Day 1 of every 3 month(s)'. The second is 'The First Monday of every 1 month(s)'. At the bottom, there is a 'Start time' field set to '00:00'.

Job frequency of the scheduler job

For the SAR, organ dose and alerts calculations the user can choose the grace period. The grace is needed to avoid calculations on studies that are still not completely analyzed.

For the limits calculation, the user can also choose the time interval between the limits calculations.

14.2. Engine monitor

14.2.1. Engine Monitor

This tab shows the user general statistics for the system, such as the number of analyzed studies and the overall cumulative dose.

Engine Monitor

Total of analyzed studies : 83907


CT: 17016 DR: 2790 CR: 1899 DX: 1198 XA: 331 NM: 0

MG: 46476 US: 1333 CBCT: 2533 MR: 2616 RF: 7715

Engine motor: details of analyzed studies

14.2.2. Study analysis

This tool shows the longest analysis duration per study, filtered per modality.



Study Analysis

14.2.3. Log file browser

Here the user has a list of the latest log files created. They can click on them to see the details and download the file itself.

```

▶ TQM log files
tqm_2020_01_20.log
tqm_2020_01_19.log
tqm_2020_01_18.log
tqm_2020_01_17.log.121243179
tqm_2020_01_17.log.121003863
tqm_2020_01_17.log
tqm_2020_01_16.log
tqm_2020_01_15.log.125320052
tqm_2020_01_15.log
tqm_2020_01_14.log
tqm_2020_01_13.log
tqm_2020_01_12.log
tqm_2020_01_11.log
tqm_2020_01_10.log
tqm_2020_01_09.log.145209313
tqm_2020_01_09.log
tqm_2020_01_08.log
tqm_2020_01_07.log
tqm_2020_01_06.log
tqm_2020_01_05.log
tqm_2020_01_04.log
tqm_2020_01_03.log.153113566
tqm_2020_01_03.log
tqm_2020_01_02.log
tqm_2020_01_01.log

Session created: Mon Jan 20 09:35:02 CET 2020 ID: 104yp9qn9idg18gcu4ui2n9c3 Total Sessions: 1
Session destroyed: Mon Jan 20 09:35:32 CET 2020 ID: 104yp9qn9idg18gcu4ui2n9c3 Total Sessions: 0
Session created: Mon Jan 20 09:53:19 CET 2020 ID: dc0uwgouked1h00bjg2j3d7 Total Sessions: 1
Username: qaelum
Session destroyed: Mon Jan 20 09:53:32 CET 2020 ID: dc0uwgouked1h00bjg2j3d7 Total Sessions: 0
Session created: Mon Jan 20 09:53:32 CET 2020 ID: pdnwx22bysb912oolc6ayuvcj Total Sessions: 1
first query: 0
FunctionalityList: 31
Security framework of XStream not initialized, XStream is probably vulnerable.
Sites: 157
Roles: 0
Done: 0
Security framework of XStream not initialized, XStream is probably vulnerable.
== UPDATE DATASET START PORTAL ==
collectionId=DEVO
doseIndicator=DAP | dGy.cm²
startDate=Wed Jan 20 00:00:00 CET 2020
endDate=Mon Jan 20 23:59:59 CET 2020
SUMMARY SQL=SELECT ST.STUDY_UID, ST.STUDY_ID, ST.STUDY_DATETIME, ST.ACCESSION_NUMBER, ST.BODYPART_EXAMINED, ST.COLLECTIONID,
ST.EVENT_COUNT, ST.MANUFACTURER, ST.MANUFACTURER_MODELNAME, ST.SOFTWARE_VERSION, ST.ID, ST.MODALITY, ST.OPERATOR_NAME, ST.IGNORE,
ST.PATIENT_AGE, ST.PATIENT_AGE_TYPE, ST.PATIENT_GENDER, ST.PATIENT_DOB, ST.PATIENT_WEIGHT, ST.PATIENT_LENGTH, ST.PATIENT_ID, ST.SITEID,
ST.DBO_DECRYPTPATIENTNAME(ST.PATIENT_NAME ?) AS patientname, ST.PATIENT_AGE_DAYS, ST.PERFORMING_PHYSICIAN_NAME, ST.REQUESTING_SERVICE,
ST.PROTOCOL_NAME, ST.REQUESTING_PHYSICIAN_NAME, ST.STUDY_DESCRIPTION, ST.SERIES_COMBINATION, ST.STUDY_GROUP_DEFINED,
ST.ENTRANCE_DOSE, ST.FLUORODAPTOTAL, ST.FLUOROTIMETOTAL, ST.DAPTOTAL, ST.PSD, ST.VENDOR_MAX_SKIN_DOSE, ST.LEGAL_CATEGORY,
ST.MEDICAL_ALERTS, ST.WARD, ST.PSD_PHANTOM_TYPE, ST.PSD_TARGET_REGION, ST.PSD_PATIENT_WEIGHT, ST.PSD_PATIENT_LENGTH,
ST.PSD_PATIENT_AGE, ST.DOSERTOTAL, ST.IMAGE_AMOUNT, ST.PERF_STUDY_START_DATETIME, ST.PERF_STUDY_END_DATETIME, ST.EFF_DOSE_ICRP103,
ST.DAPTOTAL FROM SUMMARY_DOSE_GEN_RAD ST WHERE ST.COLLECTIONID=? AND ST.STUDY_DATETIME=? AND ST.STUDY_DATETIME=? AND
(ST.IGNORE IS NULL OR ST.IGNORE=? OR ST.IGNORE=?2)
Loaded alerts in 273ms.
=== SUMMARY FINISHED IN 445 Ms ===
== UPDATE DATASET START PORTAL ==
collectionId=Siemens_Fluorospot_Compact_FD
doseIndicator=DAP | dGy.cm²
startDate=Wed Jan 20 00:00:00 CET 2020
endDate=Mon Jan 20 23:59:59 CET 2020
SUMMARY SQL=SELECT ST.STUDY_UID, ST.STUDY_ID, ST.STUDY_DATETIME, ST.ACCESSION_NUMBER, ST.BODYPART_EXAMINED, ST.COLLECTIONID,
ST.EVENT_COUNT, ST.MANUFACTURER, ST.MANUFACTURER_MODELNAME, ST.SOFTWARE_VERSION, ST.ID, ST.MODALITY, ST.OPERATOR_NAME, ST.IGNORE,
ST.PATIENT_AGE, ST.PATIENT_AGE_TYPE, ST.PATIENT_GENDER, ST.PATIENT_DOB, ST.PATIENT_WEIGHT, ST.PATIENT_LENGTH, ST.PATIENT_ID, ST.SITEID,
ST.DBO_DECRYPTPATIENTNAME(ST.PATIENT_NAME ?) AS patientname, ST.PATIENT_AGE_DAYS, ST.PERFORMING_PHYSICIAN_NAME, ST.REQUESTING_SERVICE,
ST.PROTOCOL_NAME, ST.REQUESTING_PHYSICIAN_NAME, ST.STUDY_DESCRIPTION, ST.SERIES_COMBINATION, ST.STUDY_GROUP_DEFINED,
ST.ENTRANCE_DOSE, ST.FLUORODAPTOTAL, ST.FLUOROTIMETOTAL, ST.DAPTOTAL, ST.PSD, ST.VENDOR_MAX_SKIN_DOSE, ST.LEGAL_CATEGORY,
ST.MEDICAL_ALERTS, ST.WARD, ST.PSD_PHANTOM_TYPE, ST.PSD_TARGET_REGION, ST.PSD_PATIENT_WEIGHT, ST.PSD_PATIENT_LENGTH,
ST.PSD_PATIENT_AGE, ST.DOSERTOTAL, ST.IMAGE_AMOUNT, ST.PERF_STUDY_START_DATETIME, ST.PERF_STUDY_END_DATETIME, ST.EFF_DOSE_ICRP103,
ST.DAPTOTAL FROM SUMMARY_DOSE_GEN_RAD ST WHERE ST.COLLECTIONID=? AND ST.STUDY_DATETIME=? AND ST.STUDY_DATETIME=? AND
(ST.IGNORE IS NULL OR ST.IGNORE=? OR ST.IGNORE=?2)

```

Log file browser

14.3. Data tables

14.3.1. Conversion factors

In DOSE, effective dose is calculated for all ionizing radiation modalities. The calculation is based on conversion factors that depend on modality, age group and body part. The effective dose is calculated as the product of the dosimetric parameter of every modality (DAP, DLP, etc) with the age & body part-specific conversion factor.

$$EffDose = DoseParameter * ConversionFactor(age, bodypart)$$

For Mammography, since the scattered radiation on other organs is not available, effective dose is calculated based on the glandular dose. For Nuclear Medicine, the conversion factors depend, besides the radiopharmaceutical, on the age and study description or condition (see [Dosimetry](#)).

This table shows all the conversion factors calculated by Qaelum for the calculation of effective dose ICRP60 and ICRP103, grouped by age category, body part and modality.

Conversion Factors

General Conversion factors Nuclear Medicine Conversion factors

Qaelum | + Add Remove

Start Age	End Age	Body part	Modality	Conv ICRP60	Conv ICRP103
0	0	ABDOMEN	CBCT	0.241	0.248
1	3	ABDOMEN	CBCT	0.117	0.115
4	7	ABDOMEN	CBCT	0.075	0.072
8	12	ABDOMEN	CBCT	0.0411	0.0453
13	17	ABDOMEN	CBCT	0.0259	0.0248
18	150	ABDOMEN	CBCT	0.0211	0.0202
0	0	ABDOMENPELVIS	CBCT	0.241	0.248
1	3	ABDOMENPELVIS	CBCT	0.117	0.115
4	7	ABDOMENPELVIS	CBCT	0.075	0.072
8	12	ABDOMENPELVIS	CBCT	0.0411	0.0453
13	17	ABDOMENPELVIS	CBCT	0.0259	0.0248
18	150	ABDOMENPELVIS	CBCT	0.0211	0.0202
0	0	ANKLE	CBCT	0.0364	0.0299
1	3	ANKLE	CBCT	0.01	0.0084

1445 rows

CT = mSv/mGy*cm GenRad= mSv/dGy*cm² Mammo= mSv/mGy NM= mSv/MBq

Conversion factors by Qaelum

The conversion factors by Qaelum are derived by dosimetry tools available for patients of all different age groups (newborn, 1-year-old, 5-year-old, 10-year-old, 15-year-old and adult) and for the typical settings of an examination. For Nuclear Medicine, they are based on published data from the ICRP reports.

The Qaelum CT conversion factors are based on simulations performed with the ImpactScan dosimetry tool for the standard settings of a typical examination. The pediatric CT conversion factors are based on the ones of adults, multiplied by an age/body part-dependent factor.

Note: In the case of head conversion factors for CT, the head phantom is considered, while in the case of neck examinations, the body phantom is considered. In case that a hospital uses the head phantom for neck protocols, it is suggested to replace the CT conversion factor of the neck with the following:

Bodypart	Adult	15 yo	10 yo	5 yo	1 yo	0 yo
NECK	0.0055	0.0061	0.0069	0.0091	0.0121	0.0135

As it is not straightforward to perform simulations on all different body areas, some body parts follow the conversion factors of the closest common body area.

Example cases for CT:

- Coccyx, Lspine and Hip follow Pelvis
- Tspine and Thorax follow Chest
- Clavicle follows Shoulder
- Humerus, Scapula, Heart follow Chest
- Skull, Brain, Mouth, Jaw follow Head

Example cases for DR/RX/DX/CR/RF/XA:

- Coccyx and Lspine follow Pelvis
- Shoulder, Cspine, Tspine and Thorax follow Chest
- Abdomen, Pelvis, Abdomenpelvis have the same factors
- Spine is the average of Cspine - Tspine - Lspine

Conversion factors for CBCT currently follow projection radiography.


14.3.1.1. CONVERSION FACTOR TEMPLATES

Sometimes there is a need to use different sets of conversion factors for the effective dose for different sites. Those sets, which are called Templates in DOSE, can be created, edited and assigned to sites and devices from the web application by any user with “Physical Parameters Management” functionality in their role.

Each action which is performed at Conversion Factors level is logged for auditing purposes.

To create a template:

1. Go to Settings → Conversion factors.
2. Click on Add.



Conversion Factors

General Conversion factors | Nuclear Medicine Conversion factors

Qaelum + Add Remove

Start Age	End Age	Body part	Modality	Conv ICRP60	Conv ICRP103
0	0	ABDOMEN	CBCT	0.241	0.248
1	3	ABDOMEN	CBCT	0.117	0.115
4	7	ABDOMEN	CBCT	0.075	0.072
8	12	ABDOMEN	CBCT	0.0411	0.0453
13	17	ABDOMEN	CBCT	0.0259	0.0248
18	150	ABDOMEN	CBCT	0.0211	0.0202
0	0	ABDOMENPELVIS	CBCT	0.241	0.248
1	3	ABDOMENPELVIS	CBCT	0.117	0.115
4	7	ABDOMENPELVIS	CBCT	0.075	0.072
8	12	ABDOMENPELVIS	CBCT	0.0411	0.0453
13	17	ABDOMENPELVIS	CBCT	0.0259	0.0248
18	150	ABDOMENPELVIS	CBCT	0.0211	0.0202
0	0	ANKLE	CBCT	0.0364	0.0299
1	3	ANKLE	CBCT	0.01	0.0084

1445 rows

CT = mSv/mGy*cm GenRad = mSv/dGy*cm² Mamma = mSv/mGy NM = mSv/MBq

Toggle Filter Bar visibility | Reset filter | Actions

3. Insert the name of the template, then click on Commit.

Conversion Template ✕

Creating a template and assigning it to one or multiple devices introduces a risk on effective doses. If conversion factors are not correct, the effective dose of the studies will be incorrect too. If the template doesn't contain all the bodyparts, the effective dose of the studies may be incomplete or missing. Qaelum cannot be held responsible for any issues caused by wrong or incomplete templates.

Name

+ Commit
⊘ Cancel

Templates are created as duplicates of the predefined template provided by Qaelum, so to prevent users from having to recreate all cases.

To edit a template:

1. Go to Settings → Conversion factors.
2. Select a template.

Conversion Factors

General Conversion factors
Nuclear Medicine Conversion factors

New template

+ Add
🗑 Remove

	End Age	Body part
0	0	ABDOMEN
1	3	ABDOMEN
4	7	ABDOMEN
8	12	ABDOMEN
13	17	ABDOMEN

3. Click on Actions → Insert data.



4. Once the edit is complete, click on Commit.

To delete a template:

1. Go to Settings → Conversion factors.
2. Select a template.
3. Click on Remove.




4. Click on OK.

Users are only allowed to remove templates which are not assigned to sites or to devices.

To assign a template to a site and all the devices belonging to the site:

1. Go to Device → Site overview.
2. Click on Configuration, then on Effective dose.
3. In the dropdown Template, select a template.




4. Click on the Save button.

5. Click on Commit.

Every newly created device is assigned to the Template of the site to which the device belongs to.

To assign a template to a device:

1. Go to Device → Device overview.
2. Click on Configuration → Bodypart mapping and effective dose → Device configuration.



3. In the dropdown Template, select a template.
4. Click on the Save button.

14.3.2. Study Groups


With this tool the user can create, edit and delete population groups (based on modality, age range and gender), for which it will be possible to configure acceptable and achievable dose limits. Afterwards those study groups can be used on device level for the export and filter function. For more information, see [Compliance Configuration](#).

NAME	MODALITY	START AGE	END AGE	GENDER
Abdomen	CT	18	150	MFO
Abdomen CHILD	CT	0	17	MFO
Cervical spine ADULT	CT	19	150	MFO
Lumbar spine	CT	18	150	MFO
Thorax	CT	18	150	MFO

Study groups overview

14.3.3. Dicom bodypart mapping

Here the user can manage the bodypart mapping described in [Mapping at series level](#) in a centralized way. Select the device in the dropdown menu in the upper-right corner.



Dicom bodypart mapping

14.4. Data management


14.4.1. Remove data

When a site or a device is deleted, studies data is not removed. To delete dose and contrast data, the administrator must use this tab. By clicking on a device, the user can decide to erase contrast data or all data for it.

Both methods will erase the data without leaving any revision or retrace possibility.

14.4.2. Analysis configuration

Here some configurations for the analysis are enabled/disabled.



Analysis configuration

14.4.3. Notification Recalculation Settings

The default notification recalculation settings can be configured here. The selected configuration will be displayed in the Notification Recalculation window (see section [Compliance Configuration](#)).

Notification Recalculation Settings

Notifications with the following status will be deleted and therefore completely recalculated.
 Notifications with unselected status will be archived

	Name
<input checked="" type="checkbox"/>	Unknown
<input checked="" type="checkbox"/>	Unresolved
<input type="checkbox"/>	Under investigation
<input type="checkbox"/>	Pending
<input type="checkbox"/>	Resolved
<input type="checkbox"/>	Obsolete

Default settings for recalculation

Here the user can also configure what will happen to the information within archived notifications i.e. notifications that have been recalculated and therefore overwritten.

Actions to perform on the newly generated notifications

Severity	Collect status of archived notifications	Collect data of archived notifications	Keep archived notifications
Equal severity	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Higher severity	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Lower severity	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

In the table above, three types of severity for the newly calculated notifications are shown:

- Equal severity: the new notification has the same severity (color) as the old notification (e.g. ORANGE -> ORANGE, RED -> RED)
- Higher severity: the new notification is more severe than the old notification (e.g. RED -> ORANGE, ORANGE -> GREEN)
- Lower severity: the new notification is less severe than the old notification (e.g. ORANGE -> RED, GREEN -> ORANGE)

For these three cases the user can configure what happens to the information within old notifications as follows:

- Collect the status of archived notifications: e.g. if the old notification was “resolved”, the new notification will be automatically marked also as “resolved”.
- Collect the data of archived notifications: if anything was added in the old notification (e.g. comments in “Handle”), that data will be transferred to the new notification (except its status).


- Keep archived notifications: both old and new notifications will be kept and shown.

If the archived notifications are kept, they can be visualized in the Notifications tab by checking the **Include archived** box located above the Notifications filter.



Include archived notifications

Archived notifications can be visualized in the Study Details by checking the **Show archived notifications** box.




Show archived notifications in Study Details

14.5. Tools

14.5.1. User Action log

The User Action log allows the possibility to track specific actions user per user within a certain date range. Click on **Query**, select a user to be tracked, the specific action and the date range.



User Action Log

Query Reset

Id	Username	Time	User Action	Details
37640	qaelum	17/01/2020 10:26:08	Body part changed	Changed bodypart to CHEST for all studies with a series with series description=N/A for device GRADTOCT14-PatientPositioning
36758	qaelum	06/01/2020 10:39:12	Body part changed	Changed bodypart to ABDOMEN for study 1.2.124.113532.80.22201.5366.20131231.155111.1153419054
36354	qaelum	17/12/2019 09:26:56	Body part changed	Changed bodypart to CHEST for study Philips.Patient.Positioning.Abdomen.001 for the series with instance uid 1.3.46.670589.33.
36351	qaelum	17/12/2019 09:26:03	Body part changed	Changed bodypart to ACA for all studies with a series with series description=TEST for device GRADPHCT12-PatientPositioning
36346	qaelum	17/12/2019 08:41:00	Body part changed	Changed bodypart to THORAX for all studies with a series with series description=TEST for device GRADPHCT12-PatientPositioning
36337	qaelum	17/12/2019 08:36:27	Body part changed	Changed bodypart to CHEST for study Philips.Patient.Positioning.Abdomen.001
36334	qaelum	17/12/2019 08:35:28	Body part changed	Changed bodypart to THORAX for study Philips.Patient.Positioning.Abdomen.001

User action log

14.5.2. DICOM Tools

This tab contains the query / retrieve (Q/R) tools.

Query / Retrieve studies

Clear results | Reset configuration | C-ECHO Query Retrieve Query / Retrieve

Results Configuration

Use Save

Our Port: 104

Our AE Title: QAELUMQR

Our store AE Title: QAELUM

Their IP / Host: 172.16.0.25

Their Port: 11112

Their AE Title: DCMCHCEE

Dicomstore: C:\qaelum\dicomstore

Query level: STUDY

■ Patient related parameters

Patient Name: _____ Patient ID: QAELUM1 Issuer of Patient ID: _____ Patient birth date: _____

Patient Sex: _____

■ Study related parameters

Study date: _____ Study time: _____ Accession number: _____ Modalities in study: _____

Referring physician name: _____ Study ID: _____ Study UID: _____

■ Series related parameters

Modality: _____ Series Number: _____

■ Instance related parameters

SOP Class UID: 1.2.840.10008.5.1.4.1.1.88.67 SOP Instance UID: _____

Query and retrieve studies

The first line is for saved configurations. By clicking **Save**, the current configuration will be saved so that you will be able to reload it by clicking Use.

Our Port: 104

Our AE Title: QAELUMQR

Our store AE Title: QAELUM

Their IP / Host: 172.16.0.25

Their Port: 11112

Their AE Title: DCMCHCEE

Dicomstore: C:\qaelum\dicomstore

Query level: STUDY


All general info about the DOSE server and PACS

This section contains all information about the DOSE server and the PACS to be added to Q/R. Is it possible to query on Patient, Study, Series or Instance levels. The parameters for such levels are shown in the following sections.

It is possible to run only queries, to retrieve the studies based on query selection or to query / retrieve.

14.5.3. Notification configuration

The goal is to notify users when there is a specific problem on a study. It can be done by using predefined scripts or selecting an alert level for a specific parameter. The study's where the users want to have a notification of can be on site, modality or equipment level. Furthermore, users can make an action object whereby they can select which people will get a notification and which type of message.




Notification configuration

The steps for this procedure are:

1. Create an email notification **action**
2. Create a **configuration** (alert) and link it to the action.


14.5.3.1. CREATE AN E-MAIL NOTIFICATION ACTION


When logged in in DOSE, please select Settings, then Notification configuration. There, please select the second tab (Actions).




By clicking on Actions in the lower right side of the page, a menu will appear, select Add configuration.

The popup will allow you to create an action. In this context, an action is a group of users which need to be notified when an alert is created.

Each action requires a name and at least one user or Team. Please type the name in the correspondent field and select the users by clicking on their name and then clicking on the  symbol. As this happens, the users or teams are moved to the right table.




You can unselect a user by clicking on the name in the right table and then on .


To save, please click on **Commit**.

14.5.3.2. CREATE AN E-MAIL ALERT

E-mail alerts can be configured for exceeding static and/or dynamic limits. To add a new e-mail alert access the **Configuration** tab, click on the **Actions** button in the lower-right corner and select **Add configuration**.



The following window will be displayed:



And the following alert fields can be configured:

- alert name
- parameter type
- severity level to be alerted of (orange/red)
- limit type (static/dynamic)
- filter by hospital, modality or device to be alerted for

An action (e-mail sending) can be linked to this alert. In the above example, the user Qaelum will receive an e-mail for each time that the DLP from any CT of any site is above the static acceptable (red) level.

Remark: each alert level should be configured separately. For example, configuring an orange alert level for DLP will not automatically trigger the red alert level for this parameter. For more information, refer to the video *How to configure e-mail sending for compliance notifications* in our online training center.

15. Integrations

Several integrations, described below, are available.

- **DICOM.** This integration is mandatory, as it is needed in order to collect data. DOSE contains a DICOM SCP for C-STORE, to which the PACS and/or the modalities can send the studies. For a description of the supported DICOM objects, please refer to the DICOM Conformance Statement.
- **MPPS.** DOSE can also contain an MPPS receiver. The messages are analysed and collected similarly to other DICOM objects.
- **HL7.** DOSE uses a customized installation of Mirth connect in order to collect and send data from/to third party via HL7. The following kind of messages are supported:
 - **Inbound:**
 - **ADT.** This kind of messages is used to edit patient data.
 - **ORM.** This kind of message is used to create studies for which no DICOM data is available, or to add some details.
 - **Outbound:**
 - **ORM.** This kind of messages is used to export dose data.

Additional details can be found in the HL7 Conformance Statement.

- **Webservice.** DOSE can be integrated in third party software via webservice. This allows showing information like patient passport directly on the third party.


16. About

All information about the system and its producer are shown in the lower right corner of every page.

By clicking on the **About** button, the user can find information about the project ID, the software version, the company ID and contacts, and a link to this manual. It is also possible to report a bug and/or give feedback.

© 2012-2022 QAELUM NV. All rights reserved Licensed to Development | [About](#)



Contents of the **About** button



Project ID
Customer project ID

Software identification
DOSE - Patient radiation dose monitoring software
VERSION: v24

Please read the instructions for use carefully before using this application.

 [User manual](#)
 [Technical requirements](#)


Link to Online Instructions for Use
<https://qaelum.com/media-center/documentation>


To access our on-line Training Center and watch educational videos, please visit:
knowledge.qaelum.com


To report a bug or to give feedback, please use our ServiceDesk:
[Qaelum ServiceDesk](#)
support@qaelum.com


CE 1639

[Cancel](#)

 Qaelum NV
Kolonel Begaultlaan 1B
3012 Leuven (Wilsele)
Belgium

 2024-03-25
(YYYY-MM-DD)

 HRI =
(01)05430003132
012(8012)v24




The most up to date external user manual can be accessed using the [Qaelum Disk Station](#) website link provided in the **About** box (see above), using the following credentials, given in the F20 Installation and Configuration report (chapter 3 Installation and configuration checklist):

Username: QaelumCustomer


Password: QaelumCustomer01

The user can log in using the screen below:



Log-in screen to access the user manual

The latest applicable version of the user manual can be downloaded by clicking on the "Download".



User manual download

16.1. Online Training Center

Qaelum's online training center is a tool that aims to support DOSE users once the personalized training is finished. This center contains explanatory videos that show all the steps necessary to obtain the desired results in the application. They are sorted by

theme, with titles in the format "How to configure...", "How to use...?" and includes a search tool.

16.1.1. Sign in


To access the training center, the following link can be used:

knowledge.qaelum.com

Registration can be completed via the following link:

<https://qaelum.com/training-registration>


Following registration, Qaelum will activate the user accounts and a confirmation e-mail will be set. Once activated, the training center can be entered using your credentials.



The screenshot shows a login form titled "Log in". Below the title is the instruction "Log in to access your restricted pages." The form contains two input fields: "Username *" with the text "@qaelum.com" and "Password *" with masked characters ".....". There is a checkbox labeled "Remember me" and an orange "LOG IN" button. At the bottom, there are two links: "Forgot your password?" and "Forgot your username?".


Username and/or password can be changed by clicking on ***Forgot your password/username?***

Once logged in, the ***Knowledge*** section will appear on the left-hand side of Qaelum's website menu. Once accessed, users will be redirected to the **DOSE training** section, where the training center is located.



16.1.2. How to find a video

In the training center, the videos are grouped according to the user level (basic and advanced) and by topic. Each video is accompanied by a brief description. Click on **Watch** to access a video. Additionally, a search tool is available in the upper right corner. Keywords can be entered here to search for the desired videos according to video title and caption.




16.1.2.1. MENU BY TOPIC

To find a video, you can search the menu by training level (basic or advanced) and within the level, by topic. For example, to find a basic training video related to Device Level, the **DOSE Basic Training/Device level** section is accessed:

<p>DOSE Basic Training</p>	<p>How to use the device tree on Device Level</p>
<p>- Start Portal</p>	<p>How to configure, expand/collapse, group and filter the device tree on Device Level.</p>
<p>- Study Level</p>	<p>WATCH</p>
<p>- Patient Level</p>	<p>How to use the Live Dashboards</p>
<p>- Device Level</p>	<p>How to access and use the Live Dashboards on Device Level.</p>
<p>- Modality Level</p>	<p>WATCH</p>
<p>- Location Level</p>	<p>How to use the Live Dashboard Modality function</p>
<p>- Notifications</p>	<p>How to use the Live Dashboard Modality feature that allows to see the latest studies from multiple devices of the same modality, configure the number of studies, refresh interval, filter, and export data.</p>
<p>- Logbook</p>	<p>WATCH</p>
<p>DOSE Advanced Training</p>	
<p>- Device Level</p>	
<p>- Modality Level</p>	
<p>- Mapping to DICOM Standard</p>	


When the relevant video is located, click on *Watch*.

<p>DOSE Basic Training</p>	<p>How to use the device tree on Device Level</p>
<p>- Start Portal</p>	<p>How to configure, expand/collapse, group and filter the device tree on Device Level.</p>
<p>- Study Level</p>	
<p>- Patient Level</p>	
<p>- Device Level</p>	
<p>- Modality Level</p>	
<p>- Location Level</p>	
<p>- Notifications</p>	
<p>- Logbook</p>	
<p>DOSE Advanced Training</p>	
<p>- Device Level</p>	
<p>- Modality Level</p>	


16.1.2.2. SEARCH TOOL

A keyword search engine can be found in the upper right corner. This will return results based on video title and description.


For example, to find a video related to patient information, "patient" can be entered in the search engine:



Returning a list of results:



Clicking on any of the results will navigate directly to the video page.




The screenshot shows a training interface. On the left is a navigation menu with sections: 'DOSE Basic Training' (containing Start Portal, Study Level, Patient Level, Device Level, Modality Level, Location Level, Notifications, Logbook) and 'DOSE Advanced Training' (containing Device Level). The 'Patient Level' item is highlighted. The main content area features a video player with the title 'How to use Patient Level' and a description: 'How to search for patients, view patient history, and access the Patient Dose Passport.' The video player shows a play button, a progress bar at 00:51, and the Vimeo logo.

16.2. Service desk

16.2.1. Invitation for Service Desk

An invitation to the Service Desk is sent to the contact email provided to Qaelum.



The screenshot shows an email interface with the subject 'Welcome to Qaelum ServiceDesk'. The sender is 'Qaelum ServiceDesk <jira@qaelum.atlassian.net>'. The body of the email says: 'Hi, Filippo Miniati has invited you to the Qaelum ServiceDesk portal! Visit the portal to raise requests and get help.' The Qaelum logo is at the bottom, along with the text 'Powered by JIRA Service Desk'.


Qaelum Service Desk

16.2.2. Configure your login details

When clicking on the link “Visit the portal”, you are directed to the configuration window for setting up your log-in information.

Your email is automatically generated as your username.


Provide your full name and the password you would like to use. In case you would like to see the password, click on the check box “Show password”.



16.2.3. Reporting via Service Desk

In order to create your ticket, you can use the following URL:


<https://qaelum.atlassian.net/servicedesk/customer/portal/2>



A direct link to our Service Desk can be found in the **About** box of all Qaelum products.


Log in with your credentials to start using the Service Desk.

After log-in you are directed to the start portal of the Service Desk. This page contains links to some FAQs and basic troubleshooting that can answer the common questions or complaints. We would appreciate if you first check to see if your complaint is related to a document that can clarify and possibly solve your issue immediately. Otherwise feel free to raise a request.



FAQs in the Service Desk

Please select the type of issue that you would like to report to us.




Issue types in the Service Desk

16.2.3.1. COMPLAINT REPORTING

Reporting a complaint requires some mandatory fields to be filled in. In the “Product type” dropdown please indicate which product the issue is related to. Identify the problem with

an accurate title and add it to the “Summary” box. Provide us with a detailed description of the problem in the “Description” box. In the “Affects version”, please specify which version of our product your complaint is about. In case you would like to share a screenshot or other files, you can upload it in the “Attachment” box. Finally, please indicate your customer reference number, which is provided to you by our installation team.




Qaelum HelpCenter / Qaelum ServiceDesk / Raise a request

Qaelum ServiceDesk

Welcome to Qaelum ServiceDesk. You're invited to give us feedback and to report issues.

What can we help you with?

 Complaint

Product Type *

DOSE

Please select the Qaelum product for which your request is addressed.

Customer Reference number

DEMO

Please fill in your Customer Reference, as noted in documents

Summary *

Coffee not warm enough

Please provide a short title describing your request.

Description *

In the morning the heating takes too long and results are not as expected.

Please provide further detailed information on your request: E.g. Which is the problem that you encountered? Which is the page / dashboard you were using? What were you doing or trying to do?

Affects versions

v19.11

Please provide the software version in which you experience this issue

Attachment

Drag and drop files, paste screenshots, or browse

Browse

If you prefer, you can also add a file or an image

Send Cancel

Creation of a Service Desk complaint

To send your issue, click “Create”.


After creating the issue, you can view the issue summary. You can add extra comments or attachments.

Qaelum HelpCenter / Qaelum ServiceDesk / TCP-2019
Coffee not warm enough

Add a comment TO DO

Activity

Test Test Today 5:07 PM LATEST



Details

Created at
Today 5:07 PM

Product Type
DOSE

Customer Reference number
DEMO

Description
In the morning, the heating takes ages and results are not as expected

Don't notify me


Shared with
Test Test
Creator

Service Desk complaint

Other members of your team will be automatically added as viewers of the ticket. Please let us know if you would like to limit the viewers.

An email will be sent to you, notifying that your ticket was received.

TCP-2019 Coffee not warm enough Inbox x

 **Qaelum ServiceDesk** <servicedesk@qaelum.atlassian.net>
to me

Reply above this line:

Just confirming that we got your request. We're on it.

[View request](#) - [Turn off this request's notifications](#)


This is shared with Test Test

Powered by JIRA Service Desk

Email confirmation of Service Desk ticket receipt

Every new change, modification, development, and update of your issue triggers automatic email notifications. This way every step in the handling of your issue can be tracked.

Throughout the process of handling the issue, our Service Desk agents will provide you with feedback, and if needed your confirmation or input may be asked. We kindly ask you to reply to these requests so that we can continue handling your complaint in a timely fashion according to our SLAs. You can reply via the Service Desk webpage or by simply replying to the email.



Email from the Service Desk

After the resolution has been deployed, your approval will be asked. You can approve the resolution by simply clicking on "Accept" in the Service Desk.

Qaelum HelpCenter / Qaelum ServiceDesk / TCP-2019
Coffee not warm enough

Your approval

Approve Decline

Add a comment

Activity

- Automatic response** Today 5:24 PM **LATEST**
Request requires approval. 1 approval needed.
- Automatic response** Today 5:24 PM
Your request status changed to Resolution Approval Process.
- Filippo Miniati** Today 5:24 PM
Dear,
can you confirm the issue is solved now?
Thanks.

RESOLUTION APPROVAL PROCESS

Don't notify me

Shared with

- Test Test
Creator

Resolution Approval Process

1 approval needed

- Test Test

Resolution approval

Automatic resolution acceptance will occur after a time period of 7 working days following the deployment of the resolution.

After the issue is considered resolved, you may receive an additional email, allowing you to evaluate Qaelum’s customer care.

Filippo Miniati <servicedesk@qaelum.atlassian.net>
to me

Reply above this line.

Filippo Miniati resolved this as Fixed.

How was our service for this request?

★ ★ ★ ★ ★
Very poor Poor Neither good nor poor Good Very good

16.2.3.2. OTHER REQUESTS

Similarly to complaint reporting, you can also report the requests listed below.

16.2.3.3. NEW INSTALLATION OR CONFIGURATION REQUEST

Please use this request type for adding new modules, devices / modalities or hospitals to your existing installation.

16.2.3.4. INFORMATION OR APPLICATION REQUEST


Sometimes you may simply need an explanation of a feature: in this case, or when you need a training, this request type can be raised.

16.2.3.5. NEW FEATURE REQUEST

Please use this kind of request to suggest new tools/features. We will consider them for the following releases!

16.2.4. All requests overview

When you log in to Service Desk, in the top right corner of the start portal, you can see "Requests" with a number indicating the number of open requests.




Requests in Service Desk

When you click on it, you receive an overview of all requests you and your team have sent to us. In addition, you can also search the requests and filter by all, open or closed requests.







16.2.5. Profile

In the top right corner of the start portal, you can find an icon for your profile. By clicking on your profile, you can personalize or change the information in your profile. You can change the icon that appears, your password, name, and time zone.



Profile details in the Service Desk

16.3. Symbols

Symbol	Definition
	This symbol indicates to read the instructions for use and links directly to the user guide.
	Symbol for Manufacturer
	Symbol for Date of manufacture
	Symbol for UDI Carrier
	Symbol for Manufacturer's Catalog Number
	Symbol for Medical Device

17. Supported browsers

DOSE is a zero-footprint software web application. This means that it requires no client-side application install and has no OS dependency. The application can be accessed from anywhere within the hospital network by browsing to the URL set up for DOSE. We support the following browsers:

- ✓ Google Chrome 80+
- ✓ Mozilla Firefox 70+
- ✓ Internet Explorer 11
- ✓ Microsoft Edge

Please see the Technical Requirements document for more information

18. Minimal hardware requirements

These requirements are only applicable to small installations (10 devices / 50.000 studies per year) when redundancy is not required. For further details on bigger installations, please see the Technical Requirements document or ask Qaelum for the Project Proposal.

RAM memory	12 GB
CPU	4 CORES +2GHz
Hard Disk	60GB for OS & application[C:/] 150GB data drive [e.g. D:/] 50GB cache drive[e.g. E:/]
Network connection	LAN and static IP address
Operating System	Windows Server 2012 R2 or newer

19. Functional description

19.1. Background

A quantity of ionizing radiation is related to each examination performed in radiology. Even today, when most systems are working according to the ALARA principle (“as low as reasonably achievable”), it is important to have an overview for all the study types and given doses per X-ray device and the received doses over time for each patient. For X-ray devices, the monitoring of typical dose values is needed for optimization reasons. For patients, the tracking of these doses is needed due to the potential deterministic effects. In digital radiology this dose information is often stored together with the clinical images in the PACS (Picture Archive and Communications Software). This can be done in the form of dedicated and sometimes private DICOM header values, printed as an image in SC (Secondary Capture) images, as SR objects (Structured Reports) or in log files stored on the console of the X-Ray device. Thus, the analysis of total patient dose, the monitoring of protocol quality, and the comparison of the used settings (physical-technical and configured) of different devices needs to be addressed, preferably within a single system.

19.2. Description of the solution

We aim to overcome the above problems with the software solution DOSE.

Using this software, it is possible to:

- analyze all types of clinical images for the received dose information
- store them in a separate database, apart from the PACS, for further future access
- Allow to view these data on patient level (historical overview of patient dose) and to compare the mean doses per study of the different X-Ray devices
- Analyze the workload/flow
- Compare the modalities

In addition, it is possible to:

- Analyze and track the different commonly used modalities (DX, CR, MG, CT, XA, RF, US, MR)
- Have a solution that is integrated into the clinical workflow, without impacting the normal workflow
- Have a vendor-neutral system
- Present this data to all stakeholders in the radiology department in the most optimal form. If needed, we will develop integration modules to bring this data closer to the end user and find new ways to present the collected data

- Extract other types of data related to the exposure parameters, in addition to the dose information, in order to simplify the understanding of what went wrong if dose values are exceeding limits
- Compare the extracted dose values with national/predefined reference values

19.3. Supported modalities

- ✓ General radiology (DR, DX, RF, CR)
- ✓ Computed tomography (CT)
- ✓ Interventional radiology (XA)
- ✓ Mammography (MG)

As an additional option, we can also analyze MR (including specific absorption rate, SAR) and US workflow and workload.

19.4. DOSE: Individual patient benefits

- DOSE can access all the available radiation dose information from PACS and from the devices directly, create reports that include previous examinations with the cumulative effective dose, calculated depending on age and body parts, and thus provide patient radiation dose history. The patient dose history can be used for individual risk assessment. Although this information is not used for justification of a future exam, it can be used to optimize the workflow for the patient benefit. If, for example, a department has two scanners and scanner 1 has a chest protocol with relatively lower dose than the same protocol on scanner 2, then, a patient of high cumulative dose that needs to do a chest CT can be sent to scanner 1. This feature can be also used to identify subgroups that receive high doses and pave the way towards a better clinical follow up protocol and/or specific device protocols to optimize their exposure without jeopardizing the benefits.
- The user can assess the radiation dose to the patient skin from a specific procedure or from multiple procedures. If this is high and a skin reaction is expected, the patient can be pro-actively referred to another specialist (e.g. dermatologist) to treat the specified deterministic effects as soon as possible. Additionally, high doses increase the possibility of radiation-induced stochastic effects (e.g. neoplasm), and this is also something that users may want to follow up with the goal to optimize and potentially reduce doses for next exams. The analysis of patient skin doses can help the users to optimize their practices for future patients. In cases that the complexity of the procedure allows it, different geometry and different angulations may be used to improve the outcome with respect to patient skin dose and minimize the risk for stochastic effects for future procedures.
- DOSE can make comprehensive reporting about patient organ doses. This can ensure that potentially pregnant women, pediatric patients, and other patients at

risk receive optimal doses for the clinical goal. The reporting of organ doses helps users to see the levels of exposure and evaluate the risk for stochastic effects for patients that undergo each exam type. Therefore, this allows to optimize the procedure as much as possible and minimize the risk for radiation-induced effects.

- The Live Dashboard feature in DOSE allows radiographers to see the given dose to the patient and any raised alerts immediately after the examination. Thus, the radiographer can contact the radiologist directly and discuss any potential problem/harm to the patient resulting from the examination, as well as image quality issues.
- The analysis of image noise and of the subjective image quality allows users to combine image quality with dose data in order to improve the practice based on the ALARA principle but having in mind also the adequate diagnostic result of the exam.
- Patient radiological data can be benchmarked against:
 - o Similar examinations: this is based on dynamic limits that help radiographers understand if the examination that the patient has just undergone, had a higher or lower dose if benchmarked against similar examinations from the same device or modality, having the same study description, composition, protocol name or body part.
 - o Ages/patient sizes: this comparison takes the age and size of the patient into account, which is important for the calculation of the individual patient dose. This is especially important for the pediatric patients and some body parts (e.g. breast, thyroid, uterus...) because they are more sensitive to radiation.
 - o Thresholds set by scientific guidelines or national authorities: By providing a report with this information and by generating notifications, DOSE increases awareness for the patient to ensure that the dose level is reasonable for this study type and patient characteristics (e.g. size). This way, the patient can be followed up, while errors can be found faster and ensure avoidance for future patients.
- Legal compliance: this is based on the static limits per study category that the national authorities provide each year. The exams are benchmarked automatically against exams of the same group, and if the overall group exceeds the defined thresholds, the system highlights it. By providing a report with this information, DOSE informs users to optimize this group for future patients and ensure that the level of dose is reasonable for this study type, following the ALARA (as low as reasonably achievable) principle.

19.5. DOSE residual risks

QAELUM is performing advanced risk reduction of all risks related to the software and its specific versions. However, the following list summarizes the **residual risks**:

1. **No access to, or, no/delayed/incomplete availability of patient exam dosimetry data** that may lead to missing potentially linked alerts. This can be due to:
 - a. Temporary inaccessibility of the system; e.g., due to lack of electricity, network problems or OS updates performed by local IT departments.
 - b. VM shutdown; e.g. due to human error or hardware failure.
 - c. Insufficient disk space; e.g., due to unexpected large amounts of data incoming from the PACS. Disk space is also used to collect thumbnails and quality control images: as the amount of analyzed data grows, QAELUM may ask the IT department to add disk space. If the IT department does not timely reply, the system may stop working.
 - d. Vendors or customers (i.e., users having access to the appropriate Functionalities) can change the device-data linking identifiers (e.g. StationName 0008,1010) without notifying QAELUM of this change.

2. **Inability to provide support by QAELUM.** This can be due to:

Temporary inaccessibility of VPN. The local IT department may not grant access to the DOSE server, so the QAELUM team may not be able to solve problems quickly.

3. **Unauthorized actions performed on the physical or virtual server(s) hosting the product**, without approval of QAELUM or 3rd party Supplier, can lead to disruption and incorrect functioning of the product. This can be due to:
 - a. Processes stopped. The local IT department has the possibility to connect to the DOSE server, as they are responsible for hardware and the operative system: they may stop or remove some QAELUM related services.
 - b. Database corruption. Often DOSE databases are administered by customers: therefore, they could use their access rights to accidentally remove data, modify the structure and/or change configurations, without first consulting with the Supplier.
 - c. Post-production installation of antivirus software (at customer server, by the local IT department and without informing QAELUM proactively) can be a problem for less common tools, like DOSE, as typically only executables known by the antivirus company are trusted.

Remark: If this is done prior to production, QAELUM will report on any negative effects of the antivirus on the performance of DOSE so the necessary actions can be taken prior to production use.

4. **Emails (notifications, reports, password resets, ...) not reaching the destination user.** This can be due to:
 - a. Incorrect set or outdated email on user profile or email misconfiguration.

- b. Customer email server that considers emails as spam (e.g. because of setting up email rules that generate too many emails in a predefined timeframe).
 - c. Changes in the smtp connection details without identifying QAELUM or 3rd party Supplier.
5. **Not following QAELUM's advice on backup procedures and time frames** can, in case of data corruption, lead to loss of data or delay in data recovery/catch up. The backup guidelines are part of the DOSE technical instructions/documentation, but the customer is responsible to set the configuration and maintenance plan, corresponding to local needs and policies. Whenever the database gets corrupted and the backup is not present, then some data may be lost.
6. Using the product interface on a **monitor with resolution requirements below 1024x768 or using unsupported browsers** can lead to sub-optimal display/user experience and possibly wrong actions or conclusions. For more information, please check the DOSE technical instructions/documentation.
7. **Giving users access to Functionalities with which they are insufficiently familiar** (e.g., lack of knowledge, training, etc.), can lead to modifications that in turn can lead to data modification, corruption, errors, or loss. For example, users with administrative roles may remove data by mistake and/or may misconfigure devices by adding the wrong device-data linking identifiers (e.g. StationName 0008,1010). Please assign Functionalities with care and provide user training. We expect users to have sufficient understanding of radiology, anatomy, and radiation, to avoid incorrect action or drawing inappropriate conclusions. Dose indicators should always be used according to legal and scientific guidelines. **This product should not be used outside of its intended use.**

19.6. DOSE intended user and training

DOSE is a state-of-the-art software offering a vast variety of functionalities. This software is intended for use by medical professionals only – radiographers, radiologists, the management of radiology departments, medical physicists. Therefore, software training by an expert from QAELUM is required before users can have access and start using the software.

19.6.1. Intended Use

DOSE is intended to extract, analyze and monitor patient radiation dose of the individual patient who undergoes X-ray examinations. The intended users are the medical professionals in radiology departments or diagnostic imaging facilities. Its use is independent of patient disease or condition, frequency of imaging, patient characteristics, or the imaged anatomy. It can assist the user in evaluating dose levels and potential improvements for individual patient examinations, and thus improve diagnosis and determine optimal treatments. Additionally, it helps the user to manage all the

patient radiation related data that is available in the imaging department in a faster and more efficient way.

The software allows:

IU#	Intended use
IU1	The software extracts radiation dose and other parameters for every acquisition or examination of the patient's medical exposure. This will allow the end-user to have easy access to all the details of the acquisition/examination and evaluate the situation for the patient.
IU2	The software calculates dose and dose-related indices for every acquisition/examination taking into account as many details as possible. This helps to provide a more individualized assessment for each patient and helps identifying the shortcomings with respect to procedures, techniques, workflows, persons for further improving the patient quality of care.
IU3	The extracted and/or calculated dose and dose-related indices are used to manually and/or automatically evaluate and assess the level of exposure of the patient. Based on these numbers the end-user can decide whether or not to take actions to optimize exams for all future patients.
IU4	Calculating the local reference levels and comparing with national and international DRL values allows to monitor the exposure levels for each exam type on the specific device/hospital and optimize when needed. This will allow for optimized exams for all future patients.
IU5	Benchmarking between similar studies allows for better evaluation of the practice and optimization of the technique/dose/performance and therefore outcome for patients that undergo the same exam.
IU6	After the assessment of dose and dose-related indices, the system notifies the end-user for potential outliers or situations that can have an impact on diagnosis or patient's health
IU7	Analyses and assesses the cumulative dose and dose-related indices over all examinations of the patient.
IU8	The system collects, analyses and assesses the image quality of a patient's exam in order to ensure proper diagnostic exam quality for current and future patients
IU9	Sharing or exporting patient exam data to facilitate data access without the need to connect directly to our system.
IU10	Evaluation and assessment of parameters related to device usage and performance which will benefit the patient by facilitating the optimization of the clinical imaging pathway

19.6.2. Medical conditions

DOSE is used after the medical exposure is finished, and there are no medical conditions that do not allow the use of the product.

19.6.3. Contraindications

There are no contraindications. DOSE patient radiation dose analysis can be calculated for any patient who undergoes a medical imaging examination.

19.6.4. Performance of the device

- Accuracy of measurement(s)

The accuracy of the measurements and calculations that are performed in DOSE depends on the modality and the dosimetric parameter. In some cases, it is also linked to other parameters like the examined body part or the protocol used.

In all cases, when percentiles are calculated, the minimum value of the parameter is expected to be smaller than the maximum and the 25th, 50th (median) and 75th percentiles should be values between the minimum and maximum. The mean value is expected to be $(\text{Max}-\text{Min})/2$.

These rules should apply, unless there is only 1 study taken into account and then $\text{Minimum}=\text{Maximum}$ and all the other percentiles are the same.

Additionally, if a study is selected to be ignored from the statistics (manually, by the user), then the system should exclude it from the analysis of the particular study group/study description/protocol. In this case, the statistics/percentiles are expected to be different than before.

Below an analysis per modality follows. Most of the values mentioned are indicative and should not be evaluated for an individual study, but as a measure for the majority of studies and the expected order of magnitude.

Computed Tomography (CT)

The unit of CTDIvol in DICOM standards is mGy, and therefore is mGy in DOSE. The expected range for CTDIvol of one acquisition 0-150 mGy. Even though 150mGy is quite a high number, it can occur in a protocol that does not comply. Also, in the case of stationary acquisitions, if the operator scans using long exposure time, then the CTDIvol can be high. Thus, individual studies may reach these values, and then they are considered outliers and the user must justify them. To consider the reading of DOSE accurate, the majority of studies (excluding the aforementioned outliers), should have a CTDIvol that does not exceed the aforementioned values. The same stands for DLP (suggested range 0-2000 mGy•cm for an acquisition).

DOSE also calculates SSDE in mGy. The expected range of SSDE is $0.22 \cdot \text{CTDIvol} < \text{SSDE} < 2.76 \cdot \text{CTDIvol}$ for one acquisition, as these correction factors (0.22 for a very big patient and 2.76 for a very small patient) are the extreme values in the tables in the AAPM reports 204 and 220.

For patient positioning, the values are shown in mm. Vertical or horizontal offsets should be in the range of 0-150 mm as usually there is no space for the patient to be positioned in a larger offset.

For the in-house organ dose calculation, depending on the body part examined, the highest dose values are expected in specific organs:

Chest	Breast-Lung-Heart-Thymus-Thyroid
Abdomen	Colon-Kidneys-Liver-Spleen-Stomach
Chest-abdomen	Breast-Lung-Liver-Spleen-Thymus
Pelvis	Muscle-Prostate/Uterus-Small Intestine-Testes/Ovaries- Urinary bladder
Abdomen-pelvis	Kidneys-Liver-Spleen-Stomach-Urinary bladder
Chest-abdomen-pelvis	Breast-Lung-Liver-Spleen-Thymus

Interventional Radiology (XA)

Peak skin dose map should have colors on the dark background in order to check that a calculation is performed. Also, the colors should form a shape on the map and the dark background should be visible, at least in a small part of the map.

The peak skin dose value is reported in Gy. Thus the expected range is 0.01 – 30 Gy. A PSD of 30 Gy can happen and does not indicate a calculation problem in DOSE.

The same stands for the Dose at Reference Point, which should have the same range.

Conventional radiography (DR-DX-CR-RX)

For conventional radiology, there are two parameters that should be evaluated, DAP and exposure index. However, the exposure index is not always used by the manufacturers with the standardized calculation methodology. If the standardized exposure index is used, then it is normally expected <1000 for an image (to be checked further).

DAP for an image should normally be <100 dGy•cm².

Mammography (MG)

Organ dose which indicates the mean glandular dose of the breast is in mGy.

Organ dose for each series should be below 10 mGy. This is a large number, but it is indicative of the expected order of magnitude. Cumulative organ dose of a study should be below 50 mGy.

- Capacity performance

DOSE is designed as a modular, easy to scale-up system: each component can be replicated several times in order to increase the system capacity. Based on the expected

volumes, Qaelum provides customers with an architecture proposal: this is meant to guarantee availability and speed of the system.

DOSE is able to deal with projects for up to 20 different hospitals, 300 devices, yearly study volume of 3M, and 200 users. Bigger projects are also manageable, and scalability improvements are developed regularly.

- Speed performance

Qaelum team configures each DOSE system in order to provide customers with the best user experience. Speed is granted by using the correct architecture, hardware, and database, as discussed during the project technical meetings. Improvements are regularly deployed together with the software updates. Hardware usage is evaluated regularly, and requests for upgrades are made in case of need.

During trainings users are suggested to improve their experience by using broader date ranges (>6 months) only in case of need. A shorter date range (e.g. 3-4 weeks) is usually adequate to evaluate the current trends. As a general rule, the longer the date range, the longer the loading time: while showing graphs for one device with 2 weeks of data takes less than one second, showing several years (or tens of thousands of studies) may take up to 10 seconds.

Besides this, the loading of some components is naturally slower because of the amounts of calculations to be performed: for example, the legal compliance can usually take up to 5 seconds to be shown, depending on the amount of configured study groups.

Lastly, saving configurations that require the recalculation of data (for example the study group configuration) may take several seconds to be completed, depending on the amount of studies.

- Compatibility with browsers and operating systems

DOSE is a zero-footprint software web application. This means that it requires no client-side application install and has no OS dependency. The application can be accessed from anywhere within the hospital network by browsing to the URL set up for DOSE. We support the following browsers:

- ✓ Google Chrome 80+
- ✓ Mozilla Firefox 70+
- ✓ Internet Explorer 11
- ✓ Microsoft Edge

20. Incident reporting

In case any serious incident has occurred in relation to the device, please report this to Qaelum. You can do this by phone, mail or through our ServiceDesk portal. The details for these ways of communication are available in section 16 of this document.

In case of need Qaelum will contact the competent authority of KINGDOM OF SAUDI ARABIA. You are allowed to contact the competent authority yourself, but always first inform Qaelum so we are also aware of the incident.

21. Glossary

- **Accession number** unique identifier assigned to a study within the PACS
- **ACC** acceptable
- **ACH** achievable
- **ADT** Admission, Discharge, Transfer HL7 message
- **AEC** Automatic Exposure Control
- **AFCN** Agence Fédérale de Contrôle Nucléaire
- **AGD** Average Glandular Dose
- **CBCT** Cone Beam Computer Tomography
- **Collimation width** width of the collimation to restrict radiation to the required dimensions
- **Compression force (N)** force of breast compression (in Newtons)
- **CR** Computed Radiography
- **CSV** Comma-Separated Values
- **CT** Computer Tomography
- **CTDI** Computed Tomography Dose Index
- **DAP** Dose Area Product
- **Device** equipment unit used for imaging
- **DICOM** Digital Imaging and Communications in Medicine
- **DLP** Dose Length Product
- **Effective dose** sum of equivalent doses to all organs, accounting for their respective radiation sensitivity (in mSv)

- **Exposure time** length of time that the patient is exposed to radiation
- **FANC** Federaal Agentschap voor Nucleaire Controle
- **Flow rate** the volume of fluid that passes per unit time
- **Catheter gauge** plate consisting of holes of various size used to determine the size of a catheter
- **ICRP 106** International Commission on Radiological Protection publication 106
- **Isocentre** the point in the patient through which the central ray of the radiation beam passes
- **JPEG** Joint Photographic Experts Group
- **mA** milliamperage i.e. tube current
- **mAs** tube current-exposure time
- **MG** mammography
- **Modality** a type of equipment group used for imaging
- **MPPS** Modality Performed Procedure Step
- **MR** Magnetic Resonance
- **Organ dose** energy absorbed per unit mass in an organ
- **ORM** Order Request Message in HL7
- **PACS** Picture Archiving and Communication System
- **PDF** Portable Document Format
- **Kilovoltage peak (KVp)** peak potential applied to the x-ray tube, affecting the quantity and quality of photons
- **Pitch** table distance travelled in one rotation of the gantry
- **PNG** Portable Network Graphics
- **PSD** Peak Skin Dose, the highest dose to a single area of the skin
- **RDSR** Radiation Dose Structured Report
- **RX** general radiography
- **SAR** Specific Absorption Rate, rate of energy absorbed per unit mass by a human when exposed to a radio frequency field
- **Scanlength** length of table travel during a CT scan
- **SSDE** size-specific dose estimates
- **Study ID** study identification

- **Study UID** study unique identification
- **SVG** Scalable Vector Graphics
- **US** Ultrasound
- **XA** DICOM modality that includes x-ray angiography, interventional radiology
- **XLS** excel Spreadsheet